



AVANOS*

MIC-KEY* GJ FEEDING TUBE

USE AND CARE INSTRUCTIONS

MIC-KEY*

⌀ Diameter	↔ Stoma Length	↔ Jejunum Length	Balloon Volume	⊗ Single Use Only
STERILE R Sterilized by Gamma Irradiation	⊗ Do not use if package is damaged	⊗ Do not resterilize	Not made with natural rubber latex	
Product is NOT made with DEHP as a plasticizer	Rx Only	MR Conditional	⚠ Caution	📖 Consult instructions for use

AVANOS* MIC-KEY* Gastric-Jejunal Feeding Tube (Low-Profile GJ-Tube)

Tube Use and Care Instructions

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Disposable, for single patient use only. Do not resterilize.

Product Description

The MIC-KEY* Low-Profile Gastric-Jejunal Feeding Tube is a single unit feeding/drainage device (**Fig. A**). It enters the stomach through a gastric stoma. The tube is held in place (within the stoma tract) by means of an inflatable balloon (**Fig. A-1**) and a low-profile external bolster (**Fig. A-2**). The external bolster contains two ports; one labeled "JEJUNAL" and one labeled "GASTRIC". Incorporated within each port is a one-way valve which is opened by attaching the appropriate MIC-KEY* Extension Set. MIC-KEY* Extension Sets (**Figs. B & C**) connect to the feeding and drainage ports. The JEJUNAL port is used for feeding into the small intestine. The GASTRIC port is used to drain the stomach with the use of low intermittent suction or gravity drainage. A third port labeled "BAL" is used to inflate and deflate the balloon.

Indications for Use

The use of this tube is clinically indicated when simultaneous gastric decompression and jejunal feeding are needed.

For enteral nutrition only.

⚠ Warning

Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

Feeding and Decompression through the MIC-KEY* Gastric-Jejunal Tube

The following information explains how to properly connect the MIC-KEY* Extension Sets and to feed into the small intestine and decompress the stomach.

How to Connect the Extension Set for Jejunal Feeding

1. Remove the feeding port cover (**Fig. A-3**) from the top of the MIC-KEY* Gastric-Jejunal Feeding Tube.
2. Insert the MIC-KEY* Extension Set (**Fig. B**) into the port labeled "JEJUNAL" by aligning the lock and key connector. Align the black orientation marking on the set with the corresponding black orientation line on the feeding port.
3. Lock the set into the jejunal feeding port by pushing in and rotating the connector CLOCKWISE until you feel a slight resistance (approximately a 3/4 turn).

Note: DO NOT rotate the connector past the stop point.

4. Remove the extension set by rotating the connector COUNTER-CLOCKWISE until the black line on the set aligns with the black line on the feeding port. Remove the set and cap the MIC-KEY* Gastric and Jejunal ports with the attached feeding port cover.

Jejunal Feeding

⚠ Warning: Never connect the jejunal port to suction. Do not measure residuals from the jejunal port.

⚠ Caution: Do not administer medications through the jejunal feeding port. This will clog the tube. To help avoid clogging, use a feeding pump to deliver formula through the jejunal lumen. Never allow formula to stand in the tube.

1. Assemble the equipment: Formula, Feeding set, Catheter-tip irrigating syringe, Enteral feeding pump, Water for flushing the tube, Extension Set.
2. Wash your hands with soap and water. Rinse and dry them thoroughly.
3. Shake the formula and wipe the container tops before opening. If you do not use all the formula, cover the open container, record the date and time, and store it in the refrigerator. Throw away formula more than 24 hours old. Do not mix new formula with old formula. There is always a chance it will spoil.
4. Pour the formula into the feeding bag.
5. Connect the MIC-KEY* Extension Set (**Fig. B**) to the port labeled "JEJUNAL". Open the clamp on the extension set.
6. Inject 30 ml water into the MIC-KEY* Extension Set by using the 35 ml irrigating syringe from your kit.

7. Attach the extension set to the feeding set connector. Connect the feeding set tubing to the enteral feeding pump. Follow the manufacturer's instructions to set the pump.
8. If the feed set tubing has a clamp, open it.
9. Start the feeding pump.
10. Flush the jejunal port EVERY SIX HOURS WITH AT LEAST 30 ml WATER. DO NOT USE FORCE.
Note: If you are simultaneously decompressing the stomach while feeding and you see formula in the gastric drainage, stop the feeding and notify the specialist.
11. When the feeding is complete, clamp the feeding set tubing and disconnect the feeding set tubing from the extension set.
12. Flush the extension set and feeding tube thoroughly. Disconnect the extension set from the feeding tube port.

How to Connect the Extension Set for Gastric Decompression

1. Remove the feeding port cover (**Fig. A-3**) from the top of the MIC-KEY* Gastric-Jejunal Feeding Tube.
2. Insert the MIC-KEY* Bolus Extension Set (**Fig. C**) into the port labeled "GASTRIC" by aligning the lock and key connector. Align the black orientation marking on the set with the corresponding black orientation line on the gastric port.
3. Lock the set into the gastric decompression port by pushing in and rotating the connector CLOCKWISE until you feel a slight resistance (approximately a 3/4 turn).
Note: DO NOT rotate the connector past the stop point.
4. Remove the extension set by rotating the connector COUNTER-CLOCKWISE until the black line on the set aligns with the black line on the feeding port. Remove the set and cap the MIC-KEY* Gastric and Jejunal ports with the attached feeding port cover.

Gastric Decompression

Some specialists recommend decompressing the stomach (letting the air or fluid out) before or after feeding. FOLLOW YOUR SPECIALIST'S INSTRUCTIONS. Connect the MIC-KEY* Bolus Extension Set to the port labeled "GASTRIC". The end of the Extension Set may be left open for gravity drainage or connected to low intermittent suction. This allows stomach contents and/or gas to escape. Flush the gastric port EVERY SIX HOURS with at least 30 ml water. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.

Diarrhea

- The flow rate may be too fast. Decrease the rate, then increase it slowly until you reach the prescribed flow rate. Diluting the formula with water may help slow the diarrhea, but check with the specialist first.
- Prepare fresh formula for each feeding. Discard unused formula after 24 hours.
- Keep your equipment clean. Prepare only enough formula for eight hours. Wash, rinse and dry the equipment after each use.
- If the diarrhea continues for more than three days despite your efforts, call the specialist.
- Never attempt to bolus feed through the jejunal lumen.

Stoma Care

Cleanse the site with normal saline three times daily to remove the small amount of mucus that normally accumulates around the stoma. If the mucus dries, it may irritate the skin. Applying a small amount of water may loosen the dried material. After the stoma heals a thorough cleansing with soap and water is best. A dressing is not necessary, and contributes to skin maceration and infection. Unless recommended by your specialist, avoid medicated ointments or powders. DO NOT use mineral oil or petroleum based products. DO NOT ROTATE THE TUBE. Jejunal tubes must NOT be rotated because they will kink. Should the tube become kinked, call your health care professional.

Tube Care & Maintenance

1. **Balloon Maintenance** Check the water volume in the balloon ONCE A WEEK. While holding the external bolster in place, insert a luer slip syringe into the balloon inflation port. Withdraw the water from the balloon. The amount of water you withdraw should equal the prescribed amount. If you withdraw less water than the prescribed amount, replace the amount recommended by your specialist. Never add water before you withdraw all of the water from the balloon.

Note: The use of a luer-lok syringe may make this more difficult. A luer-slip syringe is recommended.

- Never add more than 5 ml to the balloon for pediatric sized tubes (REF numbers ending in -15, -22, or -30).
 - Never add more than 10 ml to the balloon for adult sized tubes (REF numbers ending in -45).
2. **Tube Blockage** Obstructions are usually caused by formula accumulation inside the tube. Body fluids and medications can also clog the tube. Prevent this by thoroughly flushing the tube EVERY SIX HOURS WITH A MINIMUM OF 30 ml WATER. If the flow is sluggish, gently inject the water, and let it flow back out of the tube until the obstruction has cleared. When you feel too much resistance and cannot inject water at all, the tube may need to be repositioned or replaced because it may be kinked. High pressure can rupture the tube. Do not insert foreign objects into the tube. This may puncture the tube. If none of these measures work, call your specialist.
 3. **Granulation Tissue** A small amount of extra tissue around the stoma site is simply the body's attempt to heal the wound. This is normal. It is also normal for the extra tissue to secrete a small amount of mucus. If the tissue bleeds excessively, drains, or interferes with the external bolster, call your specialist.
 4. **Unusual Symptoms** Consult your health care practitioner if any of the following symptoms occur: Abdominal pain • Abdominal discomfort • Abdominal tenderness • Abdominal distension • Dizziness or fainting • Unexplained fever • Unusual amount of bleeding through or around the tube.

MRI Safety Information

Non-clinical testing has demonstrated the Low-Profile (MIC-KEY®) Enteral Feeding Tube System is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3 Tesla;
- Maximum spatial field gradient of 1,960 G/cm (19.6 T/m) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of < 2 W/kg (Normal Operating Mode).

MRI-related heating: Under the scan conditions defined above, the Low-Profile (MIC-KEY®) Tube System is expected to produce a maximum temperature rise of less than 1.3 °C after 15 minutes of continuous scanning.

Artifact Information

In non-clinical testing, the image artifact caused by the device extends less than 45 mm from the Low-Profile (MIC-KEY®) Enteral Feeding Tube System when imaged with a gradient echo pulse sequence and a 3 T MRI system.

Sonde d'alimentation gastrique-jéjunale MIC-KEY* de AVANOS* (sonde GJ extra-plate)

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*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2020-10-08
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.
Sponsored in Australia by Avanos Medical Australia Pty Ltd, 475 Victoria Avenue,