

AVANOS*

MIC* GASTRO-ENTERIC TUBE

Patient Care Instructions

MIC*

 Diameter	 Length	 Single Use Only
STERILE R Sterilized by Gamma Irradiation	 Do not use if package is damaged	 Do not resterilize
Rx Only	 Caution	 Consult instructions for use

MIC* Gastro-Enteric Tube

Patient Care Instructions

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Patient Information

A specially-designed MIC* Gastro-Enteric Tube has been inserted into your stomach through the abdominal wall. It delivers nutrition beyond the stomach into the small intestine. This tube permits the intake of the food and water that your body requires. The tube also will remove any gas or secretions which may build up in your stomach. Your doctor has selected the proper tube for you and has given you instructions on its care and maintenance. This booklet will help you remember those instructions. You have also been given instructions about diet and medications to be administered through the tube. Follow these instructions closely and never use any other diet or medication.

Key Terms

- A Jejunal Feeding Port
- F Gastric Exit Ports
- B Gastric Decompression Port
- G Radiopaque Jejunal Lumen
- C Balloon Inflation Port
- H Jejunal Exit Ports
- D SECUR-LOK* External Retention Ring
- I Weighted Tip
- E Internal Retention Balloon

Warning

You or your caregiver should not attempt to remove or change the MIC* Gastro-Enteric Tube.

If the skin around the tube appears red or irritated, consult your physician at once.

Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

 **Caution:** Follow your physician's instructions for proper gastric maintenance. If you are on continuous feeding, the tube should be flushed at least 3–5 times daily. Avoid puncturing any part of the MIC* Gastro-Enteric Tube.

Care and Use

The MIC* Gastro-Enteric Tube must be cared for daily. Care is simple and easy, consisting mainly of keeping the tube and the skin around the tube (stoma) clean and dry.

Keep the following supplies on hand to make your cleaning task easier and remember to wash hands before touching your tube.

- Soap and water
- Cotton-tipped swabs
- Tissues or 4 x 4s
- Diet soft drink for use as a cleaning solution
- Disposable syringe

Inspect both the skin around the tube (the stoma) and the tube routinely at feedings and make sure they are clean.

Lift the edges of the ventilated SECUR-LOK* ring and gently clean the skin around the tube (stoma) with cotton-tip swabs, a soft cloth, or 4 x 4s, using soap and water. Rinse and dry the stoma after cleaning.

Inspect and Clean the Tube

Clean and dry the outside of the tube and both sides of the ventilated SECUR-LOK* ring using cotton-tipped swabs. The ring can be cleaned in place or can be moved up the tube to clean.

Flush the Tube

After every feeding, flush the jejunal port of the tube with 20–30 ml of water (adults) or 10–15 ml for children or infants. If you are on continuous feeding, the tube should be flushed at least 3–5 times daily. The tube should be clear of medicine or feeding solutions. Wipe down the feeding port and its plug using a swab, soft cloth, or 4 x 4s to remove any build-up of oils or feeding materials.

The gastric port should also be irrigated on a regular basis to prevent any accumulation of stomach secretions which could block the tube.

Problem Solving

Problem	Cause	Remedy
1. Gastric leakage.	The SECUR-LOK® Ring is not properly adjusted or the balloon is not against the stomach wall. The balloon is under-inflated.	Check position of the ventilated SECUR-LOK® ring: 1–2 mm from the stoma, and the balloon is against the stomach wall. If the SECUR-LOK® ring is properly adjusted, add water in small increments (1–2 ml) to the balloon. Water can be added as the stomach size and balloon volume will allow. Note: 20 ml maximum for standard volume balloon inflation, 5 ml maximum for low volume balloon inflation. Make sure that the stoma is clean and dry. Wait 15 minutes and check to see if the problem has been corrected. If the problem persists contact your physician.
2. Feeding tube disconnects from the feeding port.	Build-up of oily deposits inside the tube.	Keep the feeding port connections clean. Remove oily build-up with either soap and water or a diet soft drink. Use a small amount of solution on a cotton tip swab and scrub the top two inches of the inside of the tube.
3. Tube decreasing in length.	Oily build-up on tube resulting in the SECUR-LOK® ring not holding.	Clean the tube. Wash with warm, soapy water; rinse and dry the tube. If oil remains, clean the tube with a small amount of diet soft drink. After cleaning, readjust the ventilated SECUR-LOK® ring 1–2 mm away from the stoma site, and verify that the balloon is up against the stomach wall. Check in one hour to be sure the ventilated SECUR-LOK® ring is holding the MIC® Gastro-Enteric tube in place. If size continues to change, consult the physician.
4. Tube increasing in length.	Balloon under-inflated or balloon is leaking.	Withdraw water from the balloon and compare the volume against the amount originally used to fill the balloon. Fill again to the original volume. Recheck the volume periodically. If size continues to change, call your physician immediately.
5. Patient vomits and/or aspirates feeding formula.	Jejunal feeding tube may have dislodged. Formula rate may be set too high.	Contact your physician immediately.
6. Formula will not flow freely through feeding port.	Jejunal tube may be kinked or clogged.	Irrigate with 20–30 ml of warm water (adults) or 10–15 ml (children). If resistance is felt or water does not flow freely, call physician immediately.
7. Gastric lumen not vented properly.	Possible obstruction in gastric lumen.	Irrigate with 20–30 ml of warm water (adults) or 10–15 ml (children). If resistance is felt or water does not flow freely, call physician immediately. Tube may need to be replaced.

DO'S

- **DO** report unusual symptoms such as abdominal pain, abdominal discomfort, abdominal tenderness, abdominal distension, unexplained fever or an unusual amount of bleeding through or around the tube to your health care practitioner.
- **DO** wash your hands before handling the tube or your stoma.
- **DO** irrigate the tube with each use—no less than 3–5 times a day.
- **DO** keep the feeding port clean.
- **DO** keep the exterior of the tube clean.
- **DO** adjust the SECUR-LOK® ring properly (1–2 mm from the stoma site).
- **DO** use commercially prepared liquid medications.
- **DO** keep the gastric port clear.
- **DO** use a feeding pump.

DON'TS

- **DON'T** secure the tube in place with tape or other appliances.
- **DON'T** let feeding solutions stand in the tube.
- **DON'T** use any diet not prescribed by the physician.
- **DON'T** cover the SECUR-LOK® ring or healthy skin with tape or gauze.
- **DON'T** fill the balloon with air.
- **DON'T** push the SECUR-LOK® ring onto the abdomen.
- **DON'T** put pressure on the stoma.
- **DON'T** clamp the tube.
- **DON'T** use oil-based lubricants or petroleum jelly on the tube.
- **DON'T** use crushed medication.
- **DON'T** use blenderized diets.

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*Registered Trademark or Trademark of Avanos Medical, Inc.,
or its affiliates. © 2018 AVNS. All rights reserved. 2019-01-24
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc., 5405 Windward Parkway,
Sponsored in Australia by Avanos Medical Australia Pty Ltd,