NEOMED[®]

ENTERAL EXTENSION SET


WITH ENFit[®] CONNECTOR

- FR

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EN

Instructions for Use

Rx Only Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Product is not made with natural rubber latex.

Product is not made with DEHP as a plasticizer.

Product is not made with BPA.

This instructions for use applies to all NEOMED[®] Enteral Extension Sets with ENFit[®] Connector.

Notice: This device is not intended to be used as a standalone device and can only be used with ENFit[®] Feeding Tubes. Read this complete use guide and the Feeding Tube Instructions for Use (IFU) prior to using this device. For additional warnings and precautions, please see IFU of the feeding tube.

Device Description / Intended Purpose

The NEOMED[®] Enteral Extension Set is an accessory device used with NEOMED[®] Feeding Tubes to facilitate delivery of enteral nutrition, water, and medications.

Length	Medical device	Single Use Only
Sterilized using ethylene oxide	Single sterile barrier system	
Single sterile barrier system with protective packaging outside	Do not use if package is damaged	Do not resterilize
Rx Only	MR Safe	Caution
Consult instructions for use	Choking Hazard - Small Parts	Keep dry
Keep away from sunlight	Country of manufacture	Manufacturer
Unique Device Identifier	Batch code	Date of Manufacture
		Use-by date

Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway, Alpharetta, GA 30004 USA. In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com

Avanos Medical, Inc.,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.

Sponsored in Australia by Avanos Medical
Australia Pty Ltd, 475 Victoria Avenue,
Chatswood, NSW 2067, Australia.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2024-03-04 15-M1-1983-01

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Features

The Enteral Extension Set is provided in a variety of configurations:

- Standard extension sets have one proximal and one distal ENFit[®] connector. ENFit[®] connectors can be male or female.
- Bifurcated extension sets allow for medication or additional nutritional delivery without disrupting the feeding line.
- Bifurcated extension sets provide an additional proximal ENFit[®] connector.
- Large and extra large bore extension sets facilitate delivery of thicker (more viscous) nutrition.

Intended Users

This device is intended to be used by healthcare professionals.

Patient Target Group

This device is intended for use in neonatal and pediatric patients.

Indications for Use

The NEOMED[®] Enteral Extension Set is indicated for use in neonatal and pediatric patients as an extension set for nasogastric/oral gastric or gastric tube enteral feeding tubes, incorporating safety connectors which help mitigate the risk of accidental connection of an I.V. system to the enteral system or the enteral system to an I.V. system.

The Enteral Only Extension Set is a sterile disposable for single patient use only.

Contraindications

The NEOMED[®] Extension Set is intended to connect to products with ENFit[®] compatible connectors. The Extension Set will only connect to ENFit[®] enteral feeding tubes and ENFit[®] enteral syringes.

The Extension Set is not intended for individual use longer than 24 hours. Replace and maintain extension set in accordance to hospital protocol and policy/or reference AND (Academy of Nutrition and Dietetics) guidelines, as applicable.

Residual Risks and Side Effects

While side effects vary depending on the type of procedure performed, mitigations have been implemented to minimize the following residual risks related to the use of the NEOMED[®] Extension Sets and the associated NEOMED[®] Nasogastric Feeding Tubes: abnormal reproductive development, allergic/systemic toxic reaction, over/under delivery of medication (potentially life-threatening)/nutrition, skin/tissue damage/irritation (includes trachea, vocal cords, gastrointestinal tract/ulcer/necrosis), infection/fever/sepsis, nausea/vomiting/bloating/diarrhea, delayed or prolonged treatment/procedure/therapy, pneumothorax, choking, long-term carcinogenic and/or mutagenic effects, aspiration, peritonitis, gastroenteritis, laceration, pain.

⚠ Warnings/Cautions

Do not reuse devices labeled for single use only. Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.

Do not use if primary pouch has been previously opened or damaged as this may lead to device contamination or biological infection.

Failure to read entire instructions for use may result in potential patient harm, off-label use of the device and/or nutritional fluid loss.

This device is intended for one-time use. Each device should be used on only one patient.

Do not autoclave or re-use device as this may increase the risk of leakage or biological infection.

Use of the device in non-enteral applications may lead to severe patient harm or death.

Do not use if primary pouch has been previously opened or damaged as this may lead to device contamination or biological infection.

Do not use visibly damaged device as this may lead to inadequate nutrition delivery or potential patient harm.

Do not use forceps, clamps other than included with extension sets, or sharp instruments as this may damage the device and lead to inadequate nutrition delivery.

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Do not modify device as this may lead to inadequate nutrition delivery or severe patient harm.

Use of the device for longer than 24 hours may result in leakage, inadequate nutrition delivery, device contamination or biological infection.

Obstruction or inadequate clamping may lead to blockage, leakage or inadequate nutrition delivery.

Failure to secure/use distal cap or proximal tethered plug may result in leakage, inadequate nutrition delivery, device contamination or biological infection.

Bifurcated Extension Set: Failure to close slide clamp and plug female tubing when not in use may result in nutritional fluid loss.

The distal hub cap is a choking hazard. Keep out of reach of children.

MR Information

MR Safe

The NEOMED[®] Enteral Extension Set contains only nonconducting, nonmetallic and non magnetic items and poses no known hazards in all MR imaging environments.

Device Lifetime

Extension sets are validated for use up to 24 hours.

Instructions

Follow facility protocol or current AND guidelines as applicable for use of this device.

- Always follow physician's instructions for use of contents.
- Wash hands per facility protocol.
- Remove extension set from packaging.
- Verify feeding tube compatibility prior to use.
- Remove and dispose of distal cap according to facility protocol.
- Attach Enteral Only Connector (inlet) to outlet of administration set. Attach extension set outlet to enteral only feeding tube.
- Close slide clamp and plug all female hubs until ready to use or when not in use.
- Replace and maintain administration set in accordance to facility protocol and not to exceed NANN (National Association of Neonatal Nurses) or AND guidelines of beyond 4 hours, as applicable.
- Do not re-use extension sets. After use, dispose of device per facility protocols and/or local or national regulations.

Complaint and Adverse Event Reporting

If a serious incident has occurred in relation to the device, it should be reported to the manufacturer and the competent authority in which the user and/or patient is established. A serious incident means any incident that directly or indirectly led, might have led or, in case of recurrence, could lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's, fetus or other person's state of health, or a serious public health threat.

Disposal

Contact the local authorities to determine the proper method of disposal of potentially biohazardous parts and accessories. Disposal of the device and accessories should be in accordance with local, state, and national laws and regulations.

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