EN

ECG Electrodes
Please read this document and keep it safe. Ensure that any individual using this product is familiar with and understands the information contained in this document.

INTENDED USE
ECG electrodes are accessories for ECG monitoring devices. ECG electrodes are non-sterile, intended for single use on intact (uninjured) skin and are used to pick up ECG signals from the patient.

- CAUTION**
- Do not use once product has expired.
 - To minimize the risk of burns at electrode placement sites during electrosurgery always follow the procedure specified by the manufacturer of the electrosurgical or monitoring unit!
 - Do not place the electrodes on areas of skin with inflammations, lesions, or injuries of any kind, or (for example) on skin stressed or damaged due to the use of medication!
 - Do not use the electrodes with direct electrode contact on patients with a history of skin irritation, and report any incompatibility reactions to the specialist dealer!
 - Do not use electrodes if their gel has dried out!
 - Do not use the electrodes longer than the time indicated on the packaging pouch!
 - Only use each electrode once. In the event of reuse, the adhesive strength and electrical properties may be inadequate. There is also a risk of cross-infection from one patient to another.
 - Use only electrodes that are suitable and appropriately labeled for X-ray imaging. These are marked with the “Radiolucent” symbol on the packaging pouch.

MRI safety information for MR conditional electrodes
Non-clinical testing has demonstrated that ECG electrodes are MR conditional. A patient with this device applied can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 teslas and 3 teslas, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T2/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of
 - 2 W/kg (Normal Operating Mode), or
 - 4 W/kg (First Level Controlled Operating Mode)


After 15 minutes of continuous scanning, under the scanning conditions defined above, the ECG electrodes are expected to produce a maximum temperature rise of less than

- 1.7 °C (2 W/kg, 1.5 teslas) or 3.4 °C (4 W/kg, 1.5 teslas) HF-related temperature increase with a background temperature rise of approx. 1.2 °C (2 W/kg, 1.5 teslas) or 2.5 °C (4 W/kg, 1.5 teslas),
- 1.3 °C (2 W/kg, 3 teslas) or 2.7 °C (4 W/kg, 3 teslas) HF-related temperature increase with a background temperature rise of approx. 0.9 °C (2 W/kg, 3 teslas) or 1.8 °C (4 W/kg, 3 teslas).

In non-clinical testing, the image artefact caused by the device extends approximately 4.79 mm from the ECG electrodes when imaged with a gradient echo pulse sequence and a 3 tesla MRI system.


- INSTRUCTIONS FOR USE**
The electrodes must only be used by a physician or a suitably trained member of medical staff.
- Select the placement sites that are favorable for your diagnostic or monitoring application.
 - Make sure the placement sites are clean, dry, and free of hair. (1)
 - Do not use solvent-based liquids to clean the skin because such agents can lead to skin reactions when trapped under the electrode.

When using Snap Connect electrodes and electrodes with a connection for a 4 mm connector:

- Connect the leads to the electrodes (2).
- Detach the electrode from the protective cover (3) and position the electrode on the skin (4).
- With circular finger movements, affix the base material of the electrode to the skin until it adheres entirely (5).

When using electrodes with the connection with tabs:

- Detach the electrode from the protective cover (2) and position the electrode on the skin (3).
- Attach clip connectors to the electrodes. The connection with tabs should point towards the leads to prevent twisting and lift-off (4).
- With circular finger movements, affix the base material of the electrode to the skin until it adheres entirely (5).
- Start the ECG.
- After use, carefully remove the electrode using one hand whilst supporting the skin beneath it with the other. Ripping off the electrode or peeling it off quickly may damage the skin. Exercise particular caution if the skin is exceptionally sensitive, especially with children but also with elderly patients, diabetics, or patients on a prolonged course of certain medications known to cause dermatological side effects.
- Always store any electrodes that are not immediately required in the original pouch. Close the opened pouch by folding it over twice at the opening.
- Once the pouch has been opened, use the electrodes within 7 days!
- Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices): if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

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