

sigvaris

Live your life. Be yourself.



**Rolly**

Donning and  
doffing device

### Rolly

Compression garments can sometimes be difficult to put on and take off; particularly if you have limited mobility or difficulty using your hands.

The Rolly is the ideal solution to allow your compression garment to move effortlessly up your limb.

#### Multi functional

Suitable for installing garments for below knee, thigh and armsleeve.

#### Easy removal

Unlike many other donning and doffing aids, the Rolly is also superb at removing the compression garment with ease.



### Product materials


- Donner - Latex Free  
Block Copolymer
- Hydrogenated Styrene Isoprene Butadiene Block
  - Antioxidant
  - Anti-blocking agent

Fluid (soapy water)

- Polyacrylamide
- Acrylamide
  - Sodium azide
  - Water

## Loading your Rolly


Follow these simple steps to load your product on to your Rolly.



<b>Step 1</b>	Put the stocking on your arm, make sure there are no wrinkles and that the stocking is stretched up the arm.
<b>Step 2</b>	With the product stretched over your arm, put your hand through your Rolly and roll your Rolly up to the top of your stocking.
<b>Step 3</b>	Pull the stocking over your Rolly by pulling downward. Then roll back your Rolly slightly up your arm in order to tuck the top end of the stocking under the base of your Rolly.
<b>Step 4</b>	Roll your Rolly to the end of your hand, winding down the stocking in the process. Then remove your Rolly (with your stocking attached), from your arm. Once finished winding it down, the toe area of your stocking at one end will be visible.

#### How to don on to your arm?

Follow the **Loading your Rolly** section in this leaflet to load your armsleeve. Then (using your non affected arm) follow these easy steps to don your armsleeve.





<b>Step 1</b>	With the opening for wrist/handpiece pointing towards you, place your hand through your Rolly and start to unroll your Rolly over your hand and then over your wrist.
<b>Step 2</b>	Continue to unroll your Rolly right the way up your arm towards your bicep and the top of your arm.
<b>Step 3</b>	Once at the top of the arm, position the griptop/top of the armsleeve into the correct place.
<b>Step 4</b>	Roll your Rolly downwards towards your wrist, then over your hand where your Rolly is then free of your armsleeve.

#### How to don on to your leg?

Follow the **Loading your Rolly** section in this leaflet to load your garment.


Follow these simple steps to don your leg hosiery.



<b>Step 1</b>	Wind the stocking all the way down the Rolly until all that is visible is the toe area of the stocking facing you.
<b>Step 2</b>	Place your toe against the stocking's toe area. Make sure that the heel of the stocking is under the Rolly.
<b>Step 3</b>	Push your foot into the Rolly, rolling it along the foot and up over the heel. Continue this until the stocking is up your leg and position the griptop/top into the correct place.
<b>Step 4</b>	Remove your Rolly, letting it roll down the leg towards your ankle. Use the ridged surface to smooth the stocking into place.

#### How to doff from your leg?

Your Rolly can also be used to remove your garment with ease and without damaging it. Follow these simple steps to doff your leg hosiery.



<b>Step 1</b>	Put your foot through the Rolly, using your hands, move the Rolly past your heel and all the way up your leg to the top of the stocking.
<b>Step 2</b>	Pull the stocking up over the Rolly and pull the stocking over the Rolly.
<b>Step 3</b>	Tuck the end of the stocking under the base of the Rolly.
<b>Step 4</b>	Roll your Rolly downwards towards your ankle, then over your foot. Your Rolly is then free of <del>your stocking</del>

# Rolly

---

## Instructions for use

- Before starting, remove any jewellery, rings and watches (if possible) to avoid damaging your garment.
- Do not place your Rolly on a radiator or any other direct heat source.
- Hand wash only
- Do not place any heavy or sharp objects on the product.
- Do not pull over your head

## Product care

- To clean your Rolly, use soap and water.
- To disinfect your Rolly, use isopropyl alcohol.
- If your Rolly becomes damaged and water leaks out, stop using it and discard the product. Although the liquid is not toxic, avoid excessive contact with your skin and avoid getting any in your eyes.

Do not hesitate to ask your Health Care professional if you have any concerns.