
ΔVANOS\*


MIC\*  
JEJUNAL  
FEEDING TUBE  
WITH ENFit® CONNECTOR

Placement Instructions for Use

J-Tube

Fig 1

|                    |                                                |                 |                                 |                                  |
|--------------------|------------------------------------------------|-----------------|---------------------------------|----------------------------------|
| Diameter           | Length                                         | Single Use Only | Sterilized by Gamma Irradiation | Do not use if package is damaged |
| Do not resterilize | Product is NOT made with DEHP as a plasticizer | Rx Only         | MR Safe                         | Caution                          |
|                    |                                                |                 |                                 | Consult instructions for use     |

## Instructions for Use

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

## Description

The AVANOS® MIC® Jejunal (J) Feeding Tube (**Fig. 1**) allows for delivery of enteral nutrition and medication directly into the distal duodenum or proximal jejunum.

## Indications for Use

The AVANOS® MIC® J Feeding Tube is indicated for use in patients who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy.

## Contraindications

Contraindications for placement of a J feeding tube include, but are not limited to:

- Colonic introposition
- Ascites
- Portal hypertension
- Peritonitis
- Uncorrected coagulopathy
- Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)
- Lack of adherence of the stomach to the abdominal wall (replacement only)
- Lack of established gastrostomy tract (replacement only)
- Evidence of infection around stoma site (replacement only)
- Presence of multiple stoma fistulous tracts (replacement only)

## ⚠ Warning

**Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.**

## Complications

The following complications may be associated with any J feeding tube:

- Nausea, vomiting, abdominal bloating or diarrhea
- Aspiration
- Peristomal pain
- Abscess, wound infection and skin breakdown
- Pressure necrosis
- Hypergranulation tissue
- Intraperitoneal leakage
- Buried bumper syndrome
- Peristomal leakage
- Balloon failure or tube dislodgement
- Tube clog
- Gastrointestinal bleeding and/or ulceration
- Gastric outlet obstruction
- Ileus or gastroparesis
- Bowel and gastric volvulus
- Hemorrhage
- Perforation
- Intussusception

Other complications such as abdominal organ injury may be associated with the procedure to place the feeding tube.

**Note:** Verify package integrity. Do not use if package is damaged or sterile barrier compromised.

**Note:** The risk for perforation may be higher in patients weighing <10kg.

## Placement

This AVANOS® J Feeding Tube may be placed in the following ways:

- Ⓢ Surgically
- Ⓡ Percutaneously under fluoroscopic (radiologic) guidance
- ⓔ Percutaneously under endoscopic guidance
- ⓧ A replacement to an existing device using an established stoma tract

**⚠ Warning: A gastrostomy must be performed to affix the stomach to the anterior abdominal wall, the feeding tube insertion site identified and stoma tract dilated prior to initial tube insertion to ensure patient safety and comfort.**

**⚠ Warning: Do not use the retention balloon of the feeding tube as a gastrostomy device. The balloon may burst and fail to attach the stomach to the anterior abdominal wall.**

## Ⓢ Ⓡ ⓔ ⓧ Tube Preparation

**⚠ Warning: Verify package integrity. Do not use if package is damaged or sterile barrier compromised.**

1. Select the appropriate size MIC® J Feeding Tube, remove from the package and inspect for damage.
2. Adjust the length of the tube by trimming the distal end with a razor blade or scalpel. Ensure that the cut is smooth and blunt and the length sufficient to be placed 10–15 cm beyond the Ligament of Treitz.

3. Using a male Luer syringe, inflate the balloon (**Fig. 1-D**) with water through the Balloon Inflation Port (**Fig. 1-A**). Do not use air.
  - Inflate the balloon with 2–3 ml of water for low volume tubes identified by LV following the REF code number.
  - Inflate the balloon with 7–10 ml of water for standard tubes.
4. Remove the syringe and verify balloon integrity by gently squeezing the balloon to check for leaks. Visually inspect the balloon to verify symmetry. Symmetry may be achieved by gently rolling the balloon between the fingers. Reinsert the syringe and remove all the water from the balloon.
5. Check the external retention bolster (**Fig. 1-C**). The bolster should slide along the tube with moderate resistance.
6. Inspect the entire length of the tube for any irregularities.
7. Using an ENFit® syringe, flush the jejunal lumen (**Fig. 1-B**) of the tube with water to confirm tube patency.
8. Lubricate the tip of the tube with a water soluble lubricant. Do not use mineral oil. Do not use petroleum jelly.
9. Generously lubricate the jejunal lumen with water-soluble lubricant. Do not use mineral oil. Do not use petroleum jelly.

## Ⓢ Suggested Surgical Procedure (Stamm Gastrostomy)

1. Through a midline laparotomy, identify the pylorus and the superior epigastric artery in the abdominal wall.
2. The gastrostomy site should be 10–15 cm from the pylorus on the greater curvature of the stomach. The gastrostomy site should also be at least 3 cm away from the costal margin to prevent damage to the retention balloon by abrasions during movement.
3. Place two concentric purse string sutures around the site. Leave the purse string needles in place.
4. On the anterior parietal peritoneum, select an exit site that approximates the gastrostomy. Avoid the superior epigastric artery, drains, or other stomas.
5. Make a stab wound from the anterior parietal peritoneum to the extra abdominal surface. Insert the tube from the outside to the inside the abdominal cavity.

**Note:** The use of a right angle clamp may be used to facilitate placement.

**⚠ Caution:** Never use a clamp with sharp teeth or tenaculum to pull the tube into position. This will damage the tube.

6. Using two Babcock clamps on the anterior stomach surface, “tent” the stomach.
7. Use electrocautery or a scalpel to open the stomach.
8. Dilate the enterotomy with a hemostat.

## Ⓡ ⓔ Suggested Site Preparation

1. Use standard Radiologic or Endoscopic techniques to visualize and prepare for gastrostomy tube placement.
2. Verify that no abnormalities are present that could pose a contraindication to placement of the tube and place the patient in the supine position.
3. Select a gastrostomy site that is free of major vessels, viscera, and scar tissue. The site is usually one third the distance from the umbilicus to the left costal margin at the midclavicular line.
4. Prep and drape the selected insertion site per facility protocol.

## Ⓡ ⓔ Gastrostomy Placement

**⚠ Warning: It is recommended to perform a three point gastrostomy in a triangle configuration to ensure attachment of the gastric wall to the anterior abdominal wall.**

1. Place a skin mark at the tube insertion site. Define the gastrostomy pattern by placing three skin marks equidistant from the tube insertion site and in a triangle configuration.
 

**⚠ Caution:** Allow adequate distance between these insertion site and gastrostomy placement to prevent interference of the T-Fastener and inflated balloon.
2. Localize the puncture sites with 1% lidocaine and administer local anesthesia to the skin and peritoneum.
3. Place the first T-Fastener and confirm intragastric position. Repeat the procedure until all three T-Fasteners are inserted at the corners of the triangle.
4. Secure the stomach to the anterior abdominal wall and complete the procedure.

## Ⓡ Create the Stoma Tract under Fluoroscopic (Radiologic) Visualization

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastrostomy pattern. With fluoroscopic guidance, confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.
 

**⚠ Warning: Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

**⚠ Warning: Take care not to advance the puncture needle too deeply in order to avoid puncturing the posterior gastric wall, pancreas, left kidney, aorta or spleen.**
2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038” compatible introducer needle at the center of the gastrostomy pattern into the gastric lumen directed toward the pylorus.
 

**Note:** For J tube placement, the best angle of insertion is a 45 degree angle to the surface of the skin.
4. Use fluoroscopic visualization to verify correct needle placement. Additionally, to aid in verification, a water filled syringe may be attached to the needle hub and air aspirated from the gastric lumen.

**Note:** Contrast may be injected upon return of air to visualize gastric folds and confirm position.

- Advance a guidewire, up to .038", through the needle and coil in the fundus of the stomach. Confirm position.
- Remove the introducer needle, leaving the guidewire in place, and dispose of according to facility protocol.
- Advance a .038" compatible flexible catheter over the guidewire and using fluoroscopic guidance, manipulate the guidewire into the antrum of the stomach.
- Advance the guidewire and flexible catheter until the catheter tip is at the pylorus.
- Negotiate through the pylorus and advance the guidewire and catheter into the duodenum and 10–15 cm beyond the Ligament of Treitz.
- Remove the catheter, leaving the guidewire in place, and dispose of according to facility protocol.

## Ⓔ Create the Stoma Tract under Endoscopic Visualization

- Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastropexy pattern. With endoscopic guidance, confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.  
⚠️ **Warning:** Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.  
⚠️ **Warning:** Take care not to advance the puncture needle too deeply in order to avoid puncturing the posterior gastric wall, pancreas, left kidney, aorta or spleen.
- Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
- Insert a .038" compatible introducer needle at the center of the gastropexy pattern into the gastric lumen directed toward the pylorus.  
**Note:** For J tube placement, the best angle of insertion is a 45 degree angle to the surface of the skin.
- Use endoscopic visualization to verify correct needle placement.
- Advance a guidewire, up to .038", through the needle and into the stomach. Confirm position.
- Use endoscopic visualization to grasp the guidewire with atraumatic forceps.
- Remove the introducer needle, leaving the guidewire in place, and dispose of according to facility protocol.

## Ⓔ Ⓔ Dilation

- Use a #11 scalpel blade to create a skin incision that extends alongside the guidewire, downward through the subcutaneous tissue and fascia of the abdominal musculature. After the incision is made, dispose of the scalpel according to facility protocol.
- Advance a dilator over the guidewire and dilate the stoma tract at least four French sizes larger than the enteral feeding tube being placed.
- Remove the dilator over the guidewire, leaving the guidewire in place.  
**Note:** After dilation, a peel-away sheath may be used to facilitate advancement of the tube through the stoma tract.

## Ⓔ Surgical Tube Placement

- Advance the MIC\* J Feeding Tube until the balloon is in the stomach.
- Palpate the tube through the duodenum. When satisfied with the placement, check the position. The tip should lie 10–15 cm past the Ligament of Treitz.
- Using a male Luer syringe, inflate the balloon.
  - Inflate the LV balloon with 2–3 ml of sterile or distilled water.
  - Inflate the standard balloon with 7–10 ml of sterile or distilled water.
 ⚠️ **Caution:** Do not exceed 5 ml total balloon volume inside the LV balloon and 20 ml total balloon volume in the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Tie the purse string sutures around the tube.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Use the purse string sutures to attach the stomach to the peritoneum. Take care to avoid puncture of the balloon.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 2–3 mm above the skin. Do not suture the bolster to the skin.

## Ⓔ Tube Placement under Fluoroscopic (Radiologic) Visualization

- Advance the distal end of the tube over the guidewire, through the stoma tract and into the stomach.
- Rotate the MIC\* J Feeding Tube while advancing to facilitate passage of the tube through the pylorus and in the jejunum.
- Advance the tube until the tip of the tube is 10–15 cm beyond the Ligament of Treitz and the balloon is in the stomach.
- Using a male Luer syringe, inflate the balloon.
  - Inflate the LV balloon with 2–3 ml of sterile or distilled water.
  - Inflate the standard balloon with 7–10 ml of sterile or distilled water.
 ⚠️ **Caution:** Do not exceed 5 ml total balloon volume inside the LV balloon and 20 ml total balloon volume inside the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 2–3 mm above the skin. Do not suture the bolster to the skin.
- Remove the guidewire.

## Ⓔ Tube Placement under Endoscopic Visualization

- Advance the distal end of the tube over the guidewire, through the stoma tract and into the stomach.
- Grasp the suture loop or the tip of the tube with atraumatic forceps.

- Advance the MIC\* J Feeding Tube through the pylorus and upper duodenum. Continue to advance the tube using the forceps until the tip is positioned 10–15 cm beyond the Ligament of Treitz and the balloon is in the stomach.
- Release the tube and withdraw the endoscope and forceps in tandem, leaving the tube in place.
- Using a male Luer syringe, inflate the balloon.
  - Inflate the LV balloon with 2–3 ml of sterile or distilled water.
  - Inflate the standard balloon with 7–10 ml of sterile or distilled water.
 ⚠️ **Caution:** Do not exceed 5 ml total balloon volume inside the LV balloon and 20 ml total balloon volume inside the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 2–3 mm above the skin. Do not suture the bolster to the skin.
- Remove the guidewire.

## Ⓔ Ⓔ Ⓔ Verify Tube Position and Patency

- Verify proper tube placement radiologically to avoid potential complication (e.g. bowel irritation or perforation) and ensure the tube is not looped within the stomach or small bowel.
- Flush the jejunal lumen with water to verify patency.
- Check for moisture around the stoma. If there are signs of gastric leakage, check the tube position and placement of the external retention bolster. Add sterile or distilled water as needed in 1–2 ml increments. Do not exceed balloon capacity as indicated previously.
- Check to ensure that the external retention bolster is not placed too tightly against the skin and rests 2–3 mm above the abdomen for initial placement and 1–2 mm above the abdomen for a replacement tube.
- Begin feeding only after confirmation of proper patency, placement and according to physician instructions.

## Ⓔ Replacement Procedure Through an Established Gastrostomy Tract

- Cleanse the skin around the stoma site and allow the area to air dry.
- Select the appropriate size MIC\* J Feeding Tube and prepare according to the instructions in the Tube Preparation section above.
- If using endoscopic visualization, perform routine esophagogastroduodenoscopy (EGD). Once the procedure is complete and no abnormalities are identified that could pose a contraindication to placement of the tube, place the patient in the supine position and insufflate the stomach with air. Manipulate the endoscope until the indwelling gastrostomy tube is in the visual field.
- Under either fluoroscopic guidance or endoscopic guidance, insert a floppy-tipped guidewire, up to .038", through the indwelling gastrostomy tube.
- Remove the gastrostomy tube over the guidewire.
- Place the MIC\* J Feeding Tube according to the Tube Placement section above.  
⚠️ **Caution:** For a replacement feeding tube, gently slide the external retention bolster to approximately 1–2 mm above the skin. Do not suture the bolster to the skin.
- Verify tube position and patency using the Verify Tube Position and Patency section above.

## Tube Patency Guidelines

Proper tube flushing is the best way to avoid clogging and maintain tube patency. The following are guidelines to avoid clogging and maintain tube patency.

- Flush the feeding tube with water every 4–6 hours during continuous feeding, anytime the feeding is interrupted, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube before and after medication administration and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and to mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric-coated medication or mix medication with formula.
- Avoid using acidic fluids such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.

## General Flushing Guidelines

Flush the feeding tube with water using an ENFit® syringe every 4–6 hours during continuous feeding, anytime the feeding is interrupted, at least every 8 hours if the tube is not being used, or per clinician's instructions. Flush the feeding tube before, after, and in between medication administration. Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes.

- Use a 30 to 60 ml ENFit® syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
- Use room temperature water for tube flushing. Sterile water may be appropriate where the quality of municipal water supplies is of concern. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

- Document the time and amount of water used in the patient's record. This will enable all caregivers to monitor the patient's needs more accurately.

## Nutrition Administration

1. Open the cap to the access port of the MIC\* J Feeding Tube.
2. Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
3. Remove the flushing syringe from the access port.
4. Securely connect an ENFit® feed set to the access port.  
⚠️ **Caution:** Do not over-tighten the feed set connector or the syringe to the access port.
5. Complete feeding per the clinician's instructions.
6. Remove the feed set or syringe from the access port.
7. Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
8. Remove the flushing syringe from the access port.
9. Close the cap to the access port.

## Medication Administration

Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric coated medication or mix medication with formula.

1. Open the cap to the access port of the MIC\* J tube.
2. Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
3. Remove the flushing syringe from the access port.
4. Securely connect an ENFit® syringe containing the medication to the access port.  
⚠️ **Caution:** Do not over-tighten the syringe to the access port.
5. Deliver the medication by depressing the ENFit® syringe plunger.
6. Remove the syringe from the access port.
7. Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
8. Remove the flushing syringe from the access port.
9. Close the cap to the access port.

## Balloon Maintenance

Precise balloon life cannot be predicted. Silicone balloons generally last 1–8 months, but the life span of the balloon varies according to several factors. These factors may include medications, volume of water used to inflate the balloon, gastric pH, and tube care.

Check the water volume in the balloon once a week.

- Insert a male Luer syringe into the Balloon Inflation Port (BAL.) and withdraw the fluid while holding the tube in place. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. If the amount is less than recommended or prescribed, refill the balloon with the water initially removed, then draw up and add the amount needed to bring the balloon volume up to the recommended and prescribed amount of water. Be aware as you deflate the balloon there may be some gastric contents that can leak from around the tube. Document the fluid volume, the amount of volume to be replaced (if any), the date and time.
- Wait 10–20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. A deflated or ruptured balloon could cause the tube to dislodge or be displaced. If the balloon is ruptured, it will need to be replaced. Secure the tube into position using tape, then follow facility protocol and/or call the physician for instructions.  
⚠️ **Caution:** Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

## Daily Care & Maintenance Check List

- Assess the patient for any signs of pain, pressure or discomfort.
- Assess the stoma site for any signs of infection, such as redness, irritation, edema, swelling, tenderness, warmth, rashes, purulent, or gastrointestinal drainage. Assess the patient for any signs of pressure necrosis, skin breakdown, or hypergranulation tissue.
- Use warm water and mild soap.
- Use a circular motion moving from the tube outwards.
- Rinse thoroughly and dry well.
- Assess the tube for any abnormalities such as damage, clogging, or abnormal discoloration.
- Use warm water and mild soap being careful not to pull or manipulate the tube excessively.
- Rinse thoroughly and dry well.
- Clean the Jejunal and Balloon Inflation Ports. Use a cotton tip applicator or soft cloth to remove all residual formula and medication.
- Verify that the external bolster rests 1–2 mm above the skin.
- Flush the feeding tube as described in the General Flushing Guidelines section above.  
⚠️ **Caution:** Do not rotate the external retention bolster. Rotating the bolster may cause the tube to kink and possibly lose position.

## Tube Occlusion

Tube occlusion is generally caused by:

- Poor flushing techniques
- Inappropriate administration of medication
- Pill fragments
- Thick formulas, such as concentrated or enriched formulas that are generally thicker
- Formula contamination that leads to coagulation
- Reflux of gastric or intestinal contents up the tube

## To Unclog a Tube

1. Make sure that the feeding tube is not kinked or clamped off.
2. If the clog is visible above the skin surface, gently massage or milk the tube between fingers to break up the clog.
3. Connect a 30 to 60 ml ENFit® syringe filled with warm water into the Jejunal access port of the tube and gently pull back on then depress the plunger to dislodge the clog. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
4. If the clog remains, repeat step #3. Gentle suction alternating with syringe pressure will relieve most obstructions.
5. If this fails, consult with the physician. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients. If the clog is stubborn and cannot be removed, the tube will have to be replaced.  
⚠️ **Caution:** Do not insert foreign objects through the tube.

## MRI Safety Information

The MIC\* Jejunal Feeding Tube is MR Safe.

⚠️ **Warning:** For enteral nutrition and/or enteral medication only.

For more information, please call 1-844-4AVANOS (1-844-428-2667) in the United States, or visit our web site at [www.avanos.com](http://www.avanos.com).

Educational Booklets: "A Guide to Proper Care" and "A Stoma Site and Enteral Feeding Tube Troubleshooting Guide" are available upon request. Please contact your local representative or contact Customer Care.

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[illegible]

[illegible]

Distributed in the USA by Avanos Medical Sales, LLC,  
5405 Windward Parkway, Alpharetta, GA 30004 USA.  
In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com  
 Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.  
 Avanos Medical Belgium BV,  
Leonardo da Vincilaan 1, 1930 Zaventem, Belgium.  
Sponsored in Australia by Avanos Medical Australia Pty Ltd,  
475 Victoria Avenue, Chatswood, NSW 2067, Australia.

\*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.  
© 2018-2024 AVNS. All rights reserved. 2024-10-24  
15-M1-277-02 / 50002130

CE 2797