

## Explanation of Symbols



Do not use if package is damaged



## COVADERM WOUND DRESSINGS

## COVADERM® WOUND DRESSINGS

STERILE EO

Covaderm® wound dressings are sterile, absorbent, multi-layer composite wound dressings. They are available with and without a polyurethane barrier film layer. All Covaderm, Covaderm Thin, Covaderm Plus® and Covaderm Plus Thin have rounded edges, non-adherent contact layers, and easy to shape adhesive coated fabric tape.

## USE

Covaderm wound dressings absorb wound drainage, promote a moist environment helpful to healing, and serve as a barrier against water and bacteria when used as a secondary dressing.

## INDICATIONS

- Covaderm wound dressings are indicated for minor abrasions, minor lacerations, and minor scalds and burns.
- When under the care or direction of a licensed practitioner they are indicated for:
  - Dermal ulcers (e.g., pressure ulcers, leg ulcers, and other secreting lesions)
  - Second degree burns
  - Donor sites
  - More severe lacerations, cuts and abrasions
  - Other superficial wounds
  - Covaderm Plus VAD is also indicated for IV sites, and other tubular sites and incisions.

## CONTRAINDICATIONS

- Covaderm wound dressings are contraindicated for:
  - Third-degree burns
  - Heavily bleeding wounds
  - Individuals who are sensitive to or have had an allergic reaction to the dressing or its materials

## WARNINGS

- Read and understand all warnings, cautions, and directions completely and carefully before use. Failure to do so may result in serious injury.
- Do not reuse, reprocess, or re-sterilize this dressing. This device is for single use only. Reuse, reprocessing, or re-sterilization may damage the dressing, and may result in injury, infection, or death.
- Do not use if the individual package is opened or if the package or the dressing is damaged. Carefully inspect the dressing for damage or defects prior to use.
- Seek medical attention if your wound is more severe than a minor abrasion, minor laceration, or minor scald and burn. Also seek medical attention if you have signs of an infection or allergic reaction such as fever, redness, pain, or swelling.
- If you have any signs of an allergic reaction, discontinue using the dressing and rinse the wound and surrounding skin.
- If applicable, ensure the dressing adheres to skin. Dressings with weak adhesives should be replaced immediately.
- If dressing is ineffective or does not absorb, discard and replace.
- If applicable, dispose of all used dressings with care to avoid exposure to or contamination of biohazard.

- Ensure the dressing fully covers your wound.
- If the dressing rips or tears prior to or during application, discard the torn dressing and replace it.
- Prior to removal, ensure the dressing has not adhered to the wound. If the dressing has adhered to the wound, rinse with water or saline solution to loosen.
- Take proper care when removing a dressing with adhesive as it may result in skin irritation. If the adhesive is too strong or the skin is fragile, rinse with water or saline solution to loosen.

## DIRECTIONS FOR USE

1. Carefully read the **WARNINGS** and **CAUTIONS** prior to use.
2. Inspect the packaging for damage or defects. Replace the entire dressing if the packaging is damaged.
3. Clean the wound of excess drainage and fluid, and thoroughly dry skin around the wound.
4. Open the dressing package.

## ADHESIVE DRESSINGS

1. Peel the paper liner from the back of the adhesive dressing.
2. Center the contact layer (gauze) over the wound.
3. Change dressing as needed when wet, soiled, leaking, or as directed by a licensed practitioner.
4. To remove the dressing, hold the edge of the adhesive border and carefully pull it in the direction of hair growth, lifting bandage from the wound.
5. Dispose of the soiled dressing.

## STORAGE AND TRANSPORT CONDITIONS

Store at room temperature and avoid excessive heat and humidity.

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

## WARRANTY

DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**