NEOMED®

ENTERAL FEEDING TUBE

WITH ENFit® CONNECTOR

- FR

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EN

Instructions for Use

Rx Only Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.
Product is not made with natural rubber latex.
Product is not made with DEHP as a plasticizer.
Product is not made with BPA.
This instructions for use applies to the following NEOMED® Enteral Feeding Tubes with ENFit® Connector.

Material	French (Fr) Size	Length (cm)
Polyurethane (PUR)	4,0	40
	5,0, 6,5, 8,0	40, 60, 90
	5,0, 8,0	105
	10,0, 12,0	109
Polyvinyl Chloride (PVC)	5,0, 6,5, 8,0	40, 60
	6,5, 8,0	90
Silicone (SIL)	5,0, 6,5, 8,0	40, 60, 90

Device Description / Intended Purpose
The NEOMED® Enteral Feeding Tube is intended for use in patients who require intermittent or continuous tube feedings via the nasogastric or orogastric pathway. Material options include Polyurethane, Silicone, or PVC. Refer to product packaging for product specific information.

 Outer diameter		 Length	 Medical device
 Single Use Only		 Sterilized using ethylene oxide	 Single sterile barrier system
 Single sterile barrier system with protective packaging outside	 Do not use if package is damaged	 Do not resterilize	Rx Only - Prescription only
 MR Safe	 Caution	 Consult instructions for use	 Keep dry
 Keep away from sunlight	 Country of manufacture	 Manufacturer	 Catalogue number
 Unique Device Identifier	 Batch code	 Date of Manufacture	 Use-by date

Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway, Alpharetta, GA 30004 USA. In USA, 1-844-4AVANOS (1-844-428-2667) . www.avanos.com

 Avanos Medical, Inc.,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.

Sponsored in Australia by Avanos Medical
Australia Pty Ltd, 475 Victoria Avenue,
Chatswood, NSW 2067, Australia.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2024-02-0115-M1-1981-01

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This product is provided sterile (and non-sterile). Sterility on the instructions for use is not valid unless there is a sterilization method indicated on the product label.

Intended Users

This device is intended to be used by healthcare professionals.

Patient Target Group

This device is intended to be used in patients of all age groups who require nasogastric/oralgastric enteral feeding.

Indications for Use

This product is intended for use in neonatal and pediatric patients to provide nutrition via nasal and oral gastric placement and is not intended for use beyond 30 days.

Contraindications

This device is intended to connect to products with ENFit® compatible connectors.

This device is not intended for use longer than 30 days.

This device is not intended for use in patients with congenital anomalies of the GI tract above the stomach requiring surgical intervention.

Residual Risks and Side Effects

While side effects vary depending on the type of procedure performed, mitigations have been implemented to minimize the following residual risks related to the use of the NEOMED®

Nasogastric Feeding Tubes and the associated NEOMED® Extension Sets: abnormal reproductive development, allergic/systemic toxic reaction, over/under delivery of medication (potentially life-threatening)/nutrition, skin/tissue damage/irritation (includes trachea, vocal cords, gastrointestinal tract/ulcer/necrosis), infection/fever/sepsis, nausea/vomiting/bloating/diarrhea, delayed or prolonged treatment/procedure/therapy, pneumothorax, choking, long-term carcinogenic and/or mutagenic effects, aspiration, peritonitis, gastroenteritis, laceration, pain.

⚠ Warnings / Cautions

Do not reuse devices labeled for single use only. Reuse of a single use device that has come in contact with blood, bone, tissue, or other body fluids may lead to patient or user injury. Possible risks associated with the reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.

Failure to read entire instructions for use may result in severe patient harm or death.

Use of the device by personnel not well versed in anatomical landmarks, safe technique and potential complications may lead to severe patient harm or death.

This device is intended for one-time use. Each device should be used on only one patient.

Do not autoclave or re-use device as this may increase the risk of leakage or biological infection.

Surgical implantation, inadvertent placement into the trachea or inadequate insertion into the stomach may result in severe patient harm or death.

Do not use a stylet or guidewire as this may lead to severe patient harm or death.

Use of the device in non-enteral applications may lead to severe patient harm or death.

Do not use if primary pouch has been previously opened or damaged as this may lead to device contamination or biological infection.

Do not use visibly damaged device or device with illegible markings as this may lead to inadequate nutrition delivery or patient harm.

Do not use forceps, clamps or sharp instruments as this may damage device and lead to inadequate nutrition delivery or severe patient harm.

Do not modify device as this may lead to inadequate nutrition delivery or severe patient harm.

Do not use if package has been previously opened or damaged.

Follow hospital protocols for feeding tube placement and maintenance or current ADA recommendation, as applicable. Use aseptic technique.

MR Information

** MR Safe**

The NEOMED® Enteral Feeding Tube contains only nonconducting, nonmetallic and non magnetic items and poses no known hazards in all MR imaging environments.

Clinical Benefits

Nasogastric/ Orogastric Feeding Tubes shall facilitate nutrition (and fluid and medication) delivery to stomach via the nose or mouth.

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Device Lifetime

Device has been validated for use for up to 30 days.

Recommended Placement Procedure

Insertion Measurement: measure from the tip of the nose to the earlobe and then from the earlobe to midway between the xyphoid process and umbilicus. Centimeter markings are provided on the catheter to assist in the feeding tube placement.

A strip of sterile tape can be placed at the target insertion depth of the feeding tube for reference.

Insert the tube slowly and cautiously into the mouth or nose (lubricant is optional); allow the infant to swallow as the tube is advanced.

Continue to advance to the predetermined length.

⚠.Precautions: avoid inadvertent placement of the tube into the trachea or inadequate insertion such that the proximal side hole is not fully placed in the stomach.

Secure the feeding tube either by holding the tube or taping to the patient before checking for placement.

Verify Feeding Tube Placement

Finalize securement, and verify and document according to hospital protocol (e.g., chest radiograph, aspiration of gastric contents / pH testing). If the tip placement remains uncertain reposition or remove.

Recommended Maintenance

The feeding tube should be visually checked for correct position before each administration utilizing tube markings.

Verify proper tube placement before each administration according to hospital protocol.

Flush feeding tube with sterile water before and after each administration.

Plug feeding tube with tethered closure when not in use.

Recommended Feeding Tube Removal

Flush the feeding tube with 0.5cc–1cc of sterile water.
Plug with tethered closure and/or clamp the feeding tube during tube removal to minimize the risk of aspirating fluids into the trachea during withdrawal.
Gently and slowly withdraw the feeding tube and discard per standard hospital protocol and/or local or national regulations.
Document the feeding tube removal.

Complaint and Adverse Event Reporting

If a serious incident has occurred in relation to the device, it should be reported to the manufacturer and the competent authority in which the user and/or patient is established. A serious incident means any incident that directly or indirectly led, might have led or, in case of recurrence, could lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's, fetus or other person's state of health, or a serious public health threat.

Disposal

Contact the local authorities to determine the proper method of disposal of potentially biohazardous parts and accessories. Disposal of the device and accessories should be in accordance with local, state, and national laws and regulations.

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