

MIC*
GASTRO-ENTERIC
TUBE

SURGICAL AND ENDOSCOPIC PLACEMENT

Instructions for Use

 Diameter	 Length	 Single Use Only		STERILE R Sterilized by Gamma Irradiation	 Do not use if package is damaged
 Do not resterilize	Product is NOT made with DEHP (Di(2-ethylhexyl) phthalate) as a plasticizer	Rx Only	 MR Conditional	 Caution	 Consult instructions for use

Instructions for Use

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician. The indications for use of the MIC* Gastro-Enteric Tube for enteral feeding and gastric decompression are wide and varied. The MIC* Gastro-Enteric Tube can be used for gastric motility complications, gastroparesis, head injury or stroke complications, or other situations where gastric decompression and enteral feeding distal to the stomach are indicated.

Contraindications

Contraindications for placement of a GE feeding tube include, but are not limited to:

- Colonic interposition
- Ascites
- Portal hypertension
- Peritonitis
- Uncorrected coagulopathy
- Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)

⚠ Warning

Patients receiving elemental diets through the jejunal lumen of the MIC* Gastro-Enteric Tube must be placed on a feeding pump to assure precise delivery and tube patency.

Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

Complications

The following complications may be associated with any GE feeding tube:

- Nausea, vomiting, abdominal bloating or diarrhea
- Aspiration
- Peristomal pain
- Abscess, wound infection and skin breakdown
- Pressure Necrosis
- Hypergranulation Tissue
- Intraperitoneal Leakage
- Buried bumper syndrome
- Peristomal leakage
- Balloon failure or tube dislodgement
- Tube clog
- Gastrointestinal bleeding and/or ulceration
- Gastric outlet obstruction
- Ileus or gastroparesis
- Bowel and gastric volvulus
- Perforation
- Intussusception
- Hemorrhage

Other complications such as abdominal organ injury may be associated with the procedure to place the feeding tube.

Note: Verify package integrity. Do not use if package is damaged or sterile barrier compromised.

Note: The risk for perforation may be higher in patients weighing <10kg.

Key Terms

A Jejunal Feeding Port	F Gastric Exit Ports
B Gastric Decompression Port	G Radiopaque Jejunal Lumen
C Balloon Inflation Port	H Jejunal Exit Ports
D SECUR-LOK* External Retention Ring	I Weighted Tip
E Internal Retention Balloon	

Surgical Placement

⚠ Warning

Do not use mineral oil or petroleum jelly.

Never fill the balloon with air.

The Stamm procedure is a common surgical technique for the initial insertion of the MIC* Gastro-Enteric Tube. The exit site of the gastrostomy portion should be placed on the greater curvature of the stomach.

Check the MIC* Gastro-Enteric Tube prior to surgical placement. Inflate the balloon with water or saline. The balloon should not leak and should be symmetrical. Silicone material has a tendency to adhere to itself causing the balloon to inflate asymmetrically. If this occurs, roll and bend the balloon end of the catheter in your fingers until the balloon becomes symmetrical. Also check the ventilated SECUR-LOK* Ring on the tube. The ring should slide with moderate resistance. Check the entire length of the tube. It should be clean and free of defects. Check patency on both the gastric and jejunal lumens by injecting saline into the catheter prior to insertion.

Note: Before placement, coat the guidewire generously with a water soluble lubricant. Practice inserting and removing the guidewire.

The jejunal portion of the MIC* Gastro-Enteric Tube is 56.50 cm long (25 cm for LV configuration) as measured from distal aspect of the retention balloon. Both are provided with a weighted tip.

Note: Remove the guidewire prior to making any length adjustments. Make all adjustments prior to placing the tube in the patient.

During the Stamm procedure, the tube is brought through the abdominal wall from the outside in, prior to insertion into the stomach. The use of a right angle clamp works well.

⚠ **Caution:** Never use a clamp with sharp teeth or tenaculum to pull the tube into position. This will damage the tube.

The jejunal portion of the tube can be placed prior to positioning the gastric portion of the tube in the stomach. This is accomplished by inserting the jejunal portion through the stab wound, prior to guiding it through the pylorus. The surgical assistant feeds the jejunal portion of the tube through the stab wound until the gastric portion of the tube is against the abdominal wall. A right angle clamp is used to pull the gastric tube through the abdominal wall, before positioning it in the stomach. Advance the gastric portion through the stab wound into the stomach and inflate the balloon with water or saline.

Standard Volume Balloon Inflation

Inflate the balloon with approximately 7–10 ml of water or saline using a standard luer syringe.

Do not use air. Do not exceed 20 ml.

Low Volume Balloon Inflation

Inflate the balloon with approximately 3–5 ml of water or saline using a standard luer syringe.

Do not use air. Do not exceed 5 ml.

Grasp the body of the tube and withdraw it until slight tension is felt from the balloon coming up against the stomach wall. Gently slide the ventilated SECUR-LOK* ring into position. Allow a space of 1–2 mm between the stoma site and the ring.

⚠ **Caution:** Do not apply excessive tension.

⚠ **Warning:** Do not suture the ventilated SECUR-LOK* ring in position.

A dressing may be used under the ring if desired. After the balloon is inflated and the ventilated SECUR-LOK* ring is in position, the remaining length of the jejunal tube can be positioned inside the bowel. Slowly remove the guidewire being careful not to kink or twist the jejunal portion. Before closing, verify tube patency and placement. Be sure to thoroughly irrigate the jejunal lumen of the tube if a contrast medium is used.

Removal Instructions: To remove the MIC* Gastro-Enteric Tube, deflate the balloon using a luer tipped syringe. Gently pull on the tube until the tip of the gastric portion is free from the stoma site. Grasp the jejunal portion of the tube and gently pull to remove.

Note: If the MIC* Gastro-Enteric Tube is placed surgically, the tube should not be removed until a well established tract has been formed into the stomach.

Replacement Technique

Endoscopic Placement

⚠ **Warning:** Do not use mineral oil or petroleum jelly.

⚠ **Warning:** Never fill the balloon with air.

The endoscopic placement of the MIC* Gastro-Enteric Tube should only be performed on patients with a well formed tract from an existing gastrostomy.

Check the MIC* Gastro-Enteric Tube prior to endoscopic placement. Inflate the balloon with water or saline using a standard luer tipped syringe. The balloon should not leak and should be symmetrical. Silicone material has a tendency to adhere to itself causing the balloon to inflate asymmetrically. If this occurs, roll and bend the balloon end of the catheter in your fingers until the balloon becomes symmetrical. Lubricate the distal end of the gastric portion of the MIC* Gastro-Enteric Tube with a water based lubricant. Also check the ventilated SECUR-LOK* ring on the tube. The ring should slide with moderate resistance. Check the entire length of the tube. It should be clean and free of defects. Check patency of both the gastric and jejunal lumens by injecting saline into the catheter prior to insertion. The guidewire provided with the tube may be used during endoscopic placement. Practice inserting and removing the guidewire to gain familiarity with its resistance characteristics.

The guidewire should be generously lubricated with a water based lubricant.

The standard MIC* Gastro-Enteric Tube jejunal lumen is 56.50 cm long (25 cm for LV configuration) from the distal aspect of the gastric balloon and is provided with a weighted tip. To adjust the length make a blunt end cut, using #10 scalpel.

Note: Remove the guidewire prior to making any length adjustments. Make all adjustments prior to placing the tube in the patient.

To place the MIC* Gastro-Enteric Tube, perform a routine flexible gastroscopy. Insufflate the stomach to gain adequate vision. Back-up the flexible endoscope until the present gastrostomy tube is in the visual field. Remove the existing gastrostomy tube and place the jejunal portion of the MIC* Gastro-Enteric Tube through the gastrostomy site into the stomach. Grasp the distal tip of the jejunal portion with a foreign body or other style atraumatic forceps.

⚠ **Caution:** Do not use a sharp biopsy forceps. Severe damage will occur.

Feed the tube through the pylorus into the upper 1/3 of the duodenum. Release the tube and pull the endoscope and forceps back as a unit.

Continue to feed the jejunal portion of the tube into the duodenum while the assistant pushes the tubing into the gastrostomy site. Push the lubricated gastric tube into the stomach and inflate the balloon with water or saline using a standard luer syringe.

Standard Volume Balloon Inflation: Inflate the balloon with approximately 7–10 ml of water or saline using a standard luer syringe. **Do not use air. Do not exceed 20 ml.**

Low Volume Balloon Inflation: Inflate the balloon with approximately 3–5 ml of water or saline using a standard luer syringe. **Do not use air. Do not exceed 5 ml.** Grasp the body of the tube and withdraw it until slight tension is felt from the balloon coming up against the stomach wall. Gently slide the ventilated SECUR-LOK* ring into position. Allow a space of 1–2 mm between the stoma site and the ring.

⚠ Caution: Do not apply excessive tension.

⚠ Warning: Do not suture the ventilated SECUR-LOK® ring in position.

Remove the guidewire slowly while the physician continues to keep the jejunal portion of the tube in the visual field. If the jejunal portion of the tube should buckle, twist or kink, reinsert the guidewire and reposition the tube using the forceps.

After the guidewire is removed, tube patency and placement must be verified. If a contrast medium is used, thoroughly irrigate the jejunal lumen of the tube to prevent clogging.

Quick Tip: The use of silk sutures tied around the jejunal portion can assist in placing the tube through the pylorus into the jejunum. Place the guidewire through the jejunal lumen of the tube.

Then using 3-0 or 4-0 silk, tie the sutures around the jejunal tube starting at the distal end working proximally. Make the knots about 10 to 15 cm from each other leaving a 1.50 to 2 cm "tag" on each knot to assist the grasping with the forceps.

Note: Do not make the suture too tight, this may cause the tube to kink or occlude.

After preparing the tube, be sure the guidewire moves freely within the jejunal lumen to assure tube patency.

Maintenance of the MIC® Gastro-Enteric Tube

Routine maintenance of the MIC® Gastro-Enteric tube is essential in keeping the jejunal and gastric lumens clear and free of obstructions. It is important to keep the stoma site clean and dry. Keep the ventilated SECUR-LOK® ring in the proper position to prevent tube migration.

Stoma Care: Daily maintenance should include inspection of the stoma site at the time of feeding. Lift the edges of the ventilated SECUR-LOK® ring and gently clean the skin around the stoma with cotton-tip swabs or 4 x 4s. For the healthy stoma, soap and warm water is adequate.

⚠ Warning: If the stoma appears red or irritated, consult the physician at once.

Tube Care: Clean and dry the outside of the tube and both sides of the ventilated SECUR-LOK® ring using 4 x 4s or cotton-tip swabs. The ring can be cleaned in place or moved up the tube for cleaning.

Irrigation: After every use, irrigate the tube with 20–30 ml of water for adults and 10–15 ml for children and infants. The tube should be clear of feeding solutions. Wipe out the feeding port using a swab or a soft cloth to remove any build-up of oils or feeding materials. If the patient is on continuous feeding, the tube should be irrigated 3–5 times daily. Between feedings, introduce a "digesting solution" into the jejunal lumen. Prepare by mixing 20 ml of warm water with 1/2 teaspoon of meat tenderizer. The digesting solution should stand in the tube for several hours followed by irrigation with 20 ml of warm water.

⚠ Caution: Avoid puncturing any part of the MIC® Gastro-Enteric Tube.

Pump feedings only

⚠ Warning: To avoid clogging, feeding pumps must be used to deliver enteral formulas through the jejunal lumen.

Medication delivery: Commercially prepared liquid medications can be administered through the MIC® Gastro-Enteric Tube. Highly viscous solutions should be diluted to a free flowing state before administration.

⚠ Caution: Do not use crushed medications. This could clog the tube.

It is important that the tube be flushed before and after each medication administration.

⚠ Warning: Do not inject medications through balloon inflation valve.

Gastric Decompression

The MIC® Gastro-Enteric tube may be used for gastric decompression. Use gravity drainage, low continuous or low intermittent suction.

⚠ Caution: Never use high continuous or high intermittent suction. This could cause the gastric lumen to collapse.

Never use suction through the jejunal lumen of the MIC® Gastro-Enteric Tube.

MRI Safety Information

Non-clinical testing demonstrated the MIC® Gastro-Enteric (GE) Enteral Feeding Tube is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-T and 3-T, only
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode

Under the scan conditions defined, the MIC® Gastro-Enteric (GE) Enteral Feeding Tube is expected to produce a maximum temperature rise of 2.8°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

Artifact Information

In non-clinical testing, the image artifact caused by the MIC® Gastro-Enteric (GE) Enteral Feeding Tube extends approximately 2-mm from this implant when imaged using a gradient echo pulse sequence and a 3-T MR system.

Problem Solving

Problem	Cause	Remedy
1. Gastric leakage.	The ventilated ring is not properly adjusted or the balloon is not against the stomach wall. The balloon is under-inflated.	Check position of the ventilated SECUR-LOK® ring: 1–2 mm from the stoma, and the balloon is against the stomach wall. If the SECUR-LOK® ring is properly adjusted, add water in small increments (1–2 ml) to the balloon. Water can be added as the stomach size and balloon volume will allow: NOTE: 20 ml maximum for standard volume balloon inflation, 5 ml maximum for low volume balloon inflation. Make sure that the stoma is clean and dry. Wait 15 minutes and check to see if the problem has been corrected.
2. Feeding tube disconnects from the feeding port.	Build-up of oily deposits inside the tube.	Keep the feeding port connections clean. Remove oily build-up with either a diet soft drink or soap and water. Use a small amount of solution on a cotton tip swab and scrub the top two inches of the inside of the tube. ⚠ Warning: Be careful! Do not let alcohol run down the tube and into the stomach or bowel.
3. Tube decreasing in length.	Oily build-up on tube resulting in the SECUR-LOK® ring not holding.	Clean the tube. Wash with warm, soapy water; rinse and dry the tube. If oil remains, clean the tube with a small amount of isopropyl alcohol. ⚠ Warning: Do not let the alcohol contact the stoma. After cleaning, readjust the ventilated SECUR-LOK® ring 1–2 mm away from the stoma site, and verify that the balloon is up against the stomach wall. Check in one hour to be sure the ventilated SECUR-LOK® ring is holding the MIC® Gastro-Enteric Tube in place. If size continues to change, consult the physician.
4. Tube increasing in length.	Balloon under-inflated or balloon is leaking.	Withdraw water from the balloon and compare the volume against the amount originally used to fill the balloon. Fill again to the original volume. Recheck the volume periodically. If size continues to change, the tube should be removed.
5. Balloon fails to deflate.	Valve is clogged.	Cut off feeding port adapter well above the SECUR-LOK® ring. Allow the water in the balloon to empty.
6. Patient vomits and/or aspirates feeding formula.	Jejunal tube may have dislodged. Formula rate may be set too high.	Have the patient or care giver consult the physician immediately. Adjust flow rate. Reposition or replace the MIC® Gastro-Enteric tube.
7. Formula will not flow freely through the feeding port.	Jejunal tube may be kinked or clogged.	Irrigate with 20–30 ml of warm water (adults) or 10–15 ml (children). If resistance is felt or water does not flow freely, call physician immediately . Tube may have to be repositioned or replaced.
8. Gastric lumen not vented properly.	Possible obstruction in gastric lumen.	Irrigate with 20–30 ml of warm water (adults) or 10–15 ml (children). If resistance is felt or water does not flow freely, call physician immediately . Tube may need to be replaced.

DO'S

- **DO** irrigate the tube with each use; no less than 3 to 5 times a day.
- **DO** keep the gastric and jejunal ports clean and clear.
- **DO** keep the exterior of the tube clean.
- **DO** adjust the ventilated SECUR-LOK® ring properly (1–2 mm from the stoma site).
- **DO** use commercially prepared liquid medications.
- **DO** use a feeding pump.
- **DO** use low continuous or low intermittent suction through the gastric lumen.

DON'TS

- **DON'T** secure the tube in place with tape or other appliances.
- **DON'T** tape or dress a healthy stoma.
- **DON'T** let feeding solutions stand in your tube.
- **DON'T** use any diet not prescribed by your physician.
- **DON'T** cover the ventilated ring or healthy skin with tape or gauze.
- **DON'T** use a needle to inflate the balloon.
- **DON'T** push the ventilated ring onto the abdomen.
- **DON'T** put pressure on the stoma.
- **DON'T** clamp the tube.
- **DON'T** use oil-based lubricants or petroleum jelly on the tube.
- **DON'T** used crushed medications.
- **DON'T** used blenderized diets.
- **DON'T** use gravity feedings when feeding through the jejunal lumen.
- **DON'T** suture the ventilated SECUR-LOK® ring in position.
- **DON'T** use suction through the jejunal lumen.

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[illegible]

[illegible]

Distributed in the USA by Avanos Medical Sales, LLC,
5405 Windward Parkway, Alpharetta, GA 30004 USA.
In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com
 Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.
 Avanos Medical Belgium BV,
Leonardo da Vincilaan 1, 1930 Zaventem, Belgium.
Sponsored in Australia by Avanos Medical Australia Pty Ltd,
475 Victoria Avenue, Chatswood, NSW 2067, Australia.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018-2024 AVNS. All rights reserved. 2024-10-25
15-M1-155-02 / 50002132

CE 2797