
AVANOS*


MIC*
J FEEDING TUBE

WITH ENFit® CONNECTOR

**PATIENT CARE
INSTRUCTIONS**

MIC*

Patient Care Instructions

The MIC® Jejunal (J) Feeding Tube is a single unit feeding device. It enters the stomach through a gastric stoma. A balloon inside the stomach holds the tube in place on the inside and the external retention ring holds the tube in place on the outside. The access port is labeled Jejunal and is used to deliver nutrients to the small intestine. It contains an ENFit® connector that fits with ENFit® feeding sets and syringes. A second port labeled Bal. is used to inflate and deflate the balloon.

Key Terms

- | | |
|-----------------------------|----------------------------|
| (A) Jejunal Access Port | (B) Balloon Inflation Port |
| (C) External Retention Ring | (D) Retention Balloon |
| (E) Radiopaque Stripe | (F) Exit Port |

Tube Flushing

Flush the feeding tube with water every 4–6 hours during continuous feeding, anytime the feeding is interrupted, at least every 8 hours if the tube is not being used, or per clinician's instructions. Flush the feeding tube before, after, and in between medication administration. Avoid using other liquids such as cranberry juice and cola beverages to flush feeding tubes.

- Use a 30 to 60 mL ENFit® syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
- Use room temperature water for tube flushing.
- Always flush the tube with the amount of water prescribed by the clinician.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and may cause injury.

Jejunal Feeding

⚠ Warning: Never connect the Jejunal access port to suction. Do not measure residuals from the Jejunal access port. Never attempt to bolus feed through the jejunal lumen.

1. Assemble the equipment: Formula; Feeding set; ENFit® syringe; Enteral feeding pump; Water for flushing the tube.
2. Wash your hands with soap and water. Rinse and dry them thoroughly.
3. Before starting the feeding, it is important to verify that your tube is in the correct position. Look at the centimeter numbers on the tube. Compare them to the numbers recorded in the information section. If they are different, adjust the tube and the ring. Refer to the tube care section for more detail.
4. Shake the formula and wipe the container tops before opening. If you do not use all the formula, cover the open

container, record the date and time, and store it in the refrigerator. Throw away formula more than 24 hours old. Do not mix new formula with old formula. There is always a chance it will spoil.

5. Pour the formula into the feeding bag.
6. Run formula through the feed set tubing to reduce the amount of air entering the stomach. Attach the tubing to the enteral feeding pump. Follow the manufacturer's directions to set the pump.
7. Open the cap to the Jejunal access port of the MIC® J Feeding Tube.
8. Use an ENFit® syringe to flush the tube through the Jejunal access port with the prescribed amount of water (see the Tube Flushing section above).
9. Securely connect the ENFit® feed set connector from the feed set to the Jejunal access port of the MIC® J Feeding Tube. Take care not to over-tighten the feed set connector to the access port.
10. If the feed set tubing has a clamp, open it.
11. Start the feeding pump.
12. Disconnect the ENFit® feed set connector from the Jejunal access port.
13. Use an ENFit® syringe to flush the tube through the Jejunal access port with the prescribed amount of water (see the Tube Flushing section above).
14. Close the cap to the Jejunal access port of the MIC® J Feeding Tube.

Note: If you see formula in the gastric drainage, stop the feeding and notify the specialist.

Medication Administration

Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric coated medication or mix medication with formula.

1. Assemble the equipment: Medication; ENFit® syringe; Water for flushing the tube.
2. Wash your hands with soap and water. Rinse and dry them thoroughly.
3. Before administering medication, it is important to verify that your tube is in the correct position. Look at the centimeter numbers on the tube. Compare them to the numbers recorded in the information section. If they are different, adjust the tube and the ring. Refer to the tube care section for more detail.

4. Prepare the medication as prescribed by the pharmacist.
5. Fill an ENFit® syringe with the proper amount of medication.
6. Open the cap to the Jejunal access port of the MIC* J Tube.
7. Use an ENFit® syringe to flush the tube with the prescribed amount of water (see the Tube Flushing section above).
8. Securely connect the ENFit® syringe containing the medication to the access port. Take care not to over-tighten the syringe to the access port.
9. Deliver the medication by depressing the ENFit® syringe plunger.
10. Remove the syringe from the access port.
11. Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
12. Close the cap to the access port.

Stoma Care

As the stoma heals, carefully cleanse the site with normal saline three times daily to remove the small amount of mucus that normally accumulates around the stoma. If the mucus dries, it may irritate the skin. Applying a small amount of water may loosen the dried material.

After the stoma heals, a thorough cleansing with mild soap and water is best. A dressing is not necessary, and contributes to skin maceration and infection. Unless recommended by your specialist, avoid medicated ointments or powders. Do not use mineral oil or petroleum based products. DO NOT ROTATE THE TUBE. Jejunal tubes must NOT be rotated because they will kink. If a kink in the tubes occurs, consult physician.

Tube Care

1. **Balloon Maintenance:** Check the water volume in the balloon once a week. Insert a male Luer syringe into the Balloon Inflation Port (labeled BAL). Withdraw the water from the balloon. The amount of water you withdraw less water than the prescribed amount, replace the amount recommended by your specialist. Never add water before you withdraw all of the water from the balloon. Do not exceed the prescribed amount of water.

Note: *The use of a Luer-lock syringe may make this more difficult. A Luer-slip syringe is recommended.*

2. **Change in Tube Length:** Check the centimeter number on the tube length before every feeding and during tube care. Compare the position of the tube with the number recorded in the information section. If the length is too short or too long, wash the tube and both sides of the retention ring with soap and water. Adjust the ring and the tube to the original position. The ring should be 1–2 mm from the stoma to anchor the balloon against the inside stomach wall.

3. **Tube Blockage:** Obstructions are may be caused by formula or other accumulations (body fluids, medications) inside the tube. Prevent this by thoroughly flushing the tube as described in the Tube Flushing section above. If the flow is sluggish, gently inject the water, and let it flow back out of the tube until the obstruction has cleared. When you feel too much resistance and cannot inject water at all, the tube may need to be repositioned or replaced because it may be kinked. High pressure can rupture the tube. If the obstruction is above the skin line, massage the tube between your fingers to break it up. If none of these measures work, call your specialist. Inserting any foreign object into the tube may puncture it.
4. **Granulation Tissue:** A small amount of extra tissue around the stoma site is simply the body's attempt to heal the wound. This is normal. It is also normal for the extra tissue to secrete a small amount of mucus. If the tissue bleeds excessively, drains, or interferes with positioning the external retention ring, call your specialist.
5. **Unusual Symptoms:** Consult your health care practitioner if any of the following symptoms occur:
Abdominal pain • Abdominal discomfort • Abdominal tenderness • Abdominal distension • Dizziness or fainting • Unexplained fever • Unusual amount of bleeding through or around the tube.

DOs

- DO** Irrigate the Jejunal access port every six hours with water.
- DO** Keep the Jejunal access port clean.
- DO** Wash the tube and external retention ring often.
- DO** Adjust the external retention ring 1–2 mm from the skin.
- DO** Use commercially prepared formula for jejunal feeding.
- DO** Use a feeding pump to give formula.

DON'Ts

- DON'T** Let formula sit in the tube.
- DON'T** Use any diet not prescribed by your Health Care Provider.
- DON'T** Fill the balloon with air, medications, or food.
- DON'T** Push the external retention ring tightly onto the abdomen.
- DON'T** Clamp the tube with instruments that may damage the tube.
- DON'T** Use oil based lubricants and creams.
- DON'T** Use gravity or bolus feedings through the jejunal port.
- DON'T** Use suction on the jejunal port or check for residuals.
- DON'T** Rotate the tube in a circle.
- DON'T** Tape the ring to the skin.
- DON'T** Attempt to clear an obstruction with a wire or other foreign object.

Problem Solving

Problem	Cause	Remedy
1. Gastric leakage.	The External Retention Ring is not properly adjusted or the balloon is not against the stomach wall. The balloon is under-inflated.	Adjust the External Retention Ring. It should be 1–2 mm (about the thickness of a dime) from the skin. Once you have adjusted the ring, check the water volume in the balloon. Do not add more water to the balloon before you first withdraw all of the water originally used to fill the balloon. Note: <i>The maximum balloon volume is 10 ml. Do not exceed that amount.</i> Make sure that the skin around the tube (stoma) is clean and dry. Wait 15 minutes and check to see if the problem has been corrected. If the problem has not been corrected, notify your physician.
2. External Retention Ring slips.	Build-up of skin oils or other lubricants on the tube.	Clean the tube and the External Retention Ring with warm, soapy water. After cleaning, re-adjust the External Retention Ring, 1–2 mm (about the thickness of a dime) away from the skin, and make sure the balloon is up against the stomach wall. DO NOT APPLY EXCESSIVE TENSION. Check in one hour to be sure that the External Retention Ring is still holding the tube properly in place. If the problem has not been corrected, notify your physician.
3. Tube decreasing in length.	Oily build-up on tube resulting in slippage of the External Retention Ring.	Clean the tube and the External Retention Ring following the instructions in #2 above.
4. Tube increasing in length.	Balloon under-inflated or inoperative.	Withdraw water from the balloon and compare the volume against the amount originally used to fill the balloon. If the balloon is under-inflated, fill again to the original, prescribed volume. Recheck the volume periodically. If size continues to change, notify your physician.
5. Balloon leak or rupture.	Variable.	Notify your physician.
6. Patient vomits and/or aspirates stomach contents.	Jejunum lumen in the stomach instead of the jejunum. Formula flow rate too high.	Notify your physician.
7. Formula will not flow freely through jejunal tube.	Jejunum lumen kinked or clogged.	Irrigate with 20–30 ml water (adults), or 10–15 ml (children). If resistance is felt or water does not flow freely, notify your physician immediately.

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*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2019-01-25
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.
Sponsored in Australia by Avanos Medical Australia Pty Ltd, 475 Victoria Avenue,