


## ErgoFit® Limb sling

The ErgoFit® limb sling for body parts is used to raise a leg or an arm in order to help with blood circulation.



# Product	Description	Dimensions	
		in	cm
300.20005	Limb sling	A	22.8 58
		B	8.3 21

### General Characteristics

- Guarantee against manufacturing defects: 1 year
- Maximum loading capacity: 300 lb / 136 kg / 21.43 st
- Machine washable; easy to care for and clean.

### Technical Specifications

Conform to ISO 10535 & CE marking

\*Conditionally to a normal use of an appropriate BHM Medical patient lift and maintenance in accordance with the labels.