

**permobil**

THE POWER OF MOBILITY™



# READY OR NOT...

Designed to grow with a child for years to come, the Koala Miniflex is the ideal pediatric wheelchair for toddlers. While it is scaled to fit the needs of children, the Koala is built with the same rugged motor and front-wheel-drive technology as Permobil's full-size adult chairs. The low seat-to-floor height allows little users to interact at the same height as their peers, and the integrated elevation provides them with easy access to a variety of surface heights such as tables and desks. Plus, the tilt feature offers improved stability and rest, a necessity after a long day of play.

***Koala***  
MINIFLEX



KOALA012213



### PILOT+ ELECTRONICS

Programmable Pilot+ series electronics come standard



### POWER ELEVATOR

4" power seat elevator available to help reach any object



### SAFETY FIRST

Optional bumper available - the ideal safety accommodation for beginner drivers



### REMOTE STOP

Optional remote stop control gives any parent that extra layer of protection.

### AVAILABLE SHROUD COLORS



ONYX BLACK



TITAN BLUE



CHILI PEPPER RED



POP STAR PINK



PLATINUM SILVER

### SPECIFICATIONS

Maximum User Weight	110 lbs	50 kg	Seat to Floor Height	(Elevator)	15"	375 mm
Maximum Speed Options	4 MPH	6.5 km/h		(Seat Tube)	No Seat Tube Option	
Driving Range <sup>1</sup>	15 Miles	25.5 km	Tilt Options	(Power)	0° – 25°	
Base Width	23"	585 mm		(Manual)	No Manual Option	
Base Length	32"	810 mm	Recline Options	(Power)	No Power Option	
Minimum Turning Radius	21"	535 mm		(Manual)	85° – 120°	
Ground Clearance	2.5"	65 mm	Legrest Elevation Options	(Power)	No Power Option	
Optional Seat Elevator	4"	100 mm		(Manual)	90° – 150°	
Battery Type	M40		Armrest Pad Lengths	8", 10"		200, 260 mm
Drive Electronics	Pilot+		Seat Widths	10", 12", 14"		260, 310, 360 mm
Maximum Slope <sup>2</sup>	10° Up / 6° Down		Seat Depths	10", 12", 14"		260, 310, 360 mm
			Backrest Heights	14", 16", 20"		360, 410, 500 mm

<sup>1</sup> ACTUAL DRIVING RANGE WILL VARY BASED ON DRIVING AND BATTERY CONDITIONS.  
<sup>2</sup> REFER TO OWNER MANUAL FOR SAFE DRIVING PRACTICES