


PART NUMBERS		
SIZE	STANDARD	DYNAMIC
XX-SMALL	SH250XXS	SH270XXS
X-SMALL	SH250XS	SH270XS
SMALL	SH250S	SH270S
MEDIUM	SH250M	SH270M
LARGE	SH250L	SH270L
X-LARGE	SH250XL	SH270XL

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

WITH GROMMETED TOP STRAP



-B1




FLAT-MOUNT END FITTING 2x
HKA1-R-1

-B2



CINCH-MOUNT END FITTING 2x
HKA2-R-1

WITH PLAIN WEBBING TOP STRAP



-C2-B1



FLAT-MOUNT END FITTING 2x
HKA1-R-1

-C2-B2



CINCH-MOUNT END FITTING 2x
HKA2-R-1

note

description

PIVOTFIT HARNESS, FRONT-PULL, TECHNICAL SPECS

drawing no.

D13-0417-30

date originated

4/17/2013

drawn by

D. HAGAN

rev.

C



MAXIMUM LENGTH

XX-SMALL	9cm (3-1/2")
X-SMALL	
SMALL	11.5cm (4-1/2")
MEDIUM	
LARGE	14cm (5-1/2")
X-LARGE	

MINIMUM LENGTH

1cm (1/2")
ALL SIZES

PAD LENGTH

XX-SMALL	24.5cm (9-5/8")
X-SMALL	28cm (11")
SMALL	32.5cm (12-3/4")
MEDIUM	37cm (14-5/8")
LARGE	42.5cm (16-3/4")
X-LARGE	49cm (19-1/4")

LOWER STRAP LENGTH

XX-SMALL	40.5cm (16")
X-SMALL	
SMALL	51cm (20")
MEDIUM	
LARGE	61cm (24")
X-LARGE	


COMBINED PAD WIDTH
RANGE

XX-SMALL	11.5 - 19cm (4-1/2 - 7-1/2")
X-SMALL	12.5 - 20.5cm (5 - 8")
SMALL	14 - 24cm (5-1/2 - 9-1/2")
MEDIUM	15 - 25.5cm (6 - 10")
LARGE	17 - 30cm (6-3/4 - 11-3/4")
X-LARGE	19 - 32cm (7-1/2 - 12-1/2")


AVAILABLE TOP STRAPS

ALL SIZES

GROMMETTED TOP STRAP
SH025



PLAIN WEBBING TOP STRAP
SH024



note

description

PIVOTFIT HARNESS, FRONT-PULL, TECHNICAL SPECS

drawing no.

D13-0417-30

date originated

4/17/2013

drawn by

D. HAGAN

rev.

C