

Frequently Asked Questions

Written by Courtney Bed, Inc.

Do twin size sheets fit the mattress?

The mattress used is found in a standard hospital bed (80" by 36") and deep-pocket twin size sheets will fit best.

If I have a special requirement needed in the bed design... can it be modified?

Our beds are hand-crafted, custom-made. We try to accommodate individual needs; we recommend you contact us with any questions.

How do I clean a panel if it becomes soiled?

A panel can be cleaned using warm soapy water and a cloth. Some parents have also used a hand-held steam cleaner to do the job. Panels can also be removed and placed in a washing machine on normal cycle.

My child is a jumper... wills the bed stand up to such activity?

Yes, our beds have a low center of gravity and very stable footprint that can support up to 250 pounds. Our sidewalls are made from a material that has a high burst rating and are flexible yet strong.

I don't own any tools, what will I need to put the bed together?

The assembly of the Courtney Bed requires just two Allen wrenches that are included with the bed.
