

# Molift RgoSling StandUp for Quick Raiser

THE IDEAL SLING SOLUTION FOR ACTIVE TRANSFERS WITH MOLIFT QUICK RAISER

## Molift RgoSling StandUp

Molift RgoSling StandUp is a padded sling developed for active hoisting together with Molift Quick Raiser. It provides perfect lumbar support and is the ideal sling solution for hoisting and moving between different sitting positions. The RgoSling StandUp is quickly and easily placed around the user.

By virtue of the design it is ideal for use in toilet situations. Molift RgoSling StandUp has an efficient anti-slip lining to prevent the sling from sliding up. The wide back-support distributes weight and pressure evenly over a large surface which makes it more comfortable for the user. The waist belt provides additional support and stability.



## Molift RgoSling StandUp with Support

Molift RgoSling StandUp with Support provides additional pelvic support for users who need a little extra help to stand up.



### Description

Molift RgoSling StandUp

Size	Item no.
XS	1720710
S	1720720
M	1720730
L	1720740
XL	1720750

### Waist Size

Size	Waist Size (cm)	Waist Size (in)
XS	55–60 cm	(21.7"–23.6")
S	60–75 cm	(23.6"–29.5")
M	75–100 cm	(29.5"–39.4")
L	100–150 cm	(39.4"–59.1")
XL	150–175 cm	(59.1"–68.9")



### Description

Molift RgoSling StandUp with Support

Size	Item no.
XS	1720810
S	1720820
M	1720830
L	1720840
XL	1720850

### Material

Polyester

### Inspection instruction

Inspect the slings regularly, especially after washing. Carefully examine slings for wear and tear. Pay particular attention to seams, fabric, straps and loops. Any damaged slings should be replaced immediately.

