



Our kids are our passion.™

NXT
Generation NEXT

Pediatric Manual Wheelchair



Yes, you can.™

Date: _____ Buyer: _____
 Company: _____
 Marked For: _____
 Email: _____

Quote Order PO# _____
 Location: _____
 Phone: _____
 Fax: _____

Description	Part No.	Retail	Description	Part No.	Retail
<h3>Transport Optional Base</h3>			<h3>Vent & Battery Trays</h3>		
<input type="radio"/> NXT Pediatric Folding Transport Optional (E1234) NXT-30300 \$2800 Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.			Note: All vent and battery trays not with transport. Note: All vent and battery trays not with short frames. Note: All vent trays on NXT only with Dual Lever Hand Activated Tilt. Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.		
<input type="radio"/> NXT Pediatric Rigid Transport Optional (E1233) NXT-30310 \$2700 Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.			<input type="radio"/> Vent Tray - Fixed (E1029) 11-7000 \$ 439 <input type="radio"/> Battery Tray - Fixed (K0108) 11-8000 \$ 429 <input type="radio"/> Quick Release Tray Hardware (E1028) 563 \$ 272		
<h3>Transportation Option</h3>			Note: For LP10 style vent and battery tray on frames 16" wide and greater. Note: Vent tray is 4" narrower than frame width and 16" deep. Note: Battery tray is 4" narrower than frame width and 10" deep.		
<input type="radio"/> Transportation Package NXT-99998 \$ 250 Note: Used when Transport Optional frame is selected above. Note: Must be installed by Freedom Designs manufacturing. Note: Not available with all vent and battery trays. Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved. Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.			<h3>Frame Colors</h3>		
<h3>Transport Included Base</h3>			<input type="radio"/> Sugar Cookie n/c <input type="radio"/> Pacific Blue n/c <input type="radio"/> Neon Yellow n/c <input type="radio"/> Blue Blast n/c <input type="radio"/> Yellow n/c <input type="radio"/> Electric Blue n/c <input type="radio"/> Bright Pink n/c <input type="radio"/> Galaxy Blue n/c <input type="radio"/> Pink Parfait n/c <input type="radio"/> Passion Purple n/c <input type="radio"/> Sugar Plum n/c <input type="radio"/> Octane Orange n/c <input type="radio"/> Vigorous Violet n/c <input type="radio"/> Red Sunset n/c <input type="radio"/> Candy Violet n/c <input type="radio"/> Really Red n/c <input type="radio"/> Candy Raspberry n/c <input type="radio"/> Fiery Red n/c <input type="radio"/> Extreme Green n/c <input type="radio"/> Twilight n/c <input type="radio"/> Sparkle Green n/c <input type="radio"/> Stormy Grey n/c <input type="radio"/> Mystic Teal n/c <input type="radio"/> Black n/c <input type="radio"/> Blizzard Blue		
<input type="radio"/> NXT Pediatric Folding Transport Included (E1234) NXT-30301 \$2950 Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.			<h3>Frame Dimensions</h3>		
<input type="radio"/> NXT Pediatric Rigid Transport Included (E1233) NXT-30311 \$2850 Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.			Specify Frame Width Setting		
<h3>Seat Rail Options</h3>			<input type="radio"/> 10" <input type="radio"/> 11" <input type="radio"/> 12" Adjusts 10" to 12" NXT-30250 n/c <input type="radio"/> 11" <input type="radio"/> 12" <input type="radio"/> 13" Adjusts 11" to 13" NXT-30210 n/c <input type="radio"/> 12" <input type="radio"/> 13" <input type="radio"/> 14" Adjusts 12" to 14" NXT-30220 n/c <input type="radio"/> 13" <input type="radio"/> 14" <input type="radio"/> 15" Adjusts 13" to 15" NXT-30230 n/c <input type="radio"/> 14" <input type="radio"/> 15" <input type="radio"/> 16" Adjusts 14" to 16" NXT-30240 n/c <input type="radio"/> 15" <input type="radio"/> 16" <input type="radio"/> 17" Adjusts 15" to 17" NXT-30290 n/c <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" Adjusts 16" to 18" NXT-30260 n/c <input type="radio"/> 17" <input type="radio"/> 18" Adjusts 17" to 18" NXT-30261 n/c <input type="radio"/> 18" Not width adjustable NXT-30262 n/c <input type="radio"/> 19" Not adjustable, not w/ short frame NXT-30263 n/c <input type="radio"/> 20" Not adjustable, not w/ short frame (E2201) NXT-30264 \$ 400		
Note: Not available with ALL Contracture, 2D & 3D Dynamic footrest systems.			Note: Frame widths measured across and from outside of seat frame tubes. Note: Seat frame widths of 19" and 20" not available with short frame. Note: Pediatric Long frame only with 15" wide and greater frames. Note: Vent and battery trays only with 14" and wider frames.		
<input type="radio"/> NXT LEAF Seat Frame- Pair (K0108) NXT-30324 \$ 714 <input type="radio"/> NXT LEAF Seat Frame- Left (K0108) NXT-30324-L \$ 357 <input type="radio"/> NXT LEAF Seat Frame- Right (K0108) NXT-30324-R \$ 357 Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy. Note: Freedom customized seats can be ordered to fit pre-configured LEAF frame models Note: No seating system is provided with NXT frames. A seating system must be purchased separately.			<h3>NXT Rider Weight Capacity is 225 lbs.</h3>		
<input type="radio"/> NXT SRO Seat Frame- Pair (K0108) NXT-30325 \$ 306 <input type="radio"/> NXT SRO Seat Frame- Left (K0108) NXT-30325-L \$ 153 <input type="radio"/> NXT SRO Seat Frame- Right (K0108) NXT-30325-R \$ 153 Note: SRO is manufactured without standard footrest hanger attachments. Note: SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.			Page 1 of 7		
<h3>Tilt Activations</h3>			NXT Pediatric		
<input type="radio"/> Dual Lever Hand Activated Tilt NXT-30272 n/c Note: Must be selected for ALL vent and battery trays.			€		
<input type="radio"/> Foot Activated Tilt NXT-30273 n/c Note: Foot pedal installed on right rear side frame. Note: Not with ALL vent and battery trays.			WONXTP 032016		
<input type="radio"/> Combination Dual Hand Activated Tilt, Brake and Wheel Lock System NXT-30277 n/c Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" Low Profile Poly Mag Wheels					

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Back Post Locations
Specify Frame Type and Frame Length

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

- Short Frame & Short Tilt Bracket** NXT-30511 n/c
 - 12" NXT-30504 n/c
 - 13" NXT-30505 n/c
 - 14" NXT-30506 n/c
 - 15" NXT-30521 n/c
 - 16" NXT-30522 n/c
 - * 17" Must select 2" extension NXT-30523 n/c
 - * 18" Must select 2" extension NXT-30524 n/c
 - 2" Extension - Pair (K1008) NXT-30525 \$ 26

- Specify Back Post Location** (if other than specified above) _____ "
- Note:** Floor to seat frame height range is 12.5 to 16.5".
- Note:** Seat frame widths of 19" and 20" not available with short frame.
- * **Note:** Must choose 2" extension for 17" and 18" frame lengths.
- Note:** 22" & 24" wheels not available with short frame and short tilt bracket.
- Note:** All vent and battery trays not available with short frame.

- Medium Frame & Short Tilt Bracket** NXT-30512 n/c
 - 14" NXT-30506 n/c
 - 15" NXT-30521 n/c
 - 16" NXT-30522 n/c
 - 17" NXT-30523 n/c
 - 18" NXT-30524 n/c
 - * 19" Must select 2" extension NXT-30527 n/c
 - * 20" Must select 2" extension (E2203) NXT-30528 \$ 500
 - * 21" Must select 3" extension (E2203) NXT-30529 \$ 500
 - 2" Extension - Pair (K1008) NXT-30525 \$ 26
 - 3" Extension - Pair (K1008) NXT-30530 \$ 26

- Specify Back Post Location** (if other than specified above) _____ "
- Note:** Floor to seat frame height range is 12.5" to 16.5".
- * **Note:** Must choose 2" extension for 17" and 18" frame lengths.
- Note:** 22" & 24" wheels not available with medium frame and short tilt bracket.

- Medium Frame & Tall Tilt Bracket** NXT-30540 n/c
 - 14" NXT-30506TP n/c
 - 15" NXT-30521TP n/c
 - 16" NXT-30522TP n/c
 - 17" NXT-30523TP n/c
 - 18" NXT-30524TP n/c
 - * 19" Must select 2" extension NXT-30527TP n/c
 - * 20" Must select 2" extension (E2203) NXT-30528TP \$ 500
 - * 21" Must select 3" extension (E2203) NXT-30529TP \$ 500
 - * 22" Must select 4" extension (E2204) NXT-30535TP \$ 900
 - 2" Extension - Pair (K1008) NXT-30525 \$ 26
 - 3" Extension - Pair (K1008) NXT-30530 \$ 26
 - 4" Extension - Pair (K1008) NXT-30531 \$ 26

- Specify Back Post Location** (if other than specified above) _____ "
- Note:** Tall tilt bracket increases floor to seat rail height 1.5"
- Note:** Floor to seat frame height range is 14" to 18".
- * **Note:** Must choose extension for 19", 20", 21" and 22" frame lengths.

- 5° Anterior Tilt Seat Frame** NXT-30276 n/c
- Note:** Limits posterior tilt to 40°, includes split collars for adjustment.
- Note:** Medium with Tall Tilt and Long frames only.
- Note:** NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14"

Notes:

Back Post Tubes

Specify Push Handle, Back Post and Extensions

- Standard 90° Bend Push Handles and Back Post Tubes** NXT-30518 n/c
- Ergo Bend Push Handles and Back Post Tubes** NXT-30519 n/c
- Note:** All push/stroller handle extensions not recommended with Ergo Bend.
- Seating Dynamics Rocker Back** (K0108) NXT-30532 \$1734
- Note:** Angle adjustable from 72° to 108° in 6° increments.
- Note:** Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.
- Note:** Locks into a non-dynamic position when necessary.
- Note:** Does not allow the back post tubes to fold forward.

Back Post Height

- Note:** Back Posts are measured from top of seat rail to top of handle bend.
- Short Push Handle - 20" high** NXT-30515 n/c
 - Medium Push Handle - 23" high** NXT-30514 n/c
 - Tall Push Handle - 26" high** NXT-30516 n/c

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Legrests, Footrests and Extensions

UPDATED!

Specialized Footrest Systems

Note: ALL legrest measurements are from the seat rail to footplate.
Note: Not available with LEAF or SRO seat rail options.

Super Contracture Footrest Hanger (K0108)

- Short Pair - 11" to 15"** NXT-30291-2 **\$3570**
- Long Pair - 13" to 18"** NXT-30292-2 **\$3570**

Note: Super Contracture for frames 15" wide and greater.
Note: Super Contracture is both swing-away and removable.
Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.
Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflexion positioning.
Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.
Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

2D Dynamic Footrest Hanger (K0108)

- Short Pair - 14" to 16"** NXT-30293-2 **\$4080**
- Long Pair - 16" to 18"** NXT-30294-2 **\$4080**

Note: 2D Dynamic for frames 15" wide and greater.
Note: 2D Dynamic is both swing-away and removable.
Note: 2D Dynamic has adjustable knee joints that have 48 positions that adjust in 7.5° increments.
Note: 2D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.
Note: 2D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
Note: 2D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit 2D Dynamic adjustment range and/or cause footplate/front caster interference.

3D Dynamic Footrest Hanger (K0108)

- Short Pair - 13" to 15"** NXT-30295-2 **\$4590**
- Long Pair - 15" to 17"** NXT-30296-2 **\$4590**

Note: 3D Dynamic for frames 15" wide and greater.
Note: 3D Dynamic is both swing-away and removable.
Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.
Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.
Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

Armrests and Armrest Pads

- Omit Armrests** NXT-38209

Back Post Mounted Armrests

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

- 10.5" Long Flip-Back** NXT-38204 **n/c**
Note: Full length pad on 13.5" or longer armrests only.
- 13.5" Long Flip-Back** NXT-38202 **n/c**
Note: Full length pad on 13.5" or longer armrests only.
- 10.5" Long Flip-Back Locking (E0973)** NXT-38201 **\$ 245**
Note: Full length pad on 13.5" or longer armrests only.
- 13.5" Long Flip-Back Locking (E0973)** NXT-38200 **\$ 245**
Note: Full length pad on 13.5" or longer armrests only.

Armrests and Armrest Pads (continued)

- 10.5" to 13.5" Long Flip-Back, Telescoping** (E0973) NXT-38208 **\$ 240**
Note: Full length pad on 13.5" or longer armrests only.
- 10.5" to 13.5" Long Flip-Back Locking, Telescoping** (E0973) NXT-38210 **\$ 265**
Note: Full length pad on 13.5" or longer armrests only.
- 17" to 20" Long Flip-Back Locking, Telescoping** (E0973) NXT-38215 **\$ 265**

- Foam Wraps - Standard** NXT-38216 **n/c**
- Desk Length Pads** NXT-38213 **n/c**
- Full Length Pads** NXT-38214 **n/c**
Note: 8.5" armrests not available with desk and full length armrest pads.
Note: Full length pad on 13.5" or longer armrests only.

Frame Mounted Armrests

- Single Post Fixed Height, Removable - Short** (K0015) NXT-38222 **\$ 225**
- Single Post Fixed Height, Removable - Tall** (K0015) NXT-38224 **\$ 225**
Note: Single Post Fixed Height arms include a removable side guard.
Note: Adjustable in height with tools only. See matrix for height ranges.

- Single Post Height Adjustable, Removable - Short** (E0973) NXT-38217 **\$ 265**
- Single Post Height Adjustable, Removable - Tall** (E0973) NXT-38203 **\$ 265**

- Dual Post Fixed Height Removable** NXT-38205 **n/c**
Note: All dual post armrests not available with short frame.
Note: All dual post armrests not available with 22" or 24" wheels.
- Dual Post Flip-Back Height Adjustable, Removable** (E0973) NXT-38206 **\$ 275**
Note: All dual post armrests not available with short frame.
Note: All dual post armrests not available with 22" or 24" wheels.

- Desk Length Pads** NXT-38213 **n/c**
- Full Length Pads** NXT-38214 **n/c**

- Specify Seat Rail to Top of Armrest Measurement** _____ "

Front Wheels - Casters

- Omit Caster Wheels** NXT-34333
- 3" Polyurethane** NXT-34888 **\$ 56**
- 3" Polyurethane Light-Up** NXT-34111 **\$ 56**
Note: 3" casters recommended for low seat to floor heights.
- 4" Polyurethane** NXT-40000 **\$ 66**
- 4" Polyurethane Light-Up** NXT-41000 **\$ 101**
- 4" x 1.4" Poly, Aluminum Mag** NXT-43000 **\$ 92**
- 5" x 1" Poly, Plastic Mag** NXT-34555 **n/c**
- 5" Polyurethane Light-Up** NXT-35050 **\$ 101**
- 5" x 1.5" Poly, Plastic Mag** NXT-34222 **\$ 87**
- 5" x 1.4" Poly, Aluminum Mag** NXT-42000 **\$ 92**
- 6" x 1" Poly, Plastic Mag** NXT-34666 **n/c**
- 6" x 1.5" Poly, Plastic Mag** NXT-34999 **\$ 87**
- 6" x 1.4" Poly, Aluminum Mag** NXT-44000 **\$ 92**
- 6" x 2" Poly with Tread** NXT-34444 **\$ 92**
- 6" x 1" Pneumatic Plastic Mag (E2214/E2215)** NXT-34667 **\$ 77**
- 7" x 1" Poly, Plastic Mag** NXT-34777 **\$ 77**
- 7" x 1.5" Poly, Plastic Mag** NXT-34000 **\$ 118**
- 7" x 2" Polyurethane** NXT-34200 **\$ 118**
Note: All 7" casters not available with all 12" rear wheels.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheel Locks

- Push to Lock - Standard** (E0961) NXT-31501 n/c
- Pull to Lock** (E0961) NXT-31506 n/c
- Extension Handles** (E0961) NXT-31801 \$ 61
- Foot Activated Wheel Lock** (E2228) NXT-31508 n/c
 Note: Includes manual wheel locks to be used with public transportation.
 Note: Not available with 12" poly plastic mag wheels.
 Note: Not available with ALL vent and battery trays.
- Hand Activated Brake and Wheel Lock System** (E2228) NXT-30275 n/c
 Note: Includes manual wheel locks to be used with public transportation.
 Note: Not available with 12" mag wheels with low pro poly tires.
 Note: Used to control direction, speed, brake and lock wheels from strolling position.
 Note: Not available with all vent and battery trays.
- Combination Dual Hand Activated Tilt, Brake and Wheel Lock System** (E2228) NXT-30277 n/c UPDATED!
 Note: Includes manual wheel locks to be used with public transportation.
 Note: Not available with 12" poly plastic mag wheels.
 Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.
 Note: Not available with all vent and battery trays.

Rear Anti-Tippers

- Curb Assist Anti-Tip Assembly Flip-Up, Adjustable** (E0971) NXT-31600 \$ 97
 Note: Pair includes one anti-tip tube with curb assist feature.
- Standard Anti-Tip Assembly Flip-Up, Adjustable** (E0971) NXT-31603 \$ 97

Seat Pan Options

- Fixed Solid Seat Pan - Bolteron** (E2231) NXT-39010 \$ 150
 Note: Not available with folding frame.
- Fixed Solid Seat Pan - Aluminum** NEW! (E2231) NXT-39030 \$ 170
 Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".
 Note: Not available with folding frame.
- Removable Solid Seat Pan - Bolteron** (E2231) NXT-39012 \$ 170
- Removable Solid Seat Pan - Aluminum** NEW! (E2231) NXT-39032 \$ 255
 Note: Removable Seat Pan raises the effective seat to floor height by 5/8".
 Note: Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.
 Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.
 Note: All seat pans include 2" x 2" notches for growth adjustment.
 Note: All seat pans made to fit frame dimensions unless specified below.
- Seat Width _____" X Seat Depth _____"**

Accessories

- Medical Necessities Pouch** (E0950) NXT-23207 \$ 68
 Note: Black nylon.
- Footrest Hanger Pads - Pair** (E0950) SFDI-9HP \$ 92
 Note: Wrap-around pads made from black neoprene.
- Frame Pads - Pair** (E0950) 10911-2 \$ 82
 Note: Wrap-around pads made from black neoprene.
 Note: Typically used to pad exposed back post tubes.
- Flip-Down Rigidizing Bar** NEW! (E0950) NXT-23210 \$ 214
 Note: Mounts to back post tubes to add rigidity to frame.
 Note: Does NOT adjust in width to grow with frame.
- Freedom Backpack** NEW! (E0950) NXT-10900 \$ 56
 Note: Black Nylon only with Freedom Designs logo.
 Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip-Back** NEW! (E0950) NXT-10905 \$ 128
 Note: Canopy assembly includes removable mounting hardware.
 Note: Canopy is available in black only with Freedom Designs logo.
 Note: Canopy includes rear flap with zippered storage pocket.
 Note: Canopy includes clear-view window.

Anterior Support Surfaces – (Trays)

- 3/8" Lexan Tray - Channel Bracket Hardware** (E0950) 08-2000 \$ 271
- 3/8" Lexan Tray - Adjustable Knob Hardware & Ribs** (E0950) 08-6000 \$ 219
- 1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs** (E0950) 08-13000 \$ 183
- 1/4" Lexan Tray - Hook & Loop Straps** (E0950) 08-16000 \$ 106

Continue To Next Page

Date: _____ PO# _____					Marked For: _____																																																																																																																																																																																																																																																														
Description			Part No.		Retail		Description			Part No.		Retail																																																																																																																																																																																																																																																							
NXT Footrest / Extension Lengths						<p>HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.</p> <p><input type="radio"/> Choose NXT Seat to Floor Height _____"</p> <p>Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height. Note: Caster sizes vary slightly in height by type. Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8". Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".</p>																																																																																																																																																																																																																																																													
Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)																																																																																																																																																																																																																																																														
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"																																																																																																																																																																																																																																																														
High Mount	N/A	N/A	1.5"-3"	N/A	N/A																																																																																																																																																																																																																																																														
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"																																																																																																																																																																																																																																																														
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A																																																																																																																																																																																																																																																														
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A																																																																																																																																																																																																																																																														
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A																																																																																																																																																																																																																																																														
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A																																																																																																																																																																																																																																																														
NXT MAPS Footrest / Extension Lengths						<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th colspan="8" style="text-align: center;">NXT Wheel / Caster Seat to Floor Heights</th> </tr> <tr> <th>Wheel</th> <th>Caster</th> <th>Pivot</th> <th>Height</th> <th>Wheel</th> <th>Caster</th> <th>Pivot</th> <th>Height</th> </tr> </thead> <tbody> <tr><td>12"</td><td>3"</td><td>Short</td><td>12.5-13"</td><td>20"</td><td>3"</td><td>Tall</td><td>N/A</td></tr> <tr><td>12"</td><td>3"</td><td>Tall</td><td>14-14.5"</td><td>20"</td><td>4"</td><td>Short</td><td>14-15"</td></tr> <tr><td>12"</td><td>4"</td><td>Short</td><td>13.5"</td><td>20"</td><td>4"</td><td>Tall</td><td>15.5-16.5"</td></tr> <tr><td>12"</td><td>4"</td><td>Tall</td><td>15"</td><td>20"</td><td>5"</td><td>Short</td><td>14.5-15.5"</td></tr> <tr><td>12"</td><td>5"</td><td>Short</td><td>15.5"</td><td>20"</td><td>5"</td><td>Tall</td><td>16-17"</td></tr> <tr><td>12"</td><td>5"</td><td>Tall</td><td>17"</td><td>20"</td><td>6"</td><td>Short</td><td>15.5-16"</td></tr> <tr><td>12"</td><td>6"</td><td>Short</td><td>16"</td><td>20"</td><td>6"</td><td>Tall</td><td>17-17.5"</td></tr> <tr><td>12"</td><td>6"</td><td>Tall</td><td>17.5"</td><td>20"</td><td>7"</td><td>Short</td><td>16.5"</td></tr> <tr><td>16"</td><td>3"</td><td>Short</td><td>12.5-13"</td><td>20"</td><td>7"</td><td>Tall</td><td>18"</td></tr> <tr><td>16"</td><td>3"</td><td>Tall</td><td>14-14.5"</td><td>22"</td><td>3"</td><td>Short</td><td>N/A</td></tr> <tr><td>16"</td><td>4"</td><td>Short</td><td>13.5-14"</td><td>22"</td><td>3"</td><td>Tall</td><td>N/A</td></tr> <tr><td>16"</td><td>4"</td><td>Tall</td><td>15-15.5"</td><td>22"</td><td>4"</td><td>Short</td><td>N/A</td></tr> <tr><td>16"</td><td>5"</td><td>Short</td><td>14.5"</td><td>22"</td><td>4"</td><td>Tall</td><td>16.5"</td></tr> <tr><td>16"</td><td>5"</td><td>Tall</td><td>16"</td><td>22"</td><td>5"</td><td>Short</td><td>N/A</td></tr> <tr><td>16"</td><td>6"</td><td>Short</td><td>16</td><td>22"</td><td>5"</td><td>Tall</td><td>16.5-17"</td></tr> <tr><td>16"</td><td>6"</td><td>Tall</td><td>17.5</td><td>22"</td><td>6"</td><td>Short</td><td>N/A</td></tr> <tr><td>16"</td><td>7"</td><td>Short</td><td>16.5"</td><td>22"</td><td>6"</td><td>Tall</td><td>17-17.5"</td></tr> <tr><td>16"</td><td>7"</td><td>Tall</td><td>18"</td><td>22"</td><td>7"</td><td>Short</td><td>N/A</td></tr> <tr><td>18"</td><td>3"</td><td>Short</td><td>12.5-13"</td><td>22"</td><td>7"</td><td>Tall</td><td>18"</td></tr> <tr><td>18"</td><td>3"</td><td>Tall</td><td>14-14.5"</td><td>24"</td><td>3"</td><td>Short</td><td>N/A</td></tr> <tr><td>18"</td><td>4"</td><td>Short</td><td>13.5-14.5"</td><td>24"</td><td>3"</td><td>Tall</td><td>N/A</td></tr> <tr><td>18"</td><td>4"</td><td>Tall</td><td>15-16.5"</td><td>24"</td><td>4"</td><td>Short</td><td>N/A</td></tr> <tr><td>18"</td><td>5"</td><td>Short</td><td>14.5-15.5"</td><td>24"</td><td>4"</td><td>Tall</td><td>N/A</td></tr> <tr><td>18"</td><td>5"</td><td>Tall</td><td>16-17"</td><td>24"</td><td>5"</td><td>Short</td><td>N/A</td></tr> <tr><td>18"</td><td>6"</td><td>Short</td><td>15.5"</td><td>24"</td><td>5"</td><td>Tall</td><td>N/A</td></tr> <tr><td>18"</td><td>6"</td><td>Tall</td><td>17"</td><td>24"</td><td>6"</td><td>Short</td><td>N/A</td></tr> <tr><td>18"</td><td>7"</td><td>Short</td><td>N/A</td><td>24"</td><td>6"</td><td>Tall</td><td>17.5"</td></tr> <tr><td>18"</td><td>7"</td><td>Tall</td><td>N/A</td><td>24"</td><td>7"</td><td>Short</td><td>N/A</td></tr> <tr><td>20"</td><td>3"</td><td>Short</td><td>N/A</td><td>24"</td><td>7"</td><td>Tall</td><td>18"</td></tr> </tbody> </table>						NXT Wheel / Caster Seat to Floor Heights								Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height	12"	3"	Short	12.5-13"	20"	3"	Tall	N/A	12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"	12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"	12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"	12"	5"	Short	15.5"	20"	5"	Tall	16-17"	12"	5"	Tall	17"	20"	6"	Short	15.5-16"	12"	6"	Short	16"	20"	6"	Tall	17-17.5"	12"	6"	Tall	17.5"	20"	7"	Short	16.5"	16"	3"	Short	12.5-13"	20"	7"	Tall	18"	16"	3"	Tall	14-14.5"	22"	3"	Short	N/A	16"	4"	Short	13.5-14"	22"	3"	Tall	N/A	16"	4"	Tall	15-15.5"	22"	4"	Short	N/A	16"	5"	Short	14.5"	22"	4"	Tall	16.5"	16"	5"	Tall	16"	22"	5"	Short	N/A	16"	6"	Short	16	22"	5"	Tall	16.5-17"	16"	6"	Tall	17.5	22"	6"	Short	N/A	16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"	16"	7"	Tall	18"	22"	7"	Short	N/A	18"	3"	Short	12.5-13"	22"	7"	Tall	18"	18"	3"	Tall	14-14.5"	24"	3"	Short	N/A	18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A	18"	4"	Tall	15-16.5"	24"	4"	Short	N/A	18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A	18"	5"	Tall	16-17"	24"	5"	Short	N/A	18"	6"	Short	15.5"	24"	5"	Tall	N/A	18"	6"	Tall	17"	24"	6"	Short	N/A	18"	7"	Short	N/A	24"	6"	Tall	17.5"	18"	7"	Tall	N/A	24"	7"	Short	N/A	20"	3"	Short	N/A	24"	7"	Tall	18"
NXT Wheel / Caster Seat to Floor Heights																																																																																																																																																																																																																																																																			
Wheel	Caster	Pivot	Height	Wheel	Caster							Pivot	Height																																																																																																																																																																																																																																																						
12"	3"	Short	12.5-13"	20"	3"							Tall	N/A																																																																																																																																																																																																																																																						
12"	3"	Tall	14-14.5"	20"	4"							Short	14-15"																																																																																																																																																																																																																																																						
12"	4"	Short	13.5"	20"	4"							Tall	15.5-16.5"																																																																																																																																																																																																																																																						
12"	4"	Tall	15"	20"	5"							Short	14.5-15.5"																																																																																																																																																																																																																																																						
12"	5"	Short	15.5"	20"	5"							Tall	16-17"																																																																																																																																																																																																																																																						
12"	5"	Tall	17"	20"	6"							Short	15.5-16"																																																																																																																																																																																																																																																						
12"	6"	Short	16"	20"	6"							Tall	17-17.5"																																																																																																																																																																																																																																																						
12"	6"	Tall	17.5"	20"	7"							Short	16.5"																																																																																																																																																																																																																																																						
16"	3"	Short	12.5-13"	20"	7"							Tall	18"																																																																																																																																																																																																																																																						
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A																																																																																																																																																																																																																																																												
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A																																																																																																																																																																																																																																																												
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A																																																																																																																																																																																																																																																												
16"	5"	Short	14.5"	22"	4"	Tall	16.5"																																																																																																																																																																																																																																																												
16"	5"	Tall	16"	22"	5"	Short	N/A																																																																																																																																																																																																																																																												
16"	6"	Short	16	22"	5"	Tall	16.5-17"																																																																																																																																																																																																																																																												
16"	6"	Tall	17.5	22"	6"	Short	N/A																																																																																																																																																																																																																																																												
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"																																																																																																																																																																																																																																																												
16"	7"	Tall	18"	22"	7"	Short	N/A																																																																																																																																																																																																																																																												
18"	3"	Short	12.5-13"	22"	7"	Tall	18"																																																																																																																																																																																																																																																												
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A																																																																																																																																																																																																																																																												
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A																																																																																																																																																																																																																																																												
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A																																																																																																																																																																																																																																																												
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A																																																																																																																																																																																																																																																												
18"	5"	Tall	16-17"	24"	5"	Short	N/A																																																																																																																																																																																																																																																												
18"	6"	Short	15.5"	24"	5"	Tall	N/A																																																																																																																																																																																																																																																												
18"	6"	Tall	17"	24"	6"	Short	N/A																																																																																																																																																																																																																																																												
18"	7"	Short	N/A	24"	6"	Tall	17.5"																																																																																																																																																																																																																																																												
18"	7"	Tall	N/A	24"	7"	Short	N/A																																																																																																																																																																																																																																																												
20"	3"	Short	N/A	24"	7"	Tall	18"																																																																																																																																																																																																																																																												
NXT MAPS Footrest / Extension Lengths																																																																																																																																																																																																																																																																			
Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)																																																																																																																																																																																																																																																														
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"																																																																																																																																																																																																																																																														
High Mount	N/A	N/A	N/A	N/A	N/A																																																																																																																																																																																																																																																														
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"																																																																																																																																																																																																																																																														
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A																																																																																																																																																																																																																																																														
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A																																																																																																																																																																																																																																																														
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A																																																																																																																																																																																																																																																														
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A																																																																																																																																																																																																																																																														
<p>Note: Length from top of footrest plate to top of seat rail. Note: Without extensions bottom hanger hole not accessible for hardware. Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p>																																																																																																																																																																																																																																																																			
Freedom Designs Armrest Heights																																																																																																																																																																																																																																																																			
Armrest and Pad		Foam	Desk	Full																																																																																																																																																																																																																																																															
Flip Back		6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"																																																																																																																																																																																																																																																															
Flip Back Locking		6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"																																																																																																																																																																																																																																																															
Telescoping Flip Back		6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"																																																																																																																																																																																																																																																															
Flip Back Locking Telescoping		6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"																																																																																																																																																																																																																																																															
Single Post Height Adj. Short		N/A	8" - 14"	8" - 14"																																																																																																																																																																																																																																																															
Single Post Height Adj. Tall		N/A	10" - 16"	10" - 16"																																																																																																																																																																																																																																																															
Dual Post Fixed Height		N/A	9.5"	9.5"																																																																																																																																																																																																																																																															
Dual Post Flip Back Height Adj.		N/A	9" - 13"	9" - 13"																																																																																																																																																																																																																																																															
Note: All heights from top of seat frame to top of armrest pad.																																																																																																																																																																																																																																																																			