

Features:

- Constructed of the highest quality health care grade PVC
- Deluxe elongated extra heavy-duty, contoured open from seat
- Dual reinforcement into the caster fittings adds additional support
- Quick drying, removable mesh back support is standard in all shower chairs
- Available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors Royal Blue mesh is the default color shipped. If you desire another color or vinyl, please order via Customer Service at 1-800-347-5678
- Ergonomic curved arms provide support without sharp edges and include anti-slip handgrips for added safety
- Threaded stem casters for additional safety
- Optional telescoping footrest feature allows footrest to move in or out with ease
- Optional positioning belt available separately. (Item#89301)
- One year limited warranty on PVC frame; six month limited warranty on all other materials



with foot rest	89110	89210	89251	89351*	
without foot rest	89150		89250	89350*	
Internal Width	18"	18"	22"	26"	30"
External Width	22"	22"	26"	30"	34"
Overall Height	40"	40" - 45"	40"	40"	43"
Seat Height	21"	20" - 25"	21"	20"	23"
Seat Clearance	18"	18" - 23"	18"	16"	19"
Seat Depth	18"	18"	19 ½"	22"	21 ½"
Overall Depth	21"	21"	23"	30"	28 ¼"
Weight Capacity - EVENLY DISTRIBUTED	300 lb	250 lb	375 lb	425 lb	700 lb
Caster Style and Size	3" twin wheel	3" twin wheel	3" heavy duty	4" heavy duty	5" heavy duty

Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Mauve

*89350 and 89351 cannot be used over a toilet