

Lumex Three Position Recliner (Model 574G)

- The Standard of the Industry and the recliner that started it all! This first recliner specifically designed to improve the comfort and mobility of residents in extended care facilities. The Lumex Three Position Recliner easily adjusts into any one of three positions: sitting, TV or full recline. Recliners offer numerous therapeutic benefits including position changes that can aid in circulation and ease discomfort from respiratory, heart and other ailments.
- A position lock allows only the caregiver to secure the chair in any one of the three easily attained positions. Recliner encourages natural relaxation and improved blood circulation.
- Flex-o-Lator seat support system increases comfort by reducing the risk of “bottoming out”.
- **Warranty**
 1. Frame: Limited 2 Years
 2. Upholstery: Limited 1 Year
 3. All Other Components: Limited 1 Year
- Ergonomically designed seat and back with varied density foam construction provides correct support. Articulating headrest quickly adjusts to maximize individual comfort requirements.
- For safer, easier transport, the Lumex Three Position Recliner features two, four inch swivel casters with wheel locks and four inch fixed front wheels.
- The molded activity table can be fixed in position or left unsecured for easy exit from the recliner.
- Durable welded steel epoxy coated frame is easy to clean and maintain.
- Weight Bearing Capacity: 250 pounds.


Available Colors:

574G401	Vintage Gold	1 ea
574G409	Warm Taupe	1 ea
574G427	Blue Ridge	1 ea
574G432	Imperial Blue	1 ea
574G454	Royal Blue	1 ea
574G851	Doe Skin	1 ea
574G857	Jade	1 ea
574G863	Rosewood	1 ea

Fire Rated Colors

(Meets California Technical Bulletin 133
Flammability Standards)

FR574G401	Vintage Gold	1 ea
FR574G409	Warm Taupe	1 ea
FR574G427	Blue Ridge	1 ea
FR574G432	Imperial Blue	1 ea
FR574G454	Royal Blue	1 ea
FR574G851	Doe Skin	1 ea
FR574G857	Jade	1 ea
FR574G863	Rosewood	1 ea



- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.



© 2009 GF Health Products, Inc.

GRAHAM-FIELD

GF0900049RevA09