



# S-Ergo 125

## Ultra Lightweight Ergonomic Manual Wheelchair

- Medicare code K0004
- S-Style Ergonomic Seat
- High strength, weighs only 25 lbs. (w/o footrests)
- Silver 1/4" Anti-Bacterial cushion, removable for washing
- Ergonomic handrim
- Folding backrest/folding seat for easy travel
- Flip-Back armrest w/ wider concaved armpad
- Swing In & Away footrest allowing maximum safety while entering or exiting the chair
- "Tube-in Center" foot-plate, assures better side leg support
- Pouch for carrying small items attached to upholstery located behind backrest
- 24" Rear spoked Wheel w/polyurethane, high tread, flat free wheels
- 7x1" Polyurethane front casters
- 16"x17", 18"x17"
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black breathable mesh bottom, AEGIS top
- Weight Capacity: 250 lbs

### Features for S-Ergonomic Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slipage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/Concaved armpad, providing comfort & support



Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
S-125F16	16	17	8	19	17	36	23	12	39	25	250
S-125F18	18	17	8	19	17	36	25	12	39	25	250

*\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam Cushion, UV Resistant Canopy, Rigidifying Push-Bar, Foldable/Adjustable Headrest, Frog Leg Suspension, Magnesium Wheel, Vulcan Grip Wheel, Mag Wheel(Black Plastic)*

*\*S-115Q – Quick-release rear wheel model, same specifications as above.*



**KARMAN HEALTHCARE Inc.**