

## Seating Measurement Guide

**A** Seat Pan to Top of Shoulder: \_\_\_\_\_

**B** Chest Depth: \_\_\_\_\_

**C** Chest Width: \_\_\_\_\_

**D** Seat Depth: \_\_\_\_\_  
(Measurement from behind the knee)

**E** Seat Pan to Top of Head: \_\_\_\_\_

**F** Seat Pan to Elbow: \_\_\_\_\_

**G** Hip/Seat Width: \_\_\_\_\_


**H** Behind Knee to Heel: \_\_\_\_\_  
(Lower leg length, seat pan to footrest)

**I** Seat to Axilla/Armpit: \_\_\_\_\_

Age: \_\_\_\_\_

User Weight (kg/pounds) \_\_\_\_\_

User Height \_\_\_\_\_



Note: Please verify the measurements of the product, to ensure proper fitting for the client and adequate space availability in the user environment. medicaleshop will not be responsible for incorrect measurements/fitting and inadequate space availability in the user environment.

Need help with measurements, Please call us...

**Medicaleshop Inc.**

87 Danbury Road, Unit 1, New Milford, CT 06776, Phone : (866) 563-6812 | Fax : (860) 838-4671  
Email : [csteam@medicaleshop.com](mailto:csteam@medicaleshop.com)