



# Raz-AT600

US MSRP  
Effective July 1, 2017

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO # / Client Ref: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

## Model Selection

- Z360 **Raz-AT600** Heavy Duty Rehab Shower Commode Chair with stainless steel frame and components, 25° tilt range, four locking 6" (150mm) casters and commode pan, and the specifically engineered, Weight Assist Spring System (WASS).

Model Code	MSRP	Seat Type	Frame Width	Seat Code	"Seat Width (WW)"	"Seat Depth (DD)"	Seat Configuration (STC) and Access Opening				
<input type="checkbox"/> Z36022M	\$5,695	Molded	22"	<input type="checkbox"/> ZMS 23 18 FRO	23"	18"	<input type="checkbox"/> Front				
<input type="checkbox"/> Z36024M	\$5,920	Molded	24"	<input type="checkbox"/> ZMS 23 18 FRO	23"	18"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36022C	\$6,440	Custom	22"	<input type="checkbox"/> ZCS _ _ _ _	_	_	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36024C	\$6,665	Custom	24"	<input type="checkbox"/> ZCS _ _ _ _	_	_	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36026C	\$6,890	Custom	26"	<input type="checkbox"/> ZCS _ _ _ _	_	_	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36028C	\$7,115	Custom	28"	<input type="checkbox"/> ZCS _ _ _ _	_	_	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36030C	\$7,340	Custom	30"	<input type="checkbox"/> ZCS _ _ _ _	_	_	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear

### Notes:

- 1) For models with custom seat types, please specify seat width, depth and seat configuration to designate if the access opening is to be in the Front, Left, Right or Rear. Select Bridge if an access opening is not required. Use this information to create the seat code for the required seat.
- 2) Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for all custom seats.
- 3) Custom seat widths can be 1" or 2" wider than the frame width to provide additional thigh support without increasing overall width of the chair.
- 4) Custom seat depths must be between 16"D and 20"D.
- 5) The Ischial Pelvic Alignment System (IPAS) will allow 1" fore/aft adjustment of the aperture location and seat depth on custom seats. On the molded 23"W x 18"D seat, the IPAS will allow the aperture location to be adjusted up to 2" forward with a corresponding seat depth range from 18"D to 20"D.

## User Weight Range

<input type="checkbox"/> Z36WR450	WASS User Weight Range 275 - 450 lbs.	Standard	<input type="checkbox"/> Z36WR600	WASS User Weight Range 400 - 600 lbs	\$250
-----------------------------------	--	----------	-----------------------------------	---	-------

## Backrest Selection

<input type="checkbox"/> Z563T	Multi-Position Back with adjustable-tension vinyl upholstery 20" back height, with 1" fore/aft center of gravity adjustment, two back angles (90°/97°) and flip-up armrests. (Z710HD Not Included)	No Charge
<input type="checkbox"/> Z710HD	Locking Armrest Kit - HD	\$75
<input type="checkbox"/> ZCUSTBK	Custom Back	By Quotation
Notes: _____		

## Foot and Leg Support Options (Z709N Adjustable Footrest Included with Z360)

<input type="checkbox"/> Z709N	Adjustable Footrest HD 16" - 20.5" - (pair)	Standard	<input type="checkbox"/> Z7ELRL	Elevating Legrest - Left (includes recessed calf pad and Z7ADMFP)	\$345
<input type="checkbox"/> Z7FREXT	Footrest Extension (18" - 22.5" pair)	\$55	<input type="checkbox"/> Z7ELRR	Elevating Legrest - Right (includes recessed calf pad and Z7ADMFP)	\$345
<input type="checkbox"/> Z7ADMFP	Adjustable Angle/Depth Footplates - (pair)	\$55	<input type="checkbox"/> Z7ELROFL	Offset Elevating Legrest - Left (includes recessed calf pad and Z7ADMFP)	\$385
<input type="checkbox"/> Z708ADJ	Adjustable Heel Loops (pair)	\$95	<input type="checkbox"/> Z7ELROFR	Offset Elevating Legrest - Right (includes recessed calf pad and Z7ADMFP)	\$385
<input type="checkbox"/> Z7OFRL	Offset Footrest Receiver (moves footrest laterally 2") - Left	\$160	<input type="checkbox"/> ZCUSTFT	Custom Foot/Leg Support Construction	By Quotation
<input type="checkbox"/> Z7OFRR	Offset Footrest Receiver (moves footrest laterally 2") - Right	\$160	Notes: _____		

**Options & Accessories**

<input type="checkbox"/> Z678	Directional Caster Lock (for 6" casters)	\$145	<input type="checkbox"/> Z740AP2L	Lateral Support - with Large Depth-Adjustable Pad - Left*	\$320
<input type="checkbox"/> Z7HPL	Headrest Pad - Large (12" x 6")	\$135	<input type="checkbox"/> Z740AP2R	Lateral Support - with Large Depth-Adjustable Pad - Right*	\$320
<input type="checkbox"/> Z713	Chest Positioning Strap HD (1-piece with Velcro closure)	\$105	<input type="checkbox"/> Z744L	Lateral Mounting Hardware - Left (Required for Lateral Supports)	\$10
<input type="checkbox"/> Z714	Pelvic Positioning Strap HD (2-pieces with Fastex buckle)	\$105	<input type="checkbox"/> Z744R	Lateral Mounting Hardware - Right (Required for Lateral Supports)	\$10
<input type="checkbox"/> Z715L	Molded Flat Armpad - Left 14"L x 3.5"W	\$95	<input type="checkbox"/> Z746AP2L	Leg Adductor - with Large Depth-Adjustable Pad - Left*	\$320
<input type="checkbox"/> Z715R	Molded Flat Armpad - Right 14"L x 3.5"W	\$95	<input type="checkbox"/> Z746AP2R	Leg Adductor - with Large Depth-Adjustable Pad - Right*	\$320
<input type="checkbox"/> Z716L	Molded Arm Trough - Left	\$110	<input type="checkbox"/> Z736HDC	Chest Strap HD - Chest Bodypoint Aeromesh	\$205
<input type="checkbox"/> Z716R	Molded Arm Trough - Right	\$110	<input type="checkbox"/> Z736HDP	Pelvic Strap Strap HD Bodypoint Aeromesh	\$195
<input type="checkbox"/> Z717L	Pivoting Arm Mount - Left (for Z715L/Z716L)	\$165	<input type="checkbox"/> Z742	Lateral Adductor 1.5" Extension Plate Quantity _____	\$25 each
<input type="checkbox"/> Z717R	Pivoting Arm Mount - Right (for Z715R/Z716R)	\$165			
<input type="checkbox"/> ZSPCP1	Commode Pan Spacer Kit 1" Lower	\$200			
<input type="checkbox"/> ZSPCP2	Commode Pan Spacer Kit 2" Lower	\$215			
<input type="checkbox"/> ZSPCP3	Commode Pan Spacer Kit 3" Lower	\$230			
<input type="checkbox"/> ZSPCP4	Commode Pan Spacer Kit 4" Lower	\$245			
			Credit & Omissions		
			<input type="checkbox"/> Z905	Omit Commode Pan & Mounting Brackets	-\$20
			<input type="checkbox"/> Z908	Omit Adjustable Footrests (Z709N)	-\$45

**Comment Box**