


Safety Rail Heights – 3 Models

LOW – This is our classic *SleepSafe*® bed with safety rails ranging up to 17¼" above the mattress depending on the foundation choice.

MEDIUM – The *SleepSafe*® II provides protection ranging up to 27" above the mattress depending on the foundation choice.

TALL – The *SleepSafer*® provides protection up to 36¼" above the mattress depending on the foundation choice. The *SleepSafer*® **Extension Kit** offers an **additional 16¾" of safety rail height** when ordered or retrofitted to your existing *SleepSafer*® Bed.



| <i>SleepSafe</i> ® Bed Specifications | Low Bed <i>SleepSafe</i> ® | Medium Bed <i>SleepSafe</i> ® II | Tall Bed / Extension <i>SleepSafer</i> ® |
|---|---|--|--|
| SAFETY RAIL HEIGHT <i>Distance from the top of mattress to the top of the safety rail</i> | | | |
| Fixed Foundation with Box Spring (<i>Choice of settings</i>) | 12" or 8¾" | 17¾", 14" or 10" | 31" or 24¾" / 47¾" or 41½" |
| Fixed Foundation with Bunkie Board (<i>Choice of settings</i>) | 17¼" or 14" | 27" or 21¾" | 36¼" or 30" / 53" or 46¾" |
| Manual Articulating Foundation | 15" or 11¾" | 21", 17" or 13¼" | 34" or 27¾" / 50¾" or 44½" |
| Manual HI-LO Foundation | LO 11½" to EVEN with rail | LO 21¼" to HI 9½" or LO 16" to HI 4¼" | LO 30½" to HI 18¾" / LO 47½" to HI 35¾" |
| Articulating Foundation (<i>Choice of settings</i>) | 15" or 11¾" | 21", 17" or 13¼" | 34" or 27¾" / 50¾" or 44½" |
| Electric HI-LO Foundation (<i>LO to HI Settings</i>) | LO 11½" to HI 1¼" below rail | LO 21¼" to HI 11" or LO 16" to HI 5¾" | LO 30½" to HI 20¾" / LO 47½" to HI 37¼" |
| MATTRESS HEIGHT <i>Distance from the top of the mattress to the floor.</i> | | | |
| Fixed Foundation with Box Spring (<i>Choice of settings</i>) | 23¾" or 27" | 28", 32" or 35½" | 23¾" or 30" |
| Fixed Foundation with Bunkie Board (<i>Choice of settings</i>) | 18½" or 21¾" | 18¾" or 24" | 18½" or 24¾" |
| Manual Articulating Foundation | 20¾" or 24" | 25", 28¾" or 32½" | 20¾" or 27" |
| Manual HI-LO Foundation | LO 24¼" to HI 36" | LO 24½" to HI 36¼" or LO 29¾" to HI 41½" | LO 24¼" to HI 36" |
| Articulated Foundation (<i>Choice of settings</i>) | 20¾" or 24" | 25", 28¾" or 32½" | 20¾" or 27" |
| Electric HI-LO Foundation - (<i>LO to HI Settings</i>) | LO 24¼" to HI 34½" | LO 24½" to HI 34¾" or LO 29¾" to HI 40" | LO 24¼" to HI 34½" |
| FLOOR TO RAIL HEIGHT <i>Distance from the top of the rail to the floor.</i> | 35¾" | 45¾" | 54¾" / 72" |
| OVERALL OUTSIDE DIMENSIONS | | | |
| Twin | 42¼"W x 80½"L x 40¾"H | 42¼"W x 80½"L x 50¾"H | 42¼"W x 80¼"L x 56"H (72¾"H w/Ext.) |
| Full | 56"W x 80"L x 40¾"H | 56"W x 80"L x 50¾"H | 56"W x 80"L x 57"H (72¾"H w/Ext.) |
| Queen | 65"W x 80"L x 40¾"H | 65"W x 80"L x 50¾"H | 65"W x 80"L x 57"H (72¾"H w/Ext.) |
| MATTRESS SIZE <i>Width without pads / Width with pads / Length</i> | | | |
| Twin | 38"W / 36"W / 77"L | 38"W / 36"W / 77"L | 38"W / 36"W / 77"L |
| Full | 51"W / 50"W / 77"L | 51"W / 50"W / 77"L | 51"W / 50"W / 77"L |
| Queen | 61"W / 60"W / 77"L | 61"W / 60"W / 77"L | 61"W / 60"W / 77"L |
| UNDER BED FLOOR CLEARANCE | | | |
| Standard Models | 8" | 8" | 7¾" |
| TRANSFER HEIGHT <i>Distance from the top edge of the folded down safety rail to the floor.</i> | 22¼" with Safety Rail On 20¾" with Safety Rail Off | 32" with Safety Rail On 30½" with Safety Rail Off | 26" with Safety Rail On 24½" with Safety Rail Off |

