

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!
Please mark clearly with an X the items you wish to purchase.

Account #: _____	Ship to: _____
Date: _____	Name: _____
PO #: _____	Address: _____
Buyer: _____	City: _____
Attention: _____	Phone: () _____

Height Right™ Chair Features:



- Adjustable footrest for proper positioning
- Seat and footrest adjust allowing the chair to grow with the child
- Optional Soft-Touch® cushions easily attach for children who need lateral support
- Optional Early Pelvic Positioner, Padded Seat and Back Cushions are designed to fit the Height Right™ Chair.



Shown with Early Pelvic Positioner, and Wooden Tray with protective cover (sold separately)



Optional Accessories:







Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!
Please mark clearly with an **X** the items you wish to purchase.

Height Right™ Chair Specifications	
Seat Plate	16.5"W × 10"D
Seat Plate Height Range (from Floor)	16.5" - 24"
Foot Plate	11.5"W × 16.5"D
Foot Plate Height Range (from Floor)	4" - 15"
Tray	19"W × 8"D
Tray Height from Floor	26"
Overall Height	33"
Overall Width	17"
Overall Depth	21"
Weight Capacity	250 lbs.

Early Pelvic Positioner Specifications	
Seat Width	11"
Seat Depth	9.5"
Lateral Thigh Support Height	3.5"

Item	Part #	MSRP
Height Right™ Chair	61000000	\$ 419.25

Optional Accessories	Part #	MSRP
Wooden Tray (w/Protective Cover)	61010000	\$ 198.00
Plastic Tray Cover ONLY	61020000	\$ 103.63
Footcups (pair) - Dark Gray	61060009	\$ 313.75

Optional Cushions	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
Soft-Touch® Seat Cushion (HS)	61040001	61040002	61040004	61040009	\$ 209.11
Soft-Touch® Back Cushion (HB)	61050001	61050002	61050004	61050009	\$ 209.11
Early Pelvic Positioner (HP)	61070001	61070002	61070004	61070009	\$ 209.11

WARNING: This chair is as stable as a typical kitchen chair or baby highchair. However, if your child self-stimulates using vestibular movements (rocking back and forth), this product might not be the best choice to meet the child's needs.

WARNING: The Early Pelvic Positioner is ALWAYS required when using the Tray or Grab Rail for children UNDER the age of 3.