

# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 1 of 7

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
<b>TiLite Information:</b>	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:
More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or to view options.			



## FRAME STYLE

- |       |   |                |              |
|-------|---|----------------|--------------|
| G2FS1 | <input type="checkbox"/> <b>2GX Titanium Frame</b> (Full Adjustability. 265 Pound Weight Limit. Light Weight and High Strength. Not Available with Frame Transit Tie-Down.).....  | <b>\$4,692</b> | <b>K0005</b> |
| XFS1  | <input type="checkbox"/> <b>Swing Away</b> .....  | <b>STD</b>     |              |
| XFS3  | <input type="checkbox"/> <b>Hemi Height Frame</b> (Available with a Front Seat Height of 13.5" - 20" (34-51cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" Seat Depth or Heavy Duty Frame.).....                          | <b>NCO</b>     |              |
| XFS2  | <input type="checkbox"/> <b>Fixed Front</b> .....   | <b>NCO</b>     | <b>K0005</b> |
| XHDR1 | <input type="checkbox"/> <b>Heavy Duty Frame</b> (350 Pound Weight Limit. Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths. Not Available with 90° front end)..... | <b>NCO</b>     |              |
| XPWR1 | <input type="checkbox"/> <b>Power Adaptable Frame</b> (Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree Camber. Not available with Single-Sided Forks).....   | <b>NCO</b>     |              |
| XRFR1 | <input type="checkbox"/> <b>Reinforced Frame</b> .....  | <b>NCO</b>     |              |

## SMARTDRIVE POWER ASSIST SYSTEM



- |          |   |                |              |
|----------|---|----------------|--------------|
| XMAXSD1  | <input type="checkbox"/> <b>Permobil SmartDrive MX2+*</b> (Not available with rear seat heights less than 13.5")..... | <b>\$7,529</b> | <b>E0986</b> |
|          | XMAXSD18 <input type="checkbox"/> SwitchControl Buttons.....  | <b>NCO</b>     |              |
| XMAXSD8  | <input type="checkbox"/> SwitchControl Button with Mono Jack Option.....  | <b>NCO</b>     |              |
| XMAXSD19 | <input type="checkbox"/> <b>NEW</b> Buddy Button.....   | <b>NCO</b>     |              |
| XMAXSD4  | <input type="checkbox"/> Carrying Bag.....  | <b>\$84</b>    |              |
| XMAXSD10 | <input type="checkbox"/> Universal Padded Handle.....   | <b>\$22</b>    |              |

\*SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

## TRANSIT OPTIONS

- |      |  |              |              |
|------|--|--------------|--------------|
| XWTB | <input type="checkbox"/> <b>Transport Bracket Option</b> (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away -Flip Back Tubular Armrests. Requires a minimum seat depth of 14" or 15" seat Depth with Hemi Height Frame.)..... | <b>\$379</b> | <b>E1022</b> |
|------|--|--------------|--------------|

## FRAME FINISHES

### Natural Finishes

- |        |   |              |
|--------|---|--------------|
| XFF1   | <input type="checkbox"/> <b>Satin</b> (Available with Titanium Frame only.).....  | <b>STD</b>   |
| XSCBT1 | <input type="checkbox"/> <b>Satin Finish Buffing Pads</b> .....   | <b>\$22</b>  |
| XFF3   | <input type="checkbox"/> <b>Tattoos</b> (Tattoos on front frame tubes. Titanium frame finished in Satin. Aluminum frame tattooed in blasted aluminum with clear coat. Tattoos will be trimmed and adjusted to fit your frame.)..... | <b>\$363</b> |
|        | <input type="checkbox"/> <b>XFF45 Need for Speed</b>  |              |
|        | <input type="checkbox"/> <b>XFF44 Flower Garden</b>   |              |



### XFF43 Mosaic



- |  |   |   |
|--|---|---|
| <input type="checkbox"/> XFF8 High Gloss Black     | <input type="checkbox"/> XFF10 Red Pepper         | <input type="checkbox"/> XFF29 Neon Pink            |
| <input type="checkbox"/> XFF14 Ocean Blue Metallic | <input type="checkbox"/> XFF15 Tangerine Metallic | <input type="checkbox"/> XFF30 Acid Green Metallic  |
| <input type="checkbox"/> XFF17 Midnight Blue Pearl | <input type="checkbox"/> XFF18 White Pearl        | <input type="checkbox"/> XFF19 Black Lava           |
| <input type="checkbox"/> XFF27 Titanium Grey       | <input type="checkbox"/> XFF59 Pomegranate Red    | <input type="checkbox"/> XFF62 Blue Mist            |
| <input type="checkbox"/> XFF57 Indigo Purple       | <input type="checkbox"/> XFF40 Hyper Blue         | <input type="checkbox"/> XFF41 Flame Orange         |
| <input type="checkbox"/> XFF42 OD Green            | <input type="checkbox"/> XFF60 Silver Streak      | <input type="checkbox"/> XFF63 Lavender Dusk        |
| <input type="checkbox"/> XFF58 Black Galaxy        | <input type="checkbox"/> XFF61 Seafoam Green      | <input type="checkbox"/> XFF64 Daffodil Yellow..... |
|  |   | <b>STD</b>  |



# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 2 of 7

## FRAME MEASUREMENTS

Measurements shown below are our most popular selections. **Not seeing what you want?** Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.



**SEAT WIDTH A** Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**  
XSW2 ☐ 20" ☐ 21" ☐ 22" (21" and 22" Seat Widths A are only available with Heavy Duty Frame.).....**\$536 E2201**

**SEAT DEPTH B** Measure from 1/2" in front of back post to the front edge of seat sling. (12"-14" Seat Depths B not available with 25" (559) and 26"

(590) Rear Wheels)

XSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**  
XSD2 ☐ 20".....**\$568 E2203**  
XSD3 ☐ CustomFrame Depth: \_\_\_\_\_" Frame Depth Longer Than Seat Depth (Use 1" increments with a maximum of 2". A Custom Frame Depth with a maximum of 20" is \$546.).....**NCO**



**FRONT SEAT HEIGHT C** Measure from floor to top of seat tube at front edge of seat sling. Adjustable in 1/2" increments. (13.5" - 14" Front Seat Height C requires Hemi Height Frame Style. 20.5" - 21.5" Front Seat Height C not available with Hemi Height Frame Style.).

### FRONT SEAT HEIGHT RANGES

	Front Wheel Size	Front Seat Height C	Hemi Front Seat Height C
	3" Front Wheel	14.5" – 18.5"	13.5" – 17.5"
	3" Front Wheel w/ Glide™	16.5" – 17.5"	15.0" – 17.0"
	4" Front Wheel	15.5" – 19.0"	14.5" – 18.0"
	4" Front Wheel w/ Glide™	17.0" – 18.0"	15.5" – 17.0"
	5" Front Wheel	16.5" – 19.5"	16.0" – 18.5"
	5" Front Wheel w/ Glide™	17.5" – 18.5"	16.5" – 18.0"
	6" Front Wheel	17.5" – 20.0"	16.5" – 19.0"
	6" Front Wheel w/ Glide™	19.0"	17.5" – 18.0"
	8" Front Wheel	19.5" – 21.0"	18.5" – 19.5"
XFSH1	<input type="checkbox"/> 13.5" <input type="checkbox"/> 14" <input type="checkbox"/> 14.5"	<input type="checkbox"/> 15" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16"	<input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0"
	<input type="checkbox"/> 17.5" <input type="checkbox"/> 18.0" <input type="checkbox"/> 18.5"	<input type="checkbox"/> 19.0" <input type="checkbox"/> 19.5" <input type="checkbox"/> 20.0"	<input type="checkbox"/> 20.5" <input type="checkbox"/> 21.0"..... <b>STD</b>

**REAR SEAT HEIGHT D** Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

Rear Wheel Size	Rear Seat Height
22"	13.0" – 18.0"
24"	14.0" – 19.0"
25"	14.0" – 19.0"
26"	15.0" – 20.0"

XRSH1 ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5"  
☐ 19" ☐ 19.5" ☐ 20".....**STD**

## BACK SUPPORTS

### COMFORT ACTA-BACK® Back Support

☐ Acta-Back with 2-PT Quick Release Compass Hardware.....**\$616 E2613**  
**Shell Height** ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" .....**STD**  
☐ 20".....**\$83**  
☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....**NCO**  
☐ Color Accent (Colored piping along the edge of the cover on the backside of the backrest) ☐ Silver.....**STD**  
☐ Red ☐ Purple ☐ Pink ☐ Blue ☐ Green ☐ Orange ☐ Remove (no color).....**\$30**  
☐ Comfort Privacy Shield.....**\$44**

### COMFORT ACTA-BACK® Deep Back Support

☐ Acta-Back Deep with 2-PT Quick Release Compass Hardware.....**\$730 E2615**  
**Shell Height** ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" .....**STD**  
☐ 20".....**\$93**

# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 3 of 7

<input type="checkbox"/> Comfort-Tek Cover	<input type="checkbox"/> Stretch-Air Cover	NCO					
<input type="checkbox"/> 3" Deep Wing	<input type="checkbox"/> 5" Deep Wing	NCO					
<input type="checkbox"/> Color Accent (Colored piping along the edge of the cover on the backside of the backrest)	<input type="checkbox"/> Silver	STD					
<input type="checkbox"/> Red	<input type="checkbox"/> Purple	<input type="checkbox"/> Pink	<input type="checkbox"/> Blue	<input type="checkbox"/> Green	<input type="checkbox"/> Orange	<input type="checkbox"/> Remove (no color)	\$30
<input type="checkbox"/> Comfort Privacy Shield							\$44

## COMFORT ACTA-EMBRACE® Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS16	<input type="checkbox"/> Acta-Embrace Back with 2-PT Quick Release Compass Hardware	\$642	E2613
CPSB1	<input type="checkbox"/> Comfort Privacy Shield	\$44	

## COMFORT ACTA-RELIEF Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS13	<input type="checkbox"/> Acta-Relief Back with 2-PT Quick Release Compass Hardware	\$846	E2615
CPSB1	<input type="checkbox"/> Comfort Privacy Shield	\$44	

**SEAT BACK HEIGHT E/SEAT BACK/PUSH HANDLES** Measure from top of back post to top of seat tube at rear of frame.

### SEAT BACK TYPE

XSB13	<input type="checkbox"/> Fixed Angle, Adjustable-Height Seat Back	STD
XLB1	<input type="checkbox"/> Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected. Includes Tension Adjustable by Strap Back Upholstery.)	\$423
XSB14	<input type="checkbox"/> Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.)	NCO
XSB18	<input type="checkbox"/> Folding Adjustable-Angle and Height Seat Back (Not available with Heavy Duty Frame.)	\$250
XSB16	<input type="checkbox"/> Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" Seat Depths B, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths B of 12"-17" allow 3" of rearward adjustment. Seat Depth B of 18" allow 2" of rearward adjustment. Requires Tension Adjustable Seat Upholstery)	\$244

### PUSH HANDLES

XSB6	<input type="checkbox"/> Integrated Push Handles	NCO
XSB7	<input type="checkbox"/> Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame.)	\$125
XSB9	<input type="checkbox"/> Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests or Swing Away - Flip Back Tubular Armrests.)	\$363

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT E** (11"- 15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Height Adjustable Push Handles or Folding Stabilizer Bar.)

Push Handles (Selected on Pg. 2)	Short	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Integrated	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Bolt-On Height Adj.	-----	<input type="checkbox"/> 16" - 20"	Set at: _____

XSB18 **FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK E** (Not available with a Folding Stabilizer Bar or Heavy Duty Frame)

Push Handles (Selected on Pg. 2)	Short	Medium	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/> 11.0" - 13.0"	<input type="checkbox"/> 13.5" - 16.5"	<input type="checkbox"/> 16.5" - 19.5"	Set at: _____
Integrated	<input type="checkbox"/> 12.5" - 14.5"	<input type="checkbox"/> 13.5" - 16.5"	<input type="checkbox"/> 17.0" - 19.5"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 15.0" - 17.0"	<input type="checkbox"/> 16.0" - 19.0"	<input type="checkbox"/> 18.5" - 20.5"	Set at: _____

## FRAME MEASUREMENTS

Options shown below are our most popular selections. **Not seeing what you want?** Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.

### CAMBER F

XCBR8	<input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 4° <input type="checkbox"/> 6°	STD
-------	---	-----

**FRONT ANGLE H** Front Angle H is measured from inside of front frame tube to top of bottom frame tube.

### Swing Away Front End

XFA3	<input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 85° (Not available with 6" Front Wheels.) <input type="checkbox"/> 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.)	STD
XFA0	<input type="checkbox"/> Omit Hangers	NCO

### Fixed Front End

XFA6	<input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 85° (85° Not available with 6" Front Wheels.)	STD
------	--	-----

**SEAT TO FOOTREST I** Measure from front edge of seat sling to top rear of footrest.

**If your Front Angle H is 70°, 80°, OR 85° select one Seat to Footrest I length and one Footrest type:**

**Standard Footrests** Measurement must be at least 2" shorter than Front Seat Height C. 70° Fixed Front has a minimum seat to footrest of 14.5"

☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 4 of 7

XFTR5	<input type="checkbox"/> Flip-up Footrest.....	STD	
XFTR9	<input type="checkbox"/> Depth Adjustable Flip-up Footrest.....	\$120	K0045
XFTR10	<input type="checkbox"/> Depth & Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest I only.).....	\$293	K0040, K0045
XFTR11	<input type="checkbox"/> Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest I only.).....	\$185	K0040, K0045

## High Mount Footrests (Not Available with Luggage Carriers, or Impact Guards.)

	<input type="checkbox"/> 7" <input type="checkbox"/> 7.5" <input type="checkbox"/> 8" <input type="checkbox"/> 8.5" <input type="checkbox"/> 9" <input type="checkbox"/> 9.5" <input type="checkbox"/> 10" <input type="checkbox"/> 10.5" <input type="checkbox"/> 11" <input type="checkbox"/> 11.5" <input type="checkbox"/> 12" <input type="checkbox"/> 12.5" <input type="checkbox"/> 13" <input type="checkbox"/> 13.5" <input type="checkbox"/> 14" <input type="checkbox"/> 14.5"		
XFTR20	<input type="checkbox"/> High Mount Flip-up.....	\$202	K0037
XFTR15	<input type="checkbox"/> Angle Adjustable High Mount Flip-up Footrest.....	\$352	K0040, K0037
XFTR16	<input type="checkbox"/> Depth & Angle Adjustable High Mount Flip-up Footrest.....	\$406	K0040, K0037

If your Front Angle H is 90°, select one Seat to Footrest I length and one Footrest type:

	<input type="checkbox"/> 7" <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13"		
XFTR5	<input type="checkbox"/> Flip-up Footrest.....	STD	
XFTR9	<input type="checkbox"/> Depth Adjustable Flip-up Footrest.....	\$120	K0037
XFTR10	<input type="checkbox"/> Depth & Angle Adjustable Flip-up Footrest.....	\$293	K0040, K0037
XHLP1	<input type="checkbox"/> Heel Loops (Not available with 13" or narrower Seat Width A or 1-Piece Footrests.).....	\$55	E0951

## Elevating Leg Rest

	<input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21"		
XFA5	<input type="checkbox"/> Single Depth Adjustable Flip-up Footrest - <input type="checkbox"/> Left OR <input type="checkbox"/> Right.....	\$185	E0990
XFTR17	<input type="checkbox"/> Depth Adjustable Flip-up Footrest - Pair.....	\$363	

## CUSTOM CONFIGURATION OPTIONS

### FRONT WHEELS (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height C limitations.)

XFW24	<input type="checkbox"/> 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$206	
XFW14	<input type="checkbox"/> 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame or Heavy Duty Frame).....	NCO	
XFW3	<input type="checkbox"/> 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame or Heavy Duty Frame.).....	\$88	
XFW26	<input type="checkbox"/> 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.).....	\$112	
XFW27	<input type="checkbox"/> 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$217	
XFW6	<input type="checkbox"/> 5"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or with 90° Front Angle H.).....	NCO	
XFW28	<input type="checkbox"/> 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 90° Front Angle H.).....	\$103	
XFW30	<input type="checkbox"/> 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 90° Front Angle H.).....	\$233	
XFW9	<input type="checkbox"/> 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or 90° Front Angle H.).....	NCO	
XFW32	<input type="checkbox"/> 6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame, 85° or 90° Front Angle H.).....	\$119	
XFW33	<input type="checkbox"/> 6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle H.).....	\$244	
XFW16	<input type="checkbox"/> 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or 12"-14" Seat Depths. Available with 70° Front Angle H only.).....	NCO	

### FRONT FORKS

XFK1	<input type="checkbox"/> TiLite Standard Fork.....	STD	
XGL1	<input type="checkbox"/> Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame).....	\$526	E1015
	Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		

### ONE ARM DRIVE

XOAD1	<input type="checkbox"/> One Arm Drive (Available with 24" Mag Wheels only. Rear Seat Height, Camber or axle restrictions may apply.).....	\$1298	E0958
	Drive Handrim Side (Available with standard aluminum handrim only.) <input type="checkbox"/> Install on Left <input type="checkbox"/> Install on Right		

### STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height D limitations.)

XRW14	TiLite SHADOW™.....	STD	
	<input type="checkbox"/> 22"(501)		
	<input type="checkbox"/> 24"(540)		
	<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth B.)		
	<input type="checkbox"/> 26"(590) (Requires 15" minimum Seat Depth B.)		

### PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height D limitations.)

XRW7	Spinergy® SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Green.....	\$915	
	<input type="checkbox"/> 22"(501) - Std. Hub		
	<input type="checkbox"/> 24"(540) - Std. Hub		
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)		
	<input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)		
XRW8	Spinergy® LX (Available with Heavy Duty Frame only when 0 or 2 degree camber is selected.)		
	<input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple <input type="checkbox"/> Green.....	\$1148	
	<input type="checkbox"/> 22"(501) - Std. Hub		
	<input type="checkbox"/> 24"(540) - Std. Hub		
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)		
	<input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)		

# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 5 of 7

ZRW9 Spinergy® Design Lab LX Wheels (Available with black spokes only.)  
☐ Gold Anodize ☐ Green Anodize ☐ Blue Anodize ☐ Red Anodize ☐ Purple Anodize ☐ Orange Anodize.....\$1,364  
☐ 22"(501) - Std. Hub  
☐ 24"(540) - Std. Hub  
☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)  
☐ 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)

XRW16 Spinergy® Carbon Blade (Not available with Heavy Duty Frame.).....\$1,602  
☐ 24"(540)  
☐ 25"(559) (Requires 15" minimum Seat Depth B.)

XRW17 X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Sliver Anodized Handrims.).....NCO  
☐ 22"(501)  
☐ 24"(540)

## AXLES

XAXL1 ☐ Stainless Quick Release.....STD  
XAXL2 ☐ Titanium Quick Release.....\$277  
XAXL3 ☐ Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.).....\$169

## PNEUMATIC REAR TIRES

XRTR1 ☐ Treaded 22"(501), 24"(540), and 26"(590) (Not available with 25"(559) Rear Wheels.).....STD  
XRTR2 ☐ Airless Inserts.....\$98  
XRTR31 ☐ New MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559) and 26"(590).....STD  
XRTR30 ☐ 24"(540), 25"(559) and 26"(590) Primo Sentinel.....\$136  
XRTR8 ☐ 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....\$152  
XRTR32 ☐ 24"(540) Primo C-763 Knobby.....\$104

E2213

## SOLID REAR TIRES

XRTR22 ☐ 22"(501), 24"(540) and 25"(559) Primo Orion (Not available with Spinergy Carbon Blade Rear Wheels.).....\$120  
XRTR23 ☐ 24"(540) and 25"(559) Primo Active.....\$147

## HANDRIMS

XRIM1 ☐ Aluminum – Silver Anodized.....STD  
XRIM21 ☐ Aluminum – Black Anodized.....NCO  
XRIM12 ☐ Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....\$509  
XRIM10 ☐ The Surge (All-In-One Oval with Gription™ Strip, Larger Profile).....\$433  
XRIM11 ☐ The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile).....\$466  
XRIM20 ☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating).....\$358

**The Natural-Fit® (Original, Larger Profile)** (Not available with 20"(451) Rear Wheels.)

TRIM14 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$428  
TRIM15 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$466  
TRIM16 ☐ Without Thumb Grip.....\$390

**The Natural-Fit® LT (Reduced Weight, Smaller Profile)** (Not available with 20"(451) Rear Wheels.)

TRIM17 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$440  
TRIM18 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$514  
TRIM19 ☐ Without Thumb Grip.....\$423

## Additional Handrim Options

XRIM5 ☐ Plastic Coated – Black.....\$174  
XRIM22 ☐ Omit Handrims.....NCO

## Tab Length

XRIM7 ☐ Long Tabs.....STD  
XRIM8 ☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....\$28

## WHEEL LOCKS (Composite wheel locks not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

XWLK7 ☐ Compact Composite Push to Lock.....NCO  
XWLK8 ☐ Compact Composite Pull to Lock.....NCO  
XWLK10 ☐ Extension Handles for Compact Composite Push/Pull to Lock.....\$93

E0961

## Additional Wheel Locks

XWLK11 ☐ MBL Composite Scissor Lock (Not available on 2GX Swing Away.).....\$158  
XWLK3 Uni-Lock Set At: ☐ Push **OR** ☐ Pull.....\$179  
XWLK4 ☐ Extension Handles for Uni-Lock.....\$93  
XWLK15 ☐ Grade Aids Push to Lock (Not available with solid rear tires).....\$190

E0961

E0974



# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 6 of 7

## BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XBUP1 ☐ Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.)..... **STD**

Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)

XBUP3 ☐ Padded Nylon..... **\$369 E2611**  
XBUP5 ☐ AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.)..... **\$369 E2611**  
XUPC3 Ultrasuede Backrest Stripe ☐ Black (Not available with Padded Nylon or AIR Breathable Back Upholstery.)..... **\$115**  
XBUP0 ☐ Omit Back Upholstery..... **NCO**

## SEAT UPHOLSTERY

XSUP4 ☐ Tension Adjustable (Black only)..... **STD**  
XSUP0 ☐ Omit Seat Upholstery..... **NCO**

## SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XCSH2 ☐ Comfort Curve Cushion (Stretch-Air cover)..... **\$122 E2601**  
ZCSH23 ☐ Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion..... **\$456 E2607**  
XCSH15 ☐ ROHO® Hybrid Elite SR™ Single Compartment Cushion (Fits seat widths 14" to 19") ..... **\$552 E2622**  
XCSH14 ☐ ROHO® QUADTRO SELECT® Cushion (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High..... **\$594 E2624**  
XCSH16 ☐ ROHO® Hybrid Select Cushion (Fits seat widths 14" to 18", and 20")..... **\$571 E2624**  
XCSH17 ☐ Vicair Active O2 Low (Fits even seat widths A 14" to 20". Non-stock upcharges may apply) ..... **\$515 E2624**

## CUSTOM CONFIGURATION OPTIONS

### SIDE GUARDS

#### Aluminum

XSDG2 ☐ Rigid Removable ☐ Small **OR** ☐ Std..... **\$217**  
XSDG3 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std..... **\$233**  
XSDG6 ☐ Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.)..... **\$363**

#### 2x2 Twill Weave Carbon Fiber

XSDG7 ☐ Rigid Removable ☐ Small **OR** ☐ Std..... **\$336**  
XSDG8 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std..... **\$347**

☐ **Custom Offset Profiled Sideguards** CMD 007. Please specify desired offset to nearest  $\pm\frac{1}{2}$ " (e.g.,  $\pm 0.50$ ") and sideguard type (rigid or fold-down) in the Notes section of the order form. Sideguard material must be selected. Some limitations may apply.

## ARMRESTS (Desk Arms come with built-in Side Guards.)

XARM1 Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.)..... **\$244 E0973**  
☐ 9" - 12" **OR** ☐ 12" - 15"  
XARM4 ☐ Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.)..... **\$342 E0973**  
XARM2 Removable Height Adjustable DeskArm with Rigid Side Guard (Requires 14" minimum Seat Depth B.)..... **\$363 E0973**  
☐ 9" to 12" **OR** ☐ 12" to 15"  
☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")  
XARM3 Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth B.)..... **\$406 E0973**  
☐ 9" to 12" **OR** ☐ 12" to 15"  
☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

### Available Armrest Option

XARM5 ☐ Covered in Black Ultrasuede (Not available with Desk Arms)..... **\$120**

## SEAT BELTS

XBLT1 ☐ Velcro Adjustable, Black Nylon Webbing..... **\$61 E0978**  
XBLT2 ☐ Bodypoint® 2" Push-Button Buckle (Auto-Style.)..... **\$103 E0978**  
XBLT4 ☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... **\$196 E0978**  
XBLT5 ☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... **\$185 E0978**  
XBL10 ☐ Bodypoint® Optional Band Clamp Mounting Hardware..... **\$22 K0108**  
XBLT16 ☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button..... **\$201 E0978**

## CALF STRAP

XCLF2 Velcro Adjustable..... **\$44 K0038**  
☐ Small (Recommended for 9"- 12" Footrest Widths)  
☐ Medium (Recommended for 11"- 15" Footrest Widths)  
☐ Large (Recommended for 14"- 17" Footrest Widths)  
☐ Xtra Large (Recommended for 16"- 19" Footrest Widths)  
XCLF3 Bodypoint® Padded Velcro Adjustable..... **\$109 K0038**  
☐ Medium (Recommended for 12"- 15" Footrest Widths)  
☐ Large (Recommended for 16"- 20" Footrest Widths)  
XCLF1 ☐ Buckle Adjustable, Black Nylon Webbing..... **\$34 K0038**

# TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 7 of 7

## RESIDUAL LIMB SUPPORT

XRSL1	<input type="checkbox"/> Residual Limb Support (8" width x 10" length).....	\$428	E1020
XRSL2	<input type="checkbox"/> Residual Limb Support Pair (8" width x 10" length).....	\$855	E1020

## BACKPACK

XBCK3	<input type="checkbox"/> Permobil Medical Necessities Backpack.....	\$103	
-------	---	-------	--

## SEAT POUCH

XPCH3	<input type="checkbox"/> Permobil Medical Necessities Seat Pouch.....	\$66	
XPCH4	<input type="checkbox"/> Bodypoint® Mobility Bag.....	\$109	

## CUSTOM CONFIGURATION OPTIONS

### SPOKE GUARDS

XGRD1	<input type="checkbox"/> Clear (Not available with Mag Rear Wheels).....	\$147	K0065
-------	--	-------	-------

### ANTI-TIPS

XTIP1	<input type="checkbox"/> Rear Aluminum - Black.....	\$174	E0971
XTIP2	<input type="checkbox"/> User Friendly Flip-up.....	\$363	E0971
XTIP9	<input type="checkbox"/> Tip Assist.....	\$115	E0971

### ACCESSORIES

XUEST1	Upper Extremity Support Tray <input type="checkbox"/> 12" to 14" OR <input type="checkbox"/> 13" to 15" Seat Widths A (Requires Armrests.).....	\$315	E0950
XLUG1	<input type="checkbox"/> Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.).....	\$158	
XMPCT2	Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.) <input type="checkbox"/> Black.....	\$115	
XCH1	<input type="checkbox"/> Crutch Holder.....	\$244	
XAFP1	4" Logo Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines.....	\$115	

## GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

**Get the latest order forms at [permobilus.com](https://permobilus.com).** Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Strain® is a registered trademark of Q'Strain Inc.

## NOTES

### QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

### ORDER THIS CHAIR