

Hygiene Sling Clips Ploy

M921529

- Easy to apply padded sling
- Built in adjustable head support
- Good access for hygiene purposes
- Good upper body tone required
- Full risk assessment required for suitability
- Built in waist support with clunk-click fastening
- Universal clip fixings
- Wash temp 95°C / 203°F
- Tumble dry low heat (max 50°C)



Product Codes and Specifications:



100% Polyester		A	B	C	D	E
XS	M921529-XS	30¼ - 34¼" (77 - 87cm)	8¾" (22cm) or over	24" (61cm)	9½" (24cm)	24¾" (63cm)
S	M921529-S	34¼ - 42¼" (87 - 107cm)	10" (25cm) or over	28" (71cm)	11" (28cm)	27½" (70cm)
M	M921529-M	38¼ - 46" (97 - 117cm)	10¾" (27cm) or over	31¼" (79cm)	11½" (29cm)	30" (76cm)
L	M921529-L	44 - 52" (112 - 132cm)	11" (28cm) or over	34¾" (88cm)	12½" (32cm)	33½" (85cm)
XL	M921529-XL	49¼ - 57" (125 - 145cm)	11¾" (30cm) or over	38¼" (97cm)	13½" (34cm)	34¼" (87cm)
XXL	M921529-XXL	54¼ - 62¼" (138 - 158cm)	13½" (34cm) or over	39¾" (101cm)	13¾" (35cm)	37" (94cm)

