



# HTS (Hygiene & Toileting System)

**The Rifton HTS** is designed to improve the health of the client and the life of the caregiver.



Scan here to order, request a quote, or see detailed product information.

- Hygienic, simple and affordable
- Promotes ideal positioning for effective toileting
- Use it for toileting, showering, and bathing





# Key features

## 1 Unmatched versatility

The HTS works on, over, or off the toilet—adapting seamlessly to any setting, any toilet style, and a wide range of special needs.



**On the toilet**

**Over the toilet**

**Off the toilet**

The HTS can be securely mounted to any toilet with the appropriate mounting bar, with no tools required. Alternatively, use the mobile base to roll it over a toilet, or add a commode pan for use anywhere outside the bathroom.

## 2 Therapeutic positioning

The HTS provides the safety, security, and support children with special needs rely on, making effective, stress-free toileting possible.



## The HTS in use



Here the HTS is securely mounted on the toilet. Mounting bars allow quick, tool-free attachment of the chair to the toilet base.



The mobile base and commode pan allow the HTS to be used anywhere.

### 3 Convenient hygiene

The open seat pad design simplifies hygiene tasks, making cleaning easier for caregivers while the client remains comfortable and relaxed.



### 4 Room to grow

The HTS accommodates a wide range of sizes, body types, and needs. Adjustable frame height, seat depth, and footrest ensure a comfortable fit as children grow.



The small and medium HTS share the same frame, with seat height adjustable from 15" to 21½".

Seat depth adjusts from 8" on the small HTS to 17½" on the large model.



The natural contours of the chair and cushioned molded foam pads provide both support and comfort, while soft straps protect sensitive skin.



The HTS adjusts to fit a wide range of ages and sizes, providing years of use. Shown here: the same chair adapted for a five-year-old and an eleven-year-old.



# Key features

## 5 More than a toileting chair

The HTS offers comfort, dignity, and caregiver convenience in every toileting and hygiene task.



The small and medium HTS seats fit the tub stand for safe in-tub care, while the tilt-in-space base allows easy showering and hair washing.

## 6 Portable

The portability kit delivers the hygiene and support essential for successful toileting, giving clients and families the freedom to travel with confidence.



**Video:** See how the HTS portability kit sets up over a toilet.



## The HTS in use



The compact footprint of the tub stand fits both narrow and wide tubs.



The portability base includes a lightweight folding frame and carry bag, making toileting possible while traveling.



Scan for sizing  
and pricing.

## 7 Easier lifts and transfers

The flip-up footboard and gas-spring–assisted forward-tilting seat make transitions easier. Swing-away laterals and removable armrests provide added convenience for safe, efficient transfers.



## 8 Hygienic and easy to clean

The HTS features waterproof construction and quick disassembly, allowing every part to be washed and disinfected with a simple 10% bleach solution.



Lifts and transfers present a major challenge to successful toileting. The HTS makes them easier with a variety of features.



The HTS adapts to every stage of growth, offering a custom fit for clients of all sizes and needs.



# HTS options

## Frame size options



**Small Z110**

**Medium Z120**

**Large Z130**

## The HTS in use

---



# HTS components



Scan for sizing  
and pricing.

## Base options



**Mobile**  
(tilt-in-space with footboard)



**Mobile**  
(non-tilt with footboard)



**Stationary**  
(non-tilt with footboard)

**Note:** check clearance if using over a toilet.



**Tub base**  
(non-tilt, for medium and small HTS seats)

Use the **mounting bar** to attach the HTS seat to round-front or elongated toilets. The easy on/off design lets caregivers quickly position or remove the seat.



### Mounting bar for round front toilet

For round-front toilets with a distance of 17" or less from the seat bolt centerline to the front edge of the bowl.



### Mounting bar for elongated toilet

(compatible with small and medium HTS) For elongated toilets with a distance greater than 17" from the seat bolt centerline to the front edge of the bowl.

## Components and options



**Rear-closed**

**Rear-open**

**No-hole**



**Soft seat pad**

### Seat pads

Waterproof foam pads remove easily for cleaning. Choose a rear-open pad for hygiene care while seated, or rear-closed and no-hole pads for added support. A rigid hole filler is available for medium and large models. The soft seat pad offers extra cushioning for sensitive skin.



### Armrests

Adjustable armrests are available in three heights and remove easily for transfers.



# More HTS components



## **Push handles**

Allow caregivers to easily move a client seated in the HTS.



## **Headrest**

Provides support for limited head control; adjusts in height, depth, and angle.



## **Butterfly harness**

Offers maximum trunk support, security, and freedom of movement.



## **Anterior support/tray**

Encourages forward-leaning posture to support toileting.



## **Laterals and chest strap**

Provide trunk support for clients with limited torso control.



## **Hip guides**

Reduce seat width, align the torso, and improve trunk control; removable for transfers.



## **Large hip guides** (thick pads)

Thicker pads allow narrower seat width adjustment.



## **Pan**

Durable off-toilet pan with pour spout for easy cleaning.



## **Bowl adapter**

Fills the front gap between the HTS seat and the toilet, ensuring waste drops cleanly into the bowl. It is especially useful on wall-hung and round-front toilets with the large HTS.



Standard

Splash guard

Dynamic

### Deflectors

Choose from standard, splash guard, or dynamic deflector (with abductor) for dependable containment.



### Abductor

Separates knees, supports leg muscles, and promotes healthy hip and joint alignment.



### Calf rest and ankle straps

Enhance comfort, prevent sliding, and provide safe, moderate movement control.



### Footboard for use with mounting bar

Adjusts in height and angle to support the client's feet.



### Large casters

Improve mobility by clearing higher thresholds and obstacles.



### Portability base with carry case

Lightweight folding base and carry bag make toileting possible while traveling.



### Tub base

Tilts the seat back 5° for showering and hair washing; fits narrow and wide tubs.

## Need support?

The Rifton website offers everything you need:

- Evidence-based articles
- Sample letters of medical necessity
- Video tutorials and training resources
- Product instruction guides
- Clinical resources and best practices



Scan here to order, request a quote, or see detailed product information.