

HOGGI

MADE IN GERMANY

BINGO

EVOLUTION

Tilt in Space Wheelchair

RUGGED &
ADAPTABLE
FOR ACTIVE
FAMILIES

Sizes 1-2 • 2 XL • 3
Up to 110 lbs



Reverse Seating



Transit



Ultra Compact



Cobra Base



Indoor
Hi-Low
Full Tilt/
Recline

Coupling Kit



Reclining Backrest



45° Tilt



“Elevating Legrest”



“Growing Seat”




distributors of

HOGGI®

BINGO SWITCH

Innovative quick release connector to change 4-wheel mobility base to 3-wheel mobility base with center steer wheel attachment



◀ Hoggi Switch Video



Scan for product details

- Compatible with all current seating options, sizes and components including IUS Seat, OT Plus Seat, Respiratory Base Modification, and Tandem Coupling Unit.
- Transit approved in 4-wheel mobility configuration

Upholstery Color Selections for Bingo Evolution and Bingo Switch



Berry Soft



Black



Turquoise Soft



Gray Soft



Reed Green



Headrest Pads



Highrise Headrest Pads



Anatomic Head Support



Trunk Supports



Removable Abduction Block



Chest-shoulder Harness



5-point Harness



Lap Belt



4-point Lap Belt



Groin Strap



Hip Supports
(required component)



Armrests



Therapy Tray



Grab Rail With Upholstery



Canopy with Rain Cover



Foot Straps



Ankle Cuffs



Frame Riser



3 Piece Removable Vent Tray



Oxygen Bottle Holder

Bingo Evolution OT Plus Frame for After Market Seating



The OT Plus seat shell has all the same Tilt and Recline capabilities but allows the attachment of after market seating of your choice to the 1" round tubing.

(Seating depicted to the right is for example purposes only and purchased by dealer from a separate manufacturer)



| | Size 1 | Size 2 | Size 2XL | Size 3 |
|---|---|-------------------|-------------------|-------------------|
| MOBILITY BASE | | | | |
| Push Handle Height | 30" - 43" | 30" - 45" | 34" - 49" | 34" - 49" |
| Overall Width/Length *Switch adapter adds 10" of length | 25"/36" | 27"/36" | 27"/36" | 29"/38" |
| Wheel Diameter Front/Rear | 7.5"/11" | 7.5"/11" | 7.5"/11" | 7.5"/11" |
| Tilt-in-Space | 45° | 45° | 45° | 45° |
| Folded size w/Wheels | 32" x 25" x 15" | 32" x 27" x 15" | 32" x 27" x 15" | 32" x 29" x 17" |
| Folded Size w/ Front Center Wheel | 40.5" x 24" x 20" | 40.5" x 26" x 20" | 40.5" x 26" x 20" | 40.5" x 28" x 22" |
| Base only Weight | 19 lbs | 19 lbs | 19 lbs | 20 lbs |
| Rider Weight Limit | 110 lbs | 110 lbs | 110 lbs | 110 lbs |
| IUS SEATING SYSTEM | | | | |
| Seat Depth | 6" - 13" | 9" - 16" | 11" - 18" | 11" - 18" |
| Seat Width | 7" - 12" | 9" - 14" | 9" - 14" | 14" - 17" |
| Lower Leg Length* (without optional risers) | 8" - 15" | 8" - 16" | 8" - 16" | 8" - 18" |
| Backrest Height | 18" - 25" | 21" - 30" | 24" - 33" | 24" - 33" |
| Height of Slots for Shoulder Straps | 11" - 18" | 13" - 22" | 17" - 25" | 17" - 25" |
| Recline | 90° - 160° *with recline, backrest reaches 180° | | | 90° - 130° |

| | Size 1 | Size 2 | Size 3 |
|--|---|----------|------------|
| OT PLUS REHAB FRAME | | | |
| *Separate Back, Seat and accessories supplied by dealer. | | | |
| Seat Depth | 13.25" | 15.25" | 17.5" |
| Seat Width | 13.25" | 15.25" | 17.25" |
| * When ordering afermarket seating, max solid base width | 11.75" | 13.75" | 16.75" |
| Back Cane Tubing Diameter | 1" | 1" | 1" |
| Back Cane Height (Pan to top of cane) | 15.5" | 17.5" | 23.5" |
| Back Cane Width (Outside to outside) | 12" | 14" | 16" |
| Lower Leg Length* (without optional risers) | 8" - 15" | 8" - 16" | 8" - 18" |
| Recline | 90° - 160° *with recline, backrest reaches 180° | | 90° - 130° |
| Seat only Weight | 14 lbs | 16 lbs | 17 lbs |

| | Size 1 | Size 2 | Size 2XL | Size 3 |
|----------------------------|------------|----------|----------|----------|
| FOOTPLATE DIMENSION | | | | |
| External | 12" x 7.5" | 14" x 9" | 14" x 9" | 14" x 9" |
| Internal | 11.5"x7" | 13" x 8" | 13" x 8" | 13" x 8" |

Transit Brackets Crash tested according to ISO 7176-19 ANSU/RESNA WC / VOL. 1 - Section 19



distributors of

