



Ascent XR™ Euro-Style Rollator

The Ascent XR is a unique extra-rugged rollator that accommodates a larger range of users. It is featured in two colors: Cobalt Blue and Brandywine. The two major features of the Ascent XR are a larger weight capacity of 450 LBS and a rear shock absorbing suspension system.



- Easy push button for handle height adjustments in 1" increments



- Innovative design hides brake lines for a clean look while protecting from wear and breakage



- 10-inch Front Caster Reduces Resistance for Easy Navigation

- Aluminum brake cable conduit protects at the most frequent breakage point

- Shielded adjustable braking mechanism offers excellent durability and performance

FEATURES

- Lightweight, Aluminum frame with integrated suspension system is designed for durability and years of worry free use.
- 10-inch front and 8-inch rear heavy-duty wheels with black non-marring, soft-grip tires to traverse almost any terrain.
- Ergonomic handles – with top-down hidden braking system; Adjustable height in one-inch increments from 36" - 42.5" with audible "click" setting.
- Comes with a 24" x 13" thickly padded seat with a handle to fold and unfold side-to-side quickly and easily. This allows for the rollator to maneuver through narrow hallways and doors.
- Attractive front mesh tote attaches to the front of the frame for carrying personal belongings. It is easily attached and detached and can double as a shoulder bag when necessary.
- Folds up quickly and easily for transport or storage. Can easily fit in vehicles, storage closet, or public transport.

COLORS

Cobalt Blue



Brandywine



SPECIFICATIONS

Weight Capacity:	450 lb
Overall Width:	29" (W)
Seat Dimensions:	24" (W) x 13" (D)
Seat to Floor Height:	24"
Adjustable Handle Height:	36" – 42.5"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	23 lb
Carton Dimensions:	29.7" x 11.2" x 37.6"
Shipping Weight:	30.4 lb
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Wear Parts

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov

Item #	Description	UOM
RLASXRCB	Ascent XR Rollator, Cobalt Blue	each
RLASXRBW	Ascent XR Rollator, Brandywine	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.