

Dual Access Sling with Retention Straps

8C0X00

- Provides support from under thighs, behind the mid-to-low back, and around the torso with a hook and loop belt enclosure
- For use with patients who are seated requiring transfer to tub or toileting
- Optional side retention straps
- Divided, padded leg supports that can be applied and removed while in a seated position




Product Codes and Size Chart:



| Size | Product No. | A | C | D | E | F | G | H | I | J | K | L | M | N | Max SWL | Size Trim |
|------|-------------|----------------|----------------|----------------|-----------------|----------------|----------------|---------------|----------------|----------------|---------------|----------------|----------------|--------------|---------|-----------|
| XXS | 8C0700 | 25" (64cm) | 8½" (22cm) | 9½" (24cm) | 26¾" (68cm) | 13" (33cm) | 8¼" (22cm) | 8" (20cm) | 13" (33cm) | 6¼" (17cm) | 4¾" (12cm) | 10½" (27cm) | 18½" (47cm) | 6" (15cm) | 450lbs | Silver |
| XS | 8C0600 | 28¾" (72cm) | 8½" (22cm) | 10¼" (27cm) | 30" (76cm) | 13¾" (35cm) | 14" (36cm) | 8" (20cm) | 13¾" (35cm) | 7" (18cm) | 4¾" (12cm) | 10½" (27cm) | 18½" (47cm) | 6" (15cm) | 450lbs | Red |
| S | 8C0500 | 31½" (80cm) | 8½" (22cm) | 11¾" (30cm) | 34" (86cm) | 17¼" (44cm) | 14" (36cm) | 10" (25cm) | 17¼" (44cm) | 8" (20cm) | 4¾" (12cm) | 10½" (27cm) | 22" (56cm) | 8" (20cm) | 625lbs | Yellow |
| M | 8C0400 | 37" (94cm) | 8½" (22cm) | 11¾" (30cm) | 38" (97cm) | 20½" (52cm) | 14" (36cm) | 10" (25cm) | 20½" (51cm) | 12" (30cm) | 5½" (14cm) | 10½" (27cm) | 22" (56cm) | 8" (20cm) | 625lbs | Green |
| L | 8C0300 | 42" (107cm) | 10½" (27cm) | 13¾" (34cm) | 46" (117cm) | 21½" (55cm) | 14" (36cm) | 10" (25cm) | 21½" (55cm) | 13¼" (34cm) | 6" (15cm) | 10½" (27cm) | 22" (56cm) | 8" (20cm) | 625lbs | Blue |
| XL | 8C0200 | 46" (117cm) | 10½" (27cm) | 16½" (42cm) | 52¾" (134cm) | 24" (61cm) | 15¼" (39cm) | 10" (25cm) | 24" (61cm) | 15¼" (39cm) | 6" (15cm) | 10½" (27cm) | 22" (56cm) | 8" (10cm) | 1000lbs | Black |

Specifications:



Applications:

- From bed to commode
- From wheelchair to toilet
- Transfer to tub