



Company: <input style="width: 350px;" type="text"/>	Ship To: <input style="width: 350px;" type="text"/>
Account Number: <input style="width: 100px;" type="text"/> Date: <input style="width: 100px;" type="text"/>	Address: <input style="width: 350px;" type="text"/>
<input type="checkbox"/> Quote: <input type="checkbox"/> PO Number: <input style="width: 250px;" type="text"/>	City: <input style="width: 150px;" type="text"/> State: <input style="width: 100px;" type="text"/>
ATP: <input style="width: 350px;" type="text"/>	Zip: <input style="width: 100px;" type="text"/> Phone: <input style="width: 150px;" type="text"/>
Marked For: <input style="width: 350px;" type="text"/>	Email/Fax: <input style="width: 350px;" type="text"/>
Associated Quote/Order#: <input style="width: 250px;" type="text"/>	

Effective: 1/5/2026


Rabbit Up



A (Total Height) _____
 B (Armpit Height) _____
 C (Top of Pelvic Height) _____
 D (Buttock Height) _____
 E (Chest Width) _____
 F (Pelvic Width) _____
 Age _____
 Weight _____

Features, incl. in standard price


Stander Frame
Centre Bar
Brake on Castors
Angling Spindle




Mobile Stander Frame
Centre Bar
Driving Wheels
Angling Spindle
Closed Pushrim
0° Camber std. Angle
Spoke protectors
Central brake
Brake on Casters




Measurements



Size 1



Size 2



Size 3



Size 4

Product Specs as a standing frame

Overall width (G)	20 ¹ / ₄ "	21 ³ / ₄ "	23 ³ / ₄ "	23 ³ / ₄ "
Length (H)	27 ³ / ₄ "	31 ¹ / ₄ "	31 ¹ / ₄ "	36 ¹ / ₂ "
Total mass	19.8 lbs.	23.1 lbs.	24.2 lbs.	29.3 lbs.

Product Specs as a mobile stander

	w/26" wheels	w/32" wheels	w/32" wheels	w/36" wheels
Overall width (G)	25"	26 ³ / ₄ "	28 ³ / ₄ "	28 ¹ / ₄ "
Length (H)	24 ¹ / ₂ "	29 ¹ / ₂ "	29 ¹ / ₂ "	35 ¹ / ₂ "
Total mass	29.75 lbs.	34.1 lbs.	35.25 lbs.	40.75 lbs.

	Size 1	Size 2	Size 3	Size 4
Overall measurements				
Width, inside (G1)	12¾"	14¾"	16¾"	16¾"
Overall height (I) *	26 ¼" – 39 ¼"	33 ¾" – 46 ¾"	39" – 51 ½"	47 ¼" – 60 ¼"
Height from foot plate (I1) *	21 ¾" – 29 ¼"	29 ½" – 42 ½"	32 ½" – 46 ¾"	37 ¾" – 56"
Center bar angle (a1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Foot plate angle (a2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	6°	6°	4°
Center bar length	19¾"	27½"	31½"	39¼"

User Measurement Guide for Stander

Recommended user height	30" - 39"	39" - 54"	43" - 58"	48" - 71"
Max armpit height	29"	42"	46"	56"
User weight max	88.1 lbs.	110.2 lbs.	132.2 lbs.	154.3 lbs.
Recommended user age	1 - 3 yr./old	4 - 10 yr./old	5 - 12 yr./old	7 - 18 yr./old

* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

Rabbit Stander Options

Make a selection of either the Prone Stander or the Mobile Prone Stander

Rabbit-Up Prone Stander



- ☐ Rabbit-Up, Stander, Yellow, Size 1 (E0638) (8621220-00) \$1648
- ☐ Rabbit-Up, Stander, Blue, Size 2 (E0638) (8622220-00) \$1690
- ☐ Rabbit-Up, Stander, Grey, Size 3 (E0638) (8623220-00) \$1731
- ☐ Rabbit-Up, Stander, Silver, Size 4 (E0638) (8624220-00) \$1858

Rabbit Mobile Prone Stander



- ☐ 22" Wheels, Yellow, Size 1 (E0642) (8621222-21) \$3159
- ☐ 24" Wheels, Yellow, Size 1 (E0642) (8621223-21) \$3159
- ☐ 26" Wheels, Yellow, Size 1 (E0642) (8621224-21) \$3159
- ☐ 29" Wheels, Blue, Size 2 (E0642) (8622225-21) \$3957
- ☐ 32" Wheels, Blue, Size 2 (E0642) (8622226-21) \$3957
- ☐ 29" Wheels, Grey, Size 3 (E0642) (8623225-21) \$4000
- ☐ 32" Wheels, Grey, Size 3 (E0642) (8623226-21) \$4000
- ☐ 32" Wheels, Silver, Size 4 (E0642) (8624226-21) \$4126
- ☐ 36" Wheels, Silver, Size 4 (E0642) (8624227-21) \$4253

Chest Supports

Supports and Pad



Chest Support and Pad

- ☐ Chest Support and Pad, W: 7.75-10.75", Size 1 (861704) \$246
- ☐ Chest Support and Pad, W: 7.75-10.75", Size 2 and 3 (861702) \$246
- ☐ Chest Support and Pad, W: 9.75-12.5", Size 4 (861703) \$271



Filler Pad

Not required for Surrounding Chest Support

- ☐ Filler Pad, Small (861421) \$98
- ☐ Filler Pad, Large (861422) \$98



Extension for Center Bar

Available with Chest and Surrounding Chest Supports

- ☐ Extension for Center Bar (861480) \$98

Rabbit Up



Surrounding Chest Support and Options



Surrounding Chest Support

- ☐ Surrounding Chest Support, C:18.5-25.25", Small (861810) \$589
- ☐ Surrounding Chest Support, C: 26.5-30.75", Medium (861820) \$613
- ☐ Surrounding Chest Support, C: 34.25-44.75", Large (861830) \$649



Cover for Surrounding Chest Support

- ☐ Cover for Surrounding Chest Support, C:18.5-25.25", Small (869810-160) \$139
- ☐ Cover for Surrounding Chest Support, C: 26.5-30.75", Medium (869820-260) \$171
- ☐ Cover for Surrounding Chest Support, C: 34.25-44.75", Large (869830-360) \$182



Extra Pad for Surrounding Support

- ☐ Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225) \$75

Trays



Trays

- ☐ Tray, Plastic, Grey (862515) \$463
- ☐ Tray With Storage Box (862501) \$475
- ☐ Plastic Tray with Bowl and Cover (862525) \$672
Mounting Hardware required for Tray with Bowl

Tray Accessories

- ☐ Height Adjustment of Tray (862510) \$271
Required with Tray with Storage Box

Mounting Hardware

Required Selection for Tray with Bowl

- ☐ Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit) \$320

Abdominal Support



Abdominal Supports

- ☐ Abdominal Support (Each) (861411) \$193

Head Supports

Head Support, Head Support Cover, Connecting System for Head Support, Headrest Mounting Hardware, and Hip Support w/ Pommel are required for a complete Head Support

Head Support



- ☐ Anatomic Head Support, Small (5099001) \$285
- ☐ Anatomic Head Support, Medium (5099002) \$285
- ☐ Anatomic Head Support, Large (5099003) \$285

Head Support Cover


- ☐ Anatomic Head Support Cover, Small (9923-83) N/C
- ☐ Anatomic Head Support Cover, Medium (9945-83) N/C
- ☐ Anatomic Head Support Cover, Large (9967-83) N/C

Connecting System for Head Support




- ☐ Connecting System for Head Support, Size 1 and 2 (861424-3) \$310
- ☐ Connecting System for Head Support, Size 3 and 4 (861424-4) \$310

Headrest Mounting Hardware



- ☐ Headrest Mounting Hardware, Short (98619-25) \$250
- ☐ Headrest Mounting Hardware, Long (98619-26) \$256

Pelvic and Hip Control Options



Hip Support w/ Pommel

Front mount to center bar. Required with Head Support.

- ☐ Hip Support w/ Pommel, Size 1 (861501) \$607
- ☐ Hip Support w/ Pommel, Size 2 and 3 (861502) \$607
- ☐ Hip Support w/ Pommel, Size 4 (861503) \$672



Hip Support, With Rotation

- ☐ Hip Support, With Rotation, Size 2 and 3 (861612) \$426
- ☐ Hip Support, With Rotation, Size 4 (861613) \$463



Hip Support, Front Mount


- ☐ Hip Support, Front Mount, Size 1 (861601) \$283
- ☐ Hip Support, Front Mount, Size 2 and 3 (861602) \$319
- ☐ Hip Support, Front Mount, Size 4 (861603) \$352



Hip Support, Surrounding

Front mount to center bar

- ☐ Hip Support, Surrounding, Size 1 (861610) \$667
- ☐ Hip Support, Surrounding, Size 2 and 3 (861620) \$667
- ☐ Hip Support, Surrounding, Size 4 (861630) \$679



Pommel With Back Support


- ☐ Pommel With Back Support (861401) \$524

Knee Supports



Knee Support

- ☐ Knee Supports, Size 1 (Pair) (861301) \$336
- ☐ Knee Supports, Size 2 and 3 (Pair) (861302) \$426
- ☐ Knee Supports, Size 4 (Pair) (861303) \$426




Hinged Knee Supports

- ☐ Hinged Knee Supports, Size 1 (Pair) (861311) \$524
- ☐ Hinged Knee Supports, Size 2 and 3 (Pair) (861312) \$524
- ☐ Hinged Knee Supports, Size 4 (Pair) (861313) \$524

Footplate and Sandals Options

Flip-up Footplate required for all Sandals, except the Individual Footrest w/ Hinged Sandals

Footplates and Options




Flip-Up Footplate

- ☐ Flip-Up Footplate, Size 1 and 2 (861211) \$265
- ☐ Flip-Up Footplate, Size 2, 3 and 4 (861212) \$265



Formed Footplates

- ☐ Formed Footplates, Size 0 (Pair) (86149) \$198
- ☐ Formed Footplates, Size 1 and 2 (Pair) (86147) \$198
- ☐ Formed Footplates, Size 3 (Pair) (86148) \$198



Footplate Straps

- ☐ Foot Straps, Plastic Lock (Pair) (85096) (K0040) \$42
- ☐ Footplate Straps, Velcro (Pair) (85102) \$59



Individual Footrest w/ Hinged Sandals

- ☐ Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201) \$475
- ☐ Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202) \$475

Sandals and Options



Sandals

- ☐ Sandals for Rabbit Up/ Toucan, Small (869490-RT) \$233
- ☐ Sandals for Rabbit Up/ Toucan, Medium (869491-RT) \$244
- ☐ Sandals for Rabbit Up/ Toucan, Large (869492-RT) \$251



Sandals

- ☐ Hinged Sandals w/ Straps, Size 1 (Pair) (861221) \$162
- ☐ Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222) \$162
- ☐ Hinged Sandals w/ Straps, Size 4 (Pair) (861223) \$283



Heel Stops

- ☐ Heel Stops, Small (Pair) (861251) \$61
- ☐ Heel Stops, Large (Pair) (861252) \$71

Conversion Kits

Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit. Conversion kit is used to convert Prone Stander into Mobile Prone Stander.



Driving Wheel Conversion Kits

- ☐ Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209) \$694
- ☐ Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210) \$694
- ☐ Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211) \$694
- ☐ Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212) \$1071
- ☐ Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213) \$1071
- ☐ Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214) \$1071



Central Brake System

- ☐ Central Brake System, Size 1 (134215) \$349
- ☐ Central Brake System, Size 2 and 3 (134216) \$349
- ☐ Central Brake System, Size 4 (134218) \$349



Camber Kit Rabbit Up

- ☐ Camber Kit Rabbit Up 0° (134200) \$63
- ☐ Camber Kit Rabbit Up +3/-3° (134201) \$63
- ☐ Camber Kit Rabbit Up +6/-6° (134202) \$63