

Fusion**ONE**[™]

Model 3002-500-105
Model 3002

For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE Exam Table offers the features healthcare providers need most.

This hi-lo exam table with stirrups enables simple ADA access and ease of examination for patients of all ages and abilities.



Model 3002-500-105
includes Premium upholstery
(shown with optional armrests)

Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Large, reversible side drawer and large front drawer
- › Drain pan: plastic 1.5' deep (3.81cm)
- › Removable, seamless upholstered top with PreFixx[®] Protective Finish
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Durable, steel construction
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

FusionONE™

Model 3002-500-105
Model 3002

Specifications

- › Height: min. 19" (48.3 cm) at seat
max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x
45.7 cm D)
- › Storage Drawer (Front): 18.0" W x
12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x
21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International model: 3002-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x
64" L x 31" H (76 cm x 163 cm x 79 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket
capabilities for #48955 (patient right)



Model 3002 includes
Classic upholstery

Upholstery and Colors

Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
True Graphite	Classic Black	Fine Wine	Blue Skies
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2022 United Metal Fabricators Inc. dba UMF Medical 08.22 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1371LT Rev A.