

Universal Mesh

- Full body supportive sling
- Divided padded leg sections
- Ideal for general transfers
- Small aperture for hygiene purposes
- Color coded destination loop straps as standard
- Mesh (100% polyester open weave mesh; high density, durable, used for bathing as quick drying.)

M921501



Product Codes & Specifications:



100% Polyester Mesh		A	B	C	D
XS	M921501-XS	30 - 34" (76 - 86cm)	17½ - 19¾" (60 - 66cm)	23½" (60cm)	21¼" (54cm)
S	M921501-S	34 - 38" (86 - 96cm)	19¾ - 22½" (66 - 74cm)	27½" (70cm)	25¼" (64cm)
M	M921501-M	38 - 42½" (96 - 108cm)	22½" - 25¼" (74 - 80cm)	35½" (90cm)	29½" (75cm)
L	M921501-L	42½ - 51¼" (108 - 130cm)	25¼ - 27¼" (80 - 88cm)	39½" (100cm)	30¾" (78cm)
XL	M921501-XL	51¼ - 56¾" (130 - 144cm)	27¼ - 29¼" (88 - 100cm)	43½" (110cm)	33" (84cm)
XXL	M921501-XXL	56¾ - 63" (144 - 160cm)	29¼" - 31½" (96 - 104cm)	51¼" (130cm)	37½" (95cm)

