



Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base model

E1236 coding still applies for EIF7 with sizes 14" or less. 300 lbs weight capacity

<input type="checkbox"/>	170WM112	K0005	Quickie QS5 X	\$2,550.00
<input type="checkbox"/>	170WM113	K0005	Quickie QS5 X - Transit included	\$2,924.00

Frame Style

See figure 1. Measure from front frame to floor

<input type="checkbox"/>	170FS31		Swing Away Frame	NC
<input type="checkbox"/>	170FS43		Fixed Front 75°	NC
<input type="checkbox"/>	170FS44		Fixed Front 85°	NC

Frame Inset (Fixed Front Only)

See **FIGURE 2** Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2. N/A for 14" and 15" Wide using composite angle adj footrest. **75° only available with 0" inset.**

<input type="checkbox"/>	170FF1		0"	NC
<input type="checkbox"/>	170FF2		1"	NC

Transit Option

ISO 7176-19 Transit Approved

<input type="checkbox"/>	170TO49	K0108	Transit System	\$374.00
--------------------------	---------	-------	----------------	----------

Frame Options

Ergo Touch Points is a rubber cover on the knee bend area of the frame tubes

<input type="checkbox"/>	1709A11		Ergo Touch Points Swing Away hangers or Fixed front	\$150.00
<input type="checkbox"/>	1709A12		Frame Fold-N-Lock mechanism	NC
<input type="checkbox"/>	1709A06		Twin LED lights	\$248.00



QS5 X - Swing Away



QS5 X - Fixed Front

Power Assist

To order any Empulse Power Assist product please see the Order Form for Part #'s

Option	HCPCS	DESCRIPTION	RETAIL
<input type="checkbox"/>	209MW91	E0986 Empulse R90 Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 2.5Ahr Removable Battery	\$7,300.00
<input type="checkbox"/>	209MW92	E0986 Empulse M90 M90 Wheels, Dual Controller, Two M90 Batteries, Mounting Hardware, Quick Release Axles, and Device Charger	\$7,950.00
<input type="checkbox"/>	147MW95	E0986 Empulse F35 Includes: F35 Stick, Subframe mount, Frame clamps kit, and Charger	\$6,500.00



2. SIZE / SEATING

Frame Width

See **FIGURE 2** Measured from outside of seat tube to the outside of opposite seat tube.

<input type="checkbox"/>	170W4	13"	NC
<input type="checkbox"/>	170W5	14"	NC
<input type="checkbox"/>	170W6	15"	NC
<input type="checkbox"/>	170W7	16"	NC
<input type="checkbox"/>	170W8	17"	NC
<input type="checkbox"/>	170W9	18"	NC
<input type="checkbox"/>	170W10	19"	NC
<input type="checkbox"/>	170W11	E2201 20"	\$583.00

Seat Depth

See **Seat Squeeze Matrix** in Appendix section. Measured from back canes to front edge of upholstery. Frame Fold-N-Lock mechanism NA with 13"-14".

<input type="checkbox"/>	170SD6	13"	NC
<input type="checkbox"/>	170SD7	14"	NC
<input type="checkbox"/>	170SD8	15"	NC
<input type="checkbox"/>	170SD9	16"	NC
<input type="checkbox"/>	170SD10	17"	NC
<input type="checkbox"/>	170SD11	18"	NC
<input type="checkbox"/>	170SD12	19"	NC
<input type="checkbox"/>	170SD13	E2203 20"	\$641.00

Seat Sling

Standard seat sling includes utility bag sewn at the bottom of upholstery.


<input type="checkbox"/>	170SL1	Standard Seat Sling	NC
<input type="checkbox"/>	170SL6	Omit Seat Sling	-\$15.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	170JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____	Depth: _____



3. CASTERS / SEAT HEIGHT

Front Seat Height

<input type="checkbox"/>	170SH41	14.0"	NC	<input type="checkbox"/>	170SH18	18.0"	NC
<input type="checkbox"/>	170SH42	14.5"	NC	<input type="checkbox"/>	170SH20	18.5"	NC
<input type="checkbox"/>	170SH12	15.0"	NC	<input type="checkbox"/>	170SH21	19.0"	NC
<input type="checkbox"/>	170SH13	15.5"	NC	<input type="checkbox"/>	170SH23	19.5"	NC
<input type="checkbox"/>	170SH14	16.0"	NC	<input type="checkbox"/>	170SH24	20.0"	NC
<input type="checkbox"/>	170SH15	16.5"	NC	<input type="checkbox"/>	170SH26	20.5"	NC
<input type="checkbox"/>	170SH16	17.0"	NC	<input type="checkbox"/>	170SH27	21.0"	NC
<input type="checkbox"/>	170SH17	17.5"	NC	<input type="checkbox"/>	170SH29	21.5"	NC
				<input type="checkbox"/>	170SH36	22.0"	NC

Caster Wheels (option A)

<input type="checkbox"/>	170CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	170CW44	E2219 4" x 1.25" Semi Pneumatic	\$63.00
<input type="checkbox"/>	170CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	170CW7	E2219 5" x 1.5" Semi-Pneumatic	\$100.00
<input type="checkbox"/>	170CW34	6" Polyurethane	NC
<input type="checkbox"/>	170CW33	E2214 E2215 6" Pneumatic	\$88.00
<input type="checkbox"/>	170CW35	E2219 6" x 1.5" Semi Pneumatic	\$100.00
<input type="checkbox"/>	170CW16	7" Polyurethane	NC
<input type="checkbox"/>	170CW32	8" Polyurethane	NC


Caster Fork

<input type="checkbox"/>	170CF5	Multiple Position	NC
<input type="checkbox"/>	170CF69	Multiple Position Single Sided Fork	\$197.00
<input type="checkbox"/>	170CF53	E1015 Frog Legs 0 weight to 100 lbs.	\$507.00
<input type="checkbox"/>	170CF54	E1015 Frog Legs 101 lbs. weight to 190 lbs.	\$507.00
<input type="checkbox"/>	170CF55	E1015 Frog Legs 191 lbs. weight to 265 lbs.	\$507.00

Caster Wheels colors (option B)

Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

<input type="checkbox"/>	170CW64	4" x 1.5" Aluminum Soft Roll	\$216.00
<input type="checkbox"/>	170CW66	5" x 1.5" Aluminum Soft Roll	\$227.00
<input type="checkbox"/>	170CW68	6" x 1.5" Aluminum Soft Roll	\$239.00
	Colors:	<input type="checkbox"/> Match anodize color package	Package price
		<input type="checkbox"/> Match frame color	\$57.00
		<input type="checkbox"/> Black anodize caster only	\$58.00



Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
13"-15"	2.0"
16"-19"	2.5"
20"	3.0"

4. REAR WHEELS / SEAT HEIGHT

Rear Seat Height						
<input type="checkbox"/>	170RH13	14.0"	NC	<input type="checkbox"/>	170RH19	17.0" NC
<input type="checkbox"/>	170RH14	14.5"	NC	<input type="checkbox"/>	170RH26	17.5" NC
<input type="checkbox"/>	170RH15	15.0"	NC	<input type="checkbox"/>	170RH20	18.0" NC
<input type="checkbox"/>	170RH16	15.5"	NC	<input type="checkbox"/>	170RH27	18.5" NC
<input type="checkbox"/>	170RH17	16.0"	NC	<input type="checkbox"/>	170RH21	19.0" NC
<input type="checkbox"/>	170RH18	16.5"	NC	<input type="checkbox"/>	170RH37	19.5" NC
				<input type="checkbox"/>	170RH22	20.0" NC

Axle			
<input type="checkbox"/>	170AX1	Quick Release Steel	NC
<input type="checkbox"/>	170AX11	Quad Quick Release	\$144.00

Rear Wheel Type

Max weight capacity on Proton wheels is 220 lbs.

<input type="checkbox"/>	170RW103	20" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	170RW104	22" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	170RW102	24" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	170RW90	20" Lite Spoke Wheel	\$234.00
<input type="checkbox"/>	170RW91	22" Lite Spoke Wheel	\$234.00
<input type="checkbox"/>	170RW92	24" Lite Spoke Wheel	NC
<input type="checkbox"/>	170RW93	25" Lite Spoke Wheel	NC
<input type="checkbox"/>	170RW94	26" Lite Spoke Wheel	\$234.00
<input type="checkbox"/>	170RW95	24" Spinergy SPOX Black (silver hub)	\$985.00
<input type="checkbox"/>	170RW96	25" Spinergy SPOX Black (silver hub)	\$985.00
<input type="checkbox"/>	170RW98	24" Spinergy LX Black (silver hub)	\$1,237.00
<input type="checkbox"/>	170RW99	25" Spinergy LX Black (silver hub)	\$1,237.00
<input type="checkbox"/>	170RW148	24" Proton Wheel	\$234.00
<input type="checkbox"/>	170RW149	25" Proton Wheel	\$234.00
<input type="checkbox"/>	170RW25	Omit Rear Wheel and Axle *	-\$55.00

Spinergy Spoke Colors (Optional)

Please select one of the following Spinergy wheel options first. Includes black hub.

Spinergy Spox spokes color:

- Black Blue Red White Yellow

Spinergy LX spokes color:

- Black Blue Green Orange Pink
 Purple Red White Yellow

Rear Wheel Tire

<input type="checkbox"/>	170RT1		Pneumatic	NC
<input type="checkbox"/>	170RT24	E2211 E2213	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$114.00
<input type="checkbox"/>	170RT3		Full Poly	NC
<input type="checkbox"/>	170RT38	E2218	SHOX Performance Poly	\$110.00
<input type="checkbox"/>	170RT25	E2218	Primo Orion Poly	\$126.00
<input type="checkbox"/>	170RT39	E2211 E2212	High-Pressure Clincher Black	\$110.00
<input type="checkbox"/>	170RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$164.00

Wheel Dia.	RSH Range	
	Min Squeeze (0° seat angle)	Max Squeeze (9° seat angle)
20"	14.0"-18.0"	14.0"-17.5"
22"	14.0"-19.0"	14.0"-18.5"
24"	15.0"-19.5"	14.5"-19.0"
25"	15.5"-20.0"	15.0"-19.5"
26"	16.0"-20.0"	15.5"-20.0"

NEW! Handrim

See page 6 for matrices.

<input type="checkbox"/>	170HR1	Aluminum Anodized	NC
<input type="checkbox"/>	170HR29	Aluminum Powder coated colored	\$82.00
<input type="checkbox"/>	170HR4	Plastic Coated	\$129.00
<input type="checkbox"/>	170HR22	Natural Fit - Standard Grip	\$549.00
<input type="checkbox"/>	170HR23	Natural Fit - Super Grip	\$611.00
<input type="checkbox"/>	170HR28	Q-Grip Neoprene Coated	\$408.00
<input type="checkbox"/>	170HR39	Ellipse 3R	\$500.00
<input type="checkbox"/>	170HR40	Carbolife - Gekko	\$535.00
<input type="checkbox"/>	170HR41	Carbolife - Curve L Tetra grip	\$575.00
<input type="checkbox"/>	170HR7	Omit Handrim <i>Option may not be available on all wheels</i>	NC

Camber

<input type="checkbox"/>	170CI1	0°	NC
<input type="checkbox"/>	170CI13	2° (default)	NC
<input type="checkbox"/>	170CI3	4°	NC

Center of Gravity Preset

<input type="checkbox"/>	170CC20	1.0" (default)	NC
<input type="checkbox"/>	170CC15	1.5"	NC
<input type="checkbox"/>	170CC17	2.0"	NC
<input type="checkbox"/>	170CC18	2.5"	NC
<input type="checkbox"/>	170CC19	3.0"	NC

5. WHEEL LOCKS

Wheel Lock

No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth

<input type="checkbox"/>	170WL1	High Mount Push	NC
<input type="checkbox"/>	170WL2	High Mount Pull	NC
<input type="checkbox"/>	170WL45	Ergo Scissor	NC
<input type="checkbox"/>	170WL88	Fold Down Push (Knee lever)	\$86.00
<input type="checkbox"/>	170WL30	Compact	\$164.00
<input type="checkbox"/>	170WL10	Omit	NC

Wheel Lock Option

<input type="checkbox"/>	170WL12	Do Not Mount	NC
--------------------------	---------	--------------	----

Extension Handle

N/A w/ ergo scissor or compact wheel locks

<input type="checkbox"/>	170EH7	E0961	Extension Handle (Pair)	\$88.00
--------------------------	--------	-------	-------------------------	---------

6. HANGER & FOOTREST									
----------------------	--	--	--	--	--	--	--	--	--

Swing Away Hangers				
<i>Swing Away frame only</i>				
<input type="checkbox"/>	170H145		60° Swing In/Out Hanger	NC
<input type="checkbox"/>	170H139		70° Swing In/Out Hanger	NC
<input type="checkbox"/>	170H141		80° Swing In/Out Hanger	NC
<input type="checkbox"/>	170H160	E0990	Elevating Hanger - Ext Mount.	\$367.00
<input type="checkbox"/>	170H159	K0053	Articulating Hanger - Ext Mount	\$431.00
<input type="checkbox"/>	170H32		Omit Hanger and Footplate	-\$15.00

Footplates				
<i>Heel loop or leg strap not included, See Footplate Options Section. Platform footplates flip up (sideways).</i>				
<input type="checkbox"/>	170F1		Composite	STD
<input type="checkbox"/>	170F271	K0040 E0951	Composite Angle Adj-(16"-22"w)	\$211.00
<input type="checkbox"/>	170F5		Aluminum (ELR & ALR)	NC
<input type="checkbox"/>	170F3	K0040	Alum Adult Angle Adjustable - Ext Mour	\$221.00
<input type="checkbox"/>	170F259	K0040	Multi Angle Adustable <input type="checkbox"/> Black NC <input type="checkbox"/> Std Frame color\$51	\$221.00
<input type="checkbox"/>	170F50	K0040	Platform Footrest <input type="checkbox"/> Black NC <input type="checkbox"/> Std Frame color\$51	\$110.00
<input type="checkbox"/>	170F335	K0040	Platform Ultra Lite Footrest	\$134.00
<input type="checkbox"/>	170F336		Platform Carbon Fiber footrest	\$158.00
<input type="checkbox"/>	170F7		Omit Footplate and Ext. Tube	-\$15.00

Footrest Extension Tubes				
<i>See Extension matrix for leg length range</i>				
<input type="checkbox"/>	170E1		Short	NC
<input type="checkbox"/>	170E2		Medium	NC
<input type="checkbox"/>	170E5		Long	NC

Footrest Options				
<input type="checkbox"/>	170HL1	E0951	Heel Loop <i>Swing away only</i>	\$47.00
<input type="checkbox"/>	170HL5	K0038	Leg Strap	\$47.00
<input type="checkbox"/>	170DI04		Extra Leg Strap	\$47.00
<input type="checkbox"/>	170IG2		Impact Guard Neoprene	\$105.00
<input type="checkbox"/>	170DI06		Platform Position Plates	\$58.00
<input type="checkbox"/>	170F35		Auto-Folding Platform Carbon fiber only	\$58.00

Articulating / Elevating Hanger Cover				
<input type="checkbox"/>	170DI47		Neoprene Cover	\$397.00
<input type="checkbox"/>	170DI48		Gel Pad Cover	\$467.00

Lower Leg Length (LLL)									
------------------------	--	--	--	--	--	--	--	--	--

*Lower Leg Length is **measured from seat sling to the top rear of the footrest**. See Appendix on Page 7. A minimum 2" ground clearance is required below the footrest. If the selected configuration results in less than 2" of ground clearance the footrest will be raised to 2" ground clearance, which may alter the lower leg length measurement. Available lower leg length is dependent on frame selection including adapter spacers. **Reference LLL matrix in page 7***

- 10" 11" 12" 13" 14" 15" 16" 17" 18" 19"

7. BACKREST / SEATING									
-----------------------	--	--	--	--	--	--	--	--	--

Back Type				
<input type="checkbox"/>	170BT123		Standard Fix (90°)	NC
<input type="checkbox"/>	170BT47		Angle Adjustable Back (80° - 105°)	NC
Select desired setting <input type="checkbox"/> 80° <input type="checkbox"/> 85° <input type="checkbox"/> 90° <input type="checkbox"/> 100°				
<input type="checkbox"/>	170BT124		Fixed Half Fold down (to the front)	\$216.00
<input type="checkbox"/>	170BT43		Lumbar curve bend	\$87.00
Select desired setting <input type="checkbox"/> 5° forward <input type="checkbox"/> 5° backward				

Back Height Setting				
<i>Half Fold Down N/A with the 10" - 13" back heights. Lumbar curve backrest with integral push handles available with 12" thru 20" back heights. Lumbar curve backrest with Fold Down or Omit push handles available with 14" thru 20" back heights.</i>				
<input type="checkbox"/>	170B161	Setting	10" - 12" Back Height <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12"	NC
<input type="checkbox"/>	170B162	Setting	12" - 14" Back Height <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"	NC
<input type="checkbox"/>	170B163	Setting	14" - 16" Back Height <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16"	NC
<input type="checkbox"/>	170B164	Setting	16" - 18" Back Height <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18"	NC
<input type="checkbox"/>	170B165	Setting	18" - 20" Back Height <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20"	NC

Push Handles				
<i>Std EXO Uph N/A with Fold Down. Ht Adj Push Handles only available with tension adjustable EXO / EXO Pro upholstery.</i>				
<input type="checkbox"/>	170PH6		Integral Push Handles	NC
<input type="checkbox"/>	170PH5		Fold Down Push Handles <i>Available w/12"-20" back height</i>	\$116.00
<input type="checkbox"/>	170PH21		Height adjustable push Handles <i>Available w/ 12"-20" back height. N/A w/Half fold down back</i>	\$175.00
<input type="checkbox"/>	170PH25		Height adjustable push Handles for JAY Backs <i>N/A w/Half fold down back</i>	\$175.00
<input type="checkbox"/>	170PH22		Omit Push Handles - Straight Post	NC

Backrest Upholstery				
<i>EXO standard and Tension adjustable upholstery binding is always black.</i>				
<input type="checkbox"/>	170BU95		Standard Upholstery EXO	NC
Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00				
<input type="checkbox"/>	170BU70	E2611	Backrest Tension Adjust Ballistic <i>*Reference 186BU70 for billing</i>	\$367.00
<input type="checkbox"/>	170BU89	E2611	Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i>	\$387.00
Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00				
<input type="checkbox"/>	170BU90	E2611	Backrest Tension Adj EXO Pro 3DX <i>*Reference 186BU71 for billing</i>	\$387.00
Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00				
<input type="checkbox"/>	170BU14		Omit Back Upholstery	-\$15.00

Military Service Patches (Optional)				
<input type="checkbox"/>	170VU2		Army Patch	\$116.00
<input type="checkbox"/>	170VU4		Air Force Patch	\$116.00
<input type="checkbox"/>	170VU12		Coast Guard Patch	\$116.00
<input type="checkbox"/>	170VU1		Marines Patch	\$116.00
<input type="checkbox"/>	170VU3		Navy Patch	\$116.00
<input type="checkbox"/>	170VU5		Purple Heart Patch	\$116.00
<input type="checkbox"/>	170VU9		Purple Heart with Army Patch	\$181.00
<input type="checkbox"/>	170VU11		Purple Heart with Air Force Patch	\$181.00
<input type="checkbox"/>	170VU13		Purple Heart with Coast Guard Patch	\$181.00
<input type="checkbox"/>	170VU8		Purple Heart with Marines Patch	\$181.00
<input type="checkbox"/>	170VU10		Purple Heart with Navy Patch	\$181.00

JAY Back (Optional)				
<i>To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s</i>				
<input type="checkbox"/>	170JB01		J3 Back (See J3 Order Form)	Price List

8. ARMREST / SIDE GUARDS
Padded Swing Away

Some limitations when used w/ angle adj back in the rear most setting. 10"-12" height ranges N/A w/Swing away anti tip. Please select a setting, will default to 11" if no selection is made.

<input type="checkbox"/>	170AA162		Padded Swing-Away Fixed height (13")	NC
<input type="checkbox"/>	170AA18	E0973	Padded Swing-Away STD (9" to 14")	\$263.00
Setting <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"				

Armrest - Height Adjustable

<input type="checkbox"/>	170AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$333.00
<input type="checkbox"/>	170AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$333.00
<input type="checkbox"/>	170AA65	E0973	Flip Back Height Adj <i>LITE</i> (8.5"-13.5")	\$367.00

Armpad - For Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

<input type="checkbox"/>	170AM23		Desk Length - Classic (10")	NC
<input type="checkbox"/>	170AM24		Full Length - Classic (14")	NC
<input type="checkbox"/>	170AM27		Ergo Pad (9.5"-13.5")	NC

Omit Armrest

<input type="checkbox"/>	170AA14		Omit Armrest	NC
--------------------------	---------	--	--------------	----

9. KEY OPTIONS & ACCESSORIES
Anti-Tip

<input type="checkbox"/>	170AT1	E0971	Anti Tips (Pair)	\$175.00															
<input type="checkbox"/>	170AT	E0971	Anti Tip swing away <i>N/A Transit wheels</i>																
<table border="0" style="width:100%;"> <tr> <td style="width:10%; text-align: center;"><input type="checkbox"/></td> <td style="width:10%;">170AT88</td> <td style="width:10%;"></td> <td style="width:50%;">Right Mount</td> <td style="width:10%; text-align: right;">\$102.00</td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td>170AT89</td> <td></td> <td>Left Mount</td> <td style="text-align: right;">\$102.00</td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td>170AT90</td> <td></td> <td>Pair</td> <td style="text-align: right;">\$204.00</td> </tr> </table>					<input type="checkbox"/>	170AT88		Right Mount	\$102.00	<input type="checkbox"/>	170AT89		Left Mount	\$102.00	<input type="checkbox"/>	170AT90		Pair	\$204.00
<input type="checkbox"/>	170AT88		Right Mount	\$102.00															
<input type="checkbox"/>	170AT89		Left Mount	\$102.00															
<input type="checkbox"/>	170AT90		Pair	\$204.00															

Additional Accessories

<input type="checkbox"/>	170A1		Backpack - Adult	\$88.00
<input type="checkbox"/>	170A2		Backpack - Kids	\$70.00
<input type="checkbox"/>	1701A01		Seat Pouch Black	\$58.00

Additional Cover Accessories

<input type="checkbox"/>	1708A02		Gel Cover Hanger/Frame	\$129.00
<input type="checkbox"/>	1708A01		Gel Cover Swing Away Arm	\$175.00

Positioning Belts

<input type="checkbox"/>	170ZZ15	E0978	JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$101.00
<input type="checkbox"/>	170ZZ26	E0978	JAY 2PT Belt 2" Push Btn 9.5" Pad	\$101.00
<input type="checkbox"/>	170ZZ14	E0978	JAY 2PT Belt 1.5" Push Btn No Pad	\$48.00
<input type="checkbox"/>	170A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$139.00
<input type="checkbox"/>	170A120		Positioning Belt Slide Clamps	\$58.00

Side Guard

Side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist. Aluminum & carbon side guards some limitations with low RSH and large rear wheels

<input type="checkbox"/>	170SG33	K0108	Plastic Standard	\$169.00
<input type="checkbox"/>	170SG29	K0108	Plastic Low	\$169.00
<input type="checkbox"/>	170SG35	K0108	Aluminum	\$286.00
<input type="checkbox"/>	170SG36	K0108	Aluminum w/ Fender	\$321.00
<input type="checkbox"/>	170SG12	K0108	Carbon Fiber	\$537.00
<input type="checkbox"/>	170SG41	K0108	Carbon Fiber Long w/LG Fender	\$607.00

Side Guard - Color Option

Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber

<input type="checkbox"/>	170GE01		Black	NC
<input type="checkbox"/>	170GE02		Matching Frame- Standard Colors	NC
<input type="checkbox"/>	170GE03		Matching Frame- Personalized Color	\$58.00

Spoke Guards

Only available w/ Lite Spoke wheels.

<input type="checkbox"/>	1704A01	K0065	Quickie Black Spoke Guards	\$134.00
<input type="checkbox"/>	1704A02	K0065	Quickie Clear Spoke Guards	\$134.00

Miscellaneous

Transit wheels N/A with Swing away anti tip.

<input type="checkbox"/>	1708A03		Caddy	\$58.00															
<input type="checkbox"/>	170A2100	E2207	Crutch Holder	\$227.00															
<input type="checkbox"/>	1708A21		Tool Kit	\$63.00															
<input type="checkbox"/>	170A999		Quickie Water Bottle & Holder	\$63.00															
<input type="checkbox"/>	1708A47		Mobile phone pocket	\$30.00															
<input type="checkbox"/>	170FO4	E0971	Transit wheels (Pair)	\$122.00															
<input type="checkbox"/>	170AT	E0971	Tip Assist																
<table border="0" style="width:100%;"> <tr> <td style="width:10%; text-align: center;"><input type="checkbox"/></td> <td style="width:10%;">170AT83</td> <td style="width:10%;"></td> <td style="width:50%;">Right Mount</td> <td style="width:10%; text-align: right;">\$36.00</td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td>170AT84</td> <td></td> <td>Left Mount</td> <td style="text-align: right;">\$36.00</td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td>170AT91</td> <td></td> <td>Pair</td> <td style="text-align: right;">\$71.00</td> </tr> </table>					<input type="checkbox"/>	170AT83		Right Mount	\$36.00	<input type="checkbox"/>	170AT84		Left Mount	\$36.00	<input type="checkbox"/>	170AT91		Pair	\$71.00
<input type="checkbox"/>	170AT83		Right Mount	\$36.00															
<input type="checkbox"/>	170AT84		Left Mount	\$36.00															
<input type="checkbox"/>	170AT91		Pair	\$71.00															

10. COLOR

Personalized Crossbrace Color

Option has an extended leadtime. Frame colors only, NA with Personalized frame colors.
If no color pick crossbrace will be Gloss Black.

<input type="checkbox"/>	170GS15	Crossbrace Color package	\$250.00
--------------------------	---------	--------------------------	----------

NEW! Standard Frame Colors

Handrim color option only available if Aluminum Coated colored option was selected.

		Frame	X-brace	Handrim
170C214	Black (Gloss)	<input type="checkbox"/>	STD	<input type="checkbox"/>
170C215	Black (Matte)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C211	Celestial Sapphire	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C212	Desert Sage	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C213	Desert Tan	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C216	Metallic Shadow	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C217	Midnight Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C218	Mystic Purple	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C219	Petal Pink	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C220	Pink Power	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C221	Red Pepper	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C222	Royal Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C223	Salute Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C224	Scarlet Noir	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C225	Silver Ore	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C226	Solar Orange	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C227	Starlight White	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C228	Storm Gray	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C229	Turbo Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
170C230	Velocity Yellow	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Personalized Frame Color

Only frame will have the personalized color; Crossbrace, handrim, footplate and Aluminum side guard not included (unless selected).

<input type="checkbox"/>	170C112	Desert Camouflage	\$460.00
<input type="checkbox"/>	170C43	Mossy Oak Camouflage	\$460.00
<input type="checkbox"/>	170C42	Stars & Stripes	\$460.00
<input type="checkbox"/>	170C108	Zebra (Choose above base color) Std frame colors with White unless above base paint is chosen N/A w/ Black & Matte Black)	\$460.00

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

<input type="checkbox"/>	170D16	Black	NC
<input type="checkbox"/>	170D17	Blue	NC
<input type="checkbox"/>	170D18	Orange	NC
<input type="checkbox"/>	170D19	Red	NC
<input type="checkbox"/>	170D20	Silver	NC

Anodized Color Package

Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog suspension Legs.

<input type="checkbox"/>	170EB08	Blue	\$187.00
<input type="checkbox"/>	170EB13	Gold	\$187.00
<input type="checkbox"/>	170EB09	Lime Green	\$187.00
<input type="checkbox"/>	170EB10	Purple	\$187.00
<input type="checkbox"/>	170EB11	Red	\$187.00
<input type="checkbox"/>	170EB12	Ti Grey	\$187.00

Built For Me

To inquire about any unique personalized for Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011



11. APPENDIX

Wheel - Handrims Matrix

Handrims	Wheels			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Aluminum Ano / Painted	20,22,24	20,22,24,25,26	24,25	24,25
Plastic Coat	20,22,24	22,24	24,25	-
Natural Fit - Standard Grip	24	22,24,25,26	24,25	-
Natural Fit - Super Grip	24	22,24,25,26	24,25	-
Q Grip Neoprene Coated	-	24,25	24,25	-
Ellipse 3R	24	24,25	24,25	24,25
Carbolife - Gekko	-	24,25	24,25	24,25
Carbolife - Curve L Tetra grip	-	24,25	24,25	24,25
Omit Handrim	20,22,24	-	24,25	-

Wheel - Tire Matrix

Tire	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Pneumatic	20,22,24	20,22,24,26	24	24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	24	-
Full Poly	20,22,24	20,22,24	24	-
SHOX Performance Poly	20,22,24	20,22,24,25	24,25	-
Primo Orion Poly	20,22,24	20,22,24	24	-
High-Pressure Clincher Blk	24	22,24,26	24	24
Schwalbe Marathon	22, 24	24,25,26	24,25	24,25

20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm

11. APPENDIX

Note 1: These Front Seat Height (FSH) matrices contain rounded estimates and assume flat seat heights (0 deg seat angle). Increasing seat angle can affect FSH setting.
Note 2: Highlighted cells indicate all valid combinations the Configurator will possibly choose based on FSH and Caster selection.
Note 3: Blank spaces (denoted with a dash "-") indicate at least one of two things: A) The caster won't fit in that fork hole, and/or B) The fork hole may not exist.
Note 4: For all forks except Frog Legs, Hole Position 1+ indicates that a 10mm spacer is used above the fork bearings to achieve additional FSH. For all other hole positions, the 10mm spacer sits below the fork bearings, except for the Frog Legs which can only have the spacer on the top side of the fork.
Note 5: The lowest FSH 14.0" is only available with the 4" Low Profile Poly Caster (Multi 90 hole 2 is not obtainable with 4"x1.25" Semi Pneu or 4"x1.5" Alum Soft Roll).

FSH for 4" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1 1+	
Hemi	Multi 90	-	-	14.0	14.5	15.0
	Multi 115	-	14.5	15.0	15.5	16.0
	Multi 131	-	15.0	15.5	16.0	16.5
	Multi 158	-	16.0	16.5	17.0	17.5
	SS Low	-	-	14.0	14.5	15.0
	SS Std	-	14.5	15.0	15.5	16.0
	Frog Legs	-	15.5	16.0	16.5	-
Short	Multi 90	-	-	15.5	16.0	16.0
	Multi 115	-	16.0	16.5	17.0	17.0
	Multi 131	-	16.5	17.0	17.5	18.0
	Multi 158	-	17.5	18.0	18.5	19.0
	SS Low	-	-	15.5	16.0	16.0
	SS Std	-	16.0	16.5	17.0	17.0
	Frog Legs	-	17.0	17.5	18.0	-
Tail	Multi 90	-	-	16.5	17.0	17.5
	Multi 115	-	17.0	17.5	18.0	18.5
	Multi 131	-	17.5	18.0	18.5	19.0
	Multi 158	-	18.5	19.0	19.5	20.0
	SS Low	-	-	16.5	17.0	17.5
	SS Std	-	17.0	17.5	18.0	18.5
	Frog Legs	-	18.0	18.5	19.0	-

FSH for 5" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1 1+	
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	15.5	16.0	16.0
	Multi 131	15.0	15.5	16.0	16.5	17.0
	Multi 158	16.0	16.5	17.0	17.5	18.0
	SS Low	-	-	-	-	-
	SS Std	-	15.0	15.5	16.0	16.0
	Frog Legs	-	-	16.5	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	16.5	17.0	17.5
	Multi 131	16.5	17.0	17.5	18.0	18.0
	Multi 158	17.5	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	16.0	16.5	17.0	17.5
	Frog Legs	-	-	17.5	18.0	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	18.0	18.5	19.0
	Multi 131	17.5	18.0	18.5	19.0	19.5
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	17.5	18.0	18.5	19.0
	Frog Legs	-	-	19.0	19.5	-

FSH for 6" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1 1+	
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	16.0	16.5
	Multi 131	-	16.0	16.5	17.0	17.0
	Multi 158	16.5	17.0	17.5	18.0	18.5
	SS Low	-	-	-	-	-
	SS Std	-	-	15.5	16.0	16.5
	Frog Legs	-	-	-	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	17.5	18.0
	Multi 131	-	17.5	18.0	18.5	18.5
	Multi 158	18.0	18.5	19.0	19.5	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	17.0	17.5	18.0
	Frog Legs	-	-	-	18.5	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	19.0	19.5
	Multi 131	-	18.5	19.0	19.5	20.0
	Multi 158	19.0	19.5	20.0	20.5	21.0
	SS Low	-	-	-	-	-
	SS Std	-	-	18.5	19.0	19.5
	Frog Legs	-	-	-	20.0	-

FSH for 7" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1 1+	
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.0	17.5	18.0
	Multi 158	17.0	17.5	18.0	18.5	19.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	18.5	19.0	19.0
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.5	20.0	20.5
	Multi 158	19.5	20.0	20.5	21.0	21.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

FSH for 8" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1 1+	
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.5	18.0	18.5
	Multi 158	-	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.0	19.5	19.5
	Multi 158	-	19.5	20.0	20.5	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	20.0	20.5	21.0
	Multi 158	-	20.5	21.0	21.5	22.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

Lower Leg Length Matrix				
<i>Note for Swing Away Hangers only: subtract 1" from the max achievable LLL for long extension tubes if the Multi Angle Adjustable footplate is selected.</i>				
Hanger	Extension	LLL (inches) (measured from top of seat sling to top of footrest)		Min Allowable FSH (due to 2" ground clearance requirement)
		Min	Max	
60° Swing In/Out	Short	10	14	No lower FSH limit
	Med	11	16	
	Long	12	18	
70° Swing In/Out	Short	10	16	No lower FSH limit
	Med	12	17	
	Long	13	19	
80° Swing In/Out	Short	11	15	No lower FSH limit
	Med	12	17	
	Long	14	19	
ELR	Short	11	16	No lower FSH limit
	Med	11	18	
	Long	13	19	
ALR	Short	11	16	No lower FSH limit
	Med	13	18	
	Long	15	19	
75° Fixed Front	Short	13	18	15.5
	Med	13	20	
	Long	14	20	
85° Fixed Front 0" inset	Short	13	18	15.5
	Med	13	20	
	Long	15	20	
85° Fixed Front 1" inset	Short	13	18	15.5
	Med	15	20	
	Long	17	20	

Back Height Go/No-Go Table																
Back Type	Push Handle Type	Back Height Range														
		10-12			12-14			14-16			16-18			18-20		
		10	11	12	12	13	14	14	15	16	16	17	18	18	19	20
Standard / Angle Adj	Integral	GO														
	Fold Down	NO-GO														
	Ht Adj (reg & JAY)	NO-GO														
	Omit PH	GO														
Lumbar (FWD or REV)	Integral	NO-GO														
	Fold Down	NO-GO														
	Ht Adj (reg & JAY)	NO-GO														
	Omit PH	NO-GO														
Half Fold	Integral	NO-GO														
	Fold Down	NO-GO														
	Ht Adj (reg & JAY)	NO-GO														
	Omit PH	NO-GO														

Back Upholstery Go/No-Go Table				
Push Handle Type	Back Upholstery Type			
	STD EXO (BU95)	TA EXO (BU89 & BU90)	TA Ballistic (BU70)	Omit (BU14)
Integral	GO	GO	GO	GO
Fold Down	NO-GO	GO	GO	GO
Ht Adj reg (PH21)	NO-GO	GO	NO-GO	GO
Ht Adj JAY (PH25)	NO-GO	NO-GO	NO-GO	GO
Omit PH	GO	GO	GO	GO



Sunrise Medical (US) LLC