

# Pediatric Nebulizer Compressor

Item#

## NEB300BEAR

Comes with Disposable Neb Kit, Pediatric Aerosol Mask, & Filters  
HCPCS# E0570 + A7003 + A7015 + A7013

## NEB300BEAR-R

Comes with Disposable and Reusable Neb Kits,  
Pediatric Aerosol Mask, & Filters  
HCPCS# E0570 + A7003 + A7005 + A7015 + A7013

Neb kits include nebulizer bottle, t-piece, mouthpiece,  
6in corrugated reservoir tube, and 7ft supply tube.



Great value

**FAST  
TREATMENT  
TIMES**

## Specifications

**Weight** 2.75 lbs

**Input** 120V/60Hz

**Current** 1.5A

**Medication Capacity** 6ml

**Particle Size** Respirable Fraction 0.5 to 5um 0.5um 80%

**MMAD** 2.0um

**Average Nebulization Rate** 0.3-0.4ml/min

**Noise Level** Below 60dBA

**Compressor Pressure Range** 30-45psi

**Operation Pressure Range** 12-14.5psi (0.85-1.1 bar)

**Liter Flow Range** 8-10LPM

**Operation Temperature Range** 50°F-104°F (10°C-40°C)

**Operation Humidity Range** 10 to 95% RH

**Storage Temperature Range** 14°F-149°F (-10°C-65°C)

**Storage Humidity Range** 10 to 95% RH

Comes in an attractive retail box

