



## ORDER FORM

# Q50 R carbon

EIPW50

Customer	Ship To	
Account #	Name	
Date	Address	
PO#	City	
Buyer	State	Zip
Mark For	Phone	

Q50 R Carbon Totals		\$
Base Price		
Total of additional equipment		
Packaging and transport		
Total Price		

## Base Specifications

<b>User Weight Capacity</b>	300 lbs.
<b>Weight</b> <i>(no batteries)</i>	32 lbs.
<b>Seat Size</b>	Width 17", Depth 17"
<b>Drive Wheel Position</b>	Rear-Wheel Drive
<b>Base Size</b>	Width 22.4", Length 36.3"
<b>Drive Control</b>	nVR2 Brushless Controller
<b>Driving Range</b>	7.5 miles (15 miles with additional battery)
<b>Drive Motors</b>	Brushless motor / 3.7 mph
<b>Back Height</b>	18"
<b>STF Height</b>	19.5"
<b>Drive Wheels</b>	8" with flat-free inserts
<b>Tires / Casters</b>	7" Solid
<b>Charger</b>	2A Off-Board
<b>Battery Size</b>	10Ah Lithium Ion
<b>Color</b>	Carbon Fiber with 5 color accents
<b>Ground Clearance</b>	2.75"

## Additional Included Features and Accessories

<b>Full Carbon Fiber</b>
<b>Quick-Release Folding</b>
<b>Lap Belt</b>
<b>Transport Strap</b>
<b>1.25" Thick Seat Cushion</b>
<b>Underseat Storage Bag</b>
<b>Storage for Additional Battery</b>
<b>Cup Holder</b>
<b>Off-Chair Battery Charging</b>
<b>Rear Running Lights</b>
<b>Depth Adjustable Joystick Mount</b>
<b>Flip-Up Footplate</b>
<b>Flip-Up Armrests</b>
<b>Breathable Back Upholstery</b>

## OPTIONS

OPTION	HCPC	DESCRIPTION										RETAIL
Model												\$
500MW93	K0899	Q50 R Carbon										\$3,667.00
Colors												\$
500C200	—	Black	NC	500C202	—	Blue	NC	500C204	—	Purple	NC	
500C201	—	Red	NC	500C203	—	Yellow	NC	500C205	—	Orange	NC	
Accessories												\$
500AG43	—	Additional 10 Ah Lithium Ion Battery										\$607.00
500BJ36	—	Additional Lithium Ion Battery Charger - <i>Chair comes standard with 1 charger this would be a 2nd charger. Recommended when an extra Lithium Ion Battery is selected, so both batteries can be charged overnight.</i>										\$179.00
5008A123	—	1 Battery Storage Bag - <i>Fits 1 Lithium Ion Battery</i>										\$66.00
5008A125	—	2nd Battery Storage Bag - <i>Fits 1 Lithium Ion Battery</i> <i>If two bags are needed both battery storage bags (1 and 2nd) should be selected</i>										\$66.00
5008A126	—	Joystick Mount for Attendant Control Position - <i>Allows for primary joystick to be relocated to the attendant position</i>										\$50.00
5008A124	—	Q50 R Carbon Complete Storage Bag										\$125.00
JAY Seat Cushions												\$
500CU114	E2601	JAY® Basic - 16" x 17" Cushion										\$109.00
500CU113	E2603	JAY® Go - 16" x 17" Cushion										\$240.00
500CU124	E2622	JAY® X2 - 16" x 17" Cushion										\$501.00

Signature