


LIFT SLINGS

DELUXE SIT-TO-STAND PADDED SLINGS

Features

- Cushioned support panels under the arm and lumbar areas enhance patient comfort
- Polyester with nylon stitching and webbing
- The unique cross-over strap and non-slip inner pad keep the sling in place so that it does not ride up
- Optional stand assist buttocks strap is available for initial lift assistance
- For use with LF2020 STS lift ,
DSLSA4 for use with LF2090 STS lift
- Six (6) month limited warranty

	Size	Torso Girth		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	lb	Kg	
DSLSA1	Small	26 – 43	11.8 – 19.5	400	181.4	1 ea
DSLSA2	Medium	29 – 48	13.2 – 21.8	400	181.4	1 ea
DSLSA3	Large	36 – 59	16.3 – 26.8	400	181.4	1 ea
DSLSA4	XX-Large	50 – 79	22.7 – 35.8	600	272.2	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700075RevD23