

Chemotherapy protection that works as hard as you do



Protecting hospital staff at every step of care

Exposure risks exist in many areas of a healthcare facility, making it critical that staff have the right level of protection. Cardinal Health carries a comprehensive line of chemotherapy gowns that protect healthcare personnel from exposure to bodily fluids and 11 types of chemotherapy drugs, from preparation and handling to administration and disposal.

Safety first

Cardinal Health™ Full Coverage Chemotherapy Gowns provide comparable fluid protection to an AAMI Level 4 surgical gown while also providing a barrier against hazardous chemotherapy drugs.



Meets **USP <800> standards** for handling hazardous drugs



FDA 510(k) clearance*



AAMI Level 4 gowns, including sleeves and seams, meet strict standards for high fluid protection.



The gown meets the broad protection Part A requirements of ASTM F3267, Standard Specification for Protective Clothing for use against Liquid Chemotherapy Drugs.



The gown has been evaluated for **resistance to permeation of 11 chemotherapy drugs** per ASTM F739-20¹.



Reinforced taped seams where direct contact with infectious materials and hazardous drugs will most likely occur.



Flammability Class I per 16 CFR 1610



AAMI Level

Least protective

Most protective

1

2

3

4

This gown is classified as **Level 4** per **AAMI Standard PB70 Liquid barrier performance and classification** of protective apparel and drapes intended for use in healthcare facilities.

Meets USP <800> standards which will be federally enforceable on November 1, 2023.

As cancer treatment evolves, so do the drug touch points in the hospital, from the compounding done in the pharmacy to the different administration sites throughout the facility. Cardinal Health Full Coverage Chemotherapy Gowns meet USP <800> standards.

- Poly-coated material
- Full back style gowns that close and tie in the back
- Taped seams provide added fluid protection
- Knit cuffs

* Not intended for use in the operating room.



Non-sterile:

Full back	
9200M	Medium
9200U	Universal
9200XL	X-Large
9200XXL	XX-Large

A chemotherapy gown that works as hard as you do

Your hospital staff works hard to protect and treat patients. They deserve chemotherapy protection that works hard to protect them, too.

Contact your Cardinal Health sales representative to request product samples.

Prior to use, review all relevant information, including the label supplied with the gown.

*FDA requires gowns with chemotherapy claims and moderate or high fluid level claims, which includes AAMI Level 4, to have 510(k) clearance.

"Premarket Notification Requirements Concerning Gowns Intended for Use in Health Care Settings Guidance for Industry and Food and Drug Administration Staff DECEMBER 2015"

1. **ASTM F739** ASTM International. Standard Test Method for Permeation of Liquids and Gases Through Protective Clothing Materials Under Conditions of Continuous Contact. <https://www.astm.org/Standards/F739.htm>. Accessed October 6, 2021.org

© 2023 Cardinal Health. All Rights Reserved. CARDINAL HEALTH and the Cardinal Health LOGO are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are property of their respective owners. Patent [cardinalhealth.com/patents](https://www.cardinalhealth.com/patents). Lit. No. 2MS21-1703040-03 (10/2023)



[cardinalhealth.com](https://www.cardinalhealth.com)

Cardinal Health
3651 Birchwood Drive
Waukegan, IL 60085