


ROLLATORS

WALKABOUT LITE FOUR-WHEEL ROLLATOR

The Lumex Walkabout Lite Four-Wheel rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 14.5 lbs and supports up to a 300 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Features

- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded straight backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- Durable 6" wheel
- Basket included that can be placed underneath seat or in front for carrying of personal items
- Includes easy to read, retail-friendly user manual
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4300AQ	Aqua	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300B	Royal Blue	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300CH	Champagne	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300G	Teal Green	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300K	Black	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300L	Lavender	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300P	Pink	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4300R	Burgundy	16	40.64	24.5	62.23	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4300AQ	Aqua	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300B	Royal Blue	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300CH	Champagne	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300G	Teal Green	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300K	Black	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300L	Lavender	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300P	Pink	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea
RJ4300R	Burgundy	6	15.24	32	81.28	37	93.98	14.5	6.58	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.

graham
field

GF0900003RevF23