

The Lumex Walkabout Four-Wheel Imperial (Bariatric) is ideal for larger individuals. Similar to the Walkabout Lite, this rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 20 lbs and supports a maximum weight capacity of 450 lb - **EVENLY DISTRIBUTED**. The Walkabout folds quickly and easily into a compact unit for storage and transport.



Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded 20" wide seat for increased seating surface and added comfort
- Removable padded backbar
- Adjustable handle height to accommodate broad range of user heights
- Durable, heavy duty 8" wheel for indoor and outdoor use on even terrain
- Pouch included that can be placed underneath seat for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0149 & E1056
- Includes easy to read, retail-friendly user manual
- Five (5) year limited warranty

And now with added rear reflectors that increase visibility at night!




Straight Backbar

RJ4400K	Black	1 ea
RJ4400R	Burgundy	1 ea

Contoured Backbar

RJ4405B	Royal Blue	1 ea
RJ4405K	Black	1 ea
RJ4405R	Burgundy	1 ea



Color options

Specifications

Width Inside Handgrips	21.5"
Overall Width	28.7"
Overall Depth	26"
<i>measured wheel to wheel</i>	
Seat Width	19.5"
Seat Depth	13.8"
Seat Height	21.8"
Caster Diameter	8"
Minimum Handle Height	34"
Maximum Handle Height	39"
Weight Capacity	450 lb
EVENLY DISTRIBUTED	
Product Weight	20 lb

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.