

McKESSON


**PREVENT® R**

## Push-Button Safety IV Catheter

McKesson Prevent® R push-button safety IV catheters sample blood or administer fluids intravenously. The push-button safety IV catheters feature a blood control active check valve that inhibits blood flow when the valve is activated. Ergonomically designed to address clinicians' needs, the retractable safety catheters contain push-off tabs and contoured molding for a secure grip and to enable one-handed catheter advancement. Clinicians can use one-handed or two-handed techniques. The flashback chamber provides a clear visual confirmation of successful venipuncture. Winged hub options provide a larger surface area for easier application of tape or other catheter securement devices and are ideal for accessing small or fragile veins.

### FEATURES AND BENEFITS

- Push-button style needle retraction safety mechanism
- Blood control active check valve
- Triple facet, back cut, beveled needle
- Gently tapered catheter tip
- Push-off tab
- Flashback chamber allows for clear visual confirmation that catheter has entered the vein
- Audible click lock safety mechanism ensures the needle is safely encapsulated in the catheter body for safety through point of disposal
- Encapsulating needle shield encases the needle when safety mechanism is activated
- Color-coded catheter hub and packaging for easy identification of needle gauge size



### Satisfaction Guaranteed

For complete details, please visit  
[mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2025 McKesson Medical-Surgical Inc.

**McKESSON**



# Push-Button Safety IV Catheter

## FEATURES AND BENEFITS

- Blood control technology
- Sterile
- Non-Pyrogenic
- Radiopaque materials
- DEHP free
- Disposable
- Not made with natural latex rubber



IMAGE	E1 #	MFR#	DESCRIPTION	HUB TYPE	QTY
	<a href="#">1242604</a>	200-IVRB18116	Push-Button Safety IV Catheter Blood Control, 18G x 1.16 in	Straight Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242605</a>	200-IVRB201	Push-Button Safety IV Catheter Blood Control, 20G x 1 in	Straight Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242606</a>	200-IVRB20116	Push-Button Safety IV Catheter Blood Control, 20G x 1.16 in	Straight Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242607</a>	200-IVRB221	Push-Button Safety IV Catheter Blood Control, 22G x 1 in	Straight Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242608</a>	200-IVRB2434	Push-Button Safety IV Catheter Blood Control, 24G x 0.75 in	Straight Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242609</a>	200-IVRWB201	Push-Button Safety IV Catheter Blood Control, 20G x 1 in	Winged Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242610</a>	200-IVRWB221	Push-Button Safety IV Catheter Blood Control, 22G x 1 in	Winged Hub	1 ea 50 ea/bx 4 bx/cs
	<a href="#">1242611</a>	200-IVRWB2434	Push-Button Safety IV Catheter Blood Control, 24G x 0.75 in	Winged Hub	1 ea 50 ea/bx 4 bx/cs