

## AeroFly Ultra-Lite Wheelchair

HCPCS:



- Ultralight magnesium alloy frame weighs just 13.3 lbs. (excluding rear wheels, 19 lbs. with rear wheels)
- Includes both attendant brakes and wheel locks for assisted or independent use
- Fixed, desk length arms with molded armrest pads for comfort and durability
- Lightweight spoke wheels with treaded flat-free tires
- Quick release rear axles
- Fold down backrest
- Seat integrated carry handles
- Removable seat and back padded overlay cushions
- European style swing away footrests with ergonomic release levers
- Ships standard with rear anti-tippers for added safety

PRODUCTS	Description
AF18FDABK-SF	AeroFly,18",FixedDeskArm,SF,Black/Red
AF18FDAGY-SF	AeroFly,18",FixedDeskArm,SF,Grey/Blue
SPECIFICATIONS	Description
Warranty	Limited Lifetime (Frame) / 60 Days (Wearable Parts)
Product Weight	13.3 lbs. (Frame only weight, excludes rear wheels)   19 lbs. (Total weight inc. rear wheels)
Product Dimensions	30.3 x 26.8 x 36.6
Material	Magnesium Alloy
Seat Dimensions	18" x 16"
Weight Capacity	250 lbs.
Seat To Floor Height	19"
Wheels	22"
Carton Shipping Dimensions	31.1 x 13.4 x 30.7
Casters	6"
Folded Dimensions	30.3" x 13" x 29.7"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]