



Raz-AT600

MSRP – USA – USD
Effective May 1, 2025

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AT600 Heavy Duty Mobile Shower Commode Chair with stainless steel frame and components, 25° tilt range, Fixed Back, four dual-locking 6" (150mm) casters, flip-up arm supports, heavy duty foot supports, commode pan, head support with small pad and the proprietary Weight Assist Spring System (WASS). Also included is a 20"H Multi-Position Back with 2" fore/aft adjustment range, two back angles (90°/97°) and tension-adjustable vinyl upholstery. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 605 lbs.**

1. Raz-AT600 codes by size – Must Choose One

Z36022	22"W Seat and Back Frames	\$5,995
Z36024	24"W Seat and Back Frames	\$6,295
Z36026	26"W Seat and Back Frames	\$6,595
Z36028	28"W Seat and Back Frames	\$6,995
Z36030	30"W Seat and Back Frames – Requires custom seat	\$7,595

2. User Weight Range

Z36WR450	WASS User Weight Range 350-450 lbs.	Std	Z36WR600	WASS User Weight Range 400-605 lbs.	\$250
----------	-------------------------------------	-----	----------	-------------------------------------	-------

3. Seats – Must Choose One

Molded and Visco Foam Interface Seats – Seat Support Kit (ZSSKT) required if seat is narrower than seat frame.

ZMS23FRO	Molded 23"W x 18"D – Front Access	Fits 22" and 24"W ¹ frames only	N/C
ZMS26FRO	Molded 26"W x 18"D – Front Access	Fits 26" and 28"W ¹ frames only	N/C
ZVFI23FRO	Visco Foam Interface 23"W x 18"D – Front Access	Fits 22" and 24"W ¹ frames only	\$625
ZVFI26FRO	Visco Foam Interface 26"W x 18"D – Front Access	Fits 26" and 28"W ¹ frames only	\$745
ZSSKTHD	Seat Support Kit for Raz-AT600 – for seats that are narrower than seat frames		\$145

¹ Requires a Seat Support Kit (ZSSKTHD)

Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear.

Select Bridged if an access opening is not required. [See Custom Seat Assessment Form – Resources tab at razdesigninc.com](#)

ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if both width and depth are ≤ 28"	\$1,095
ZCSEATXHD	Custom Dual-Visco Foam Seat – XHD" (enter specs below) Select if width and/or depth is ≥ 28"	\$1,295

Base Frame Width: __ 22" / 24" / 26" / 28" / 30"

Custom Seat Width: __ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than seat frame width)

Seat Depth: __ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" (≥ 21" requires custom longer frame)

Access: __ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged

Aperture Width: __ 6.5" recommended

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom seat widths can be up to 2" wider than seat the frame width to provide additional thigh support without increasing overall width of the chair.
- The Ischial Pelvic Alignment System (IPAS) will allow 2" fore/aft adjustment of the seat position (and therefore seat depth and aperture position) on all seats.
- Access for Rear-Open seats is limited due to rear seat support tube.

4. Symphony Back (Recessed solid back that replaces the standard Tension-Adjustable Vinyl Upholstery [recommended for infection-control environments]). All Symphony Back panels are 13”H.

ZSB22	Symphony Back – 22”W	\$465
ZSB24	Symphony Back – 24”W	\$535
ZSB26	Symphony Back – 26”W	\$505
ZSB28	Symphony Back – 28”W	\$575
ZSB30	Symphony Back – 30”W	\$675

5. Arm Supports Accessories

ZALKHD	Arm Support Locks – HD (pr)	\$95
ZFAPLL	Long Molded Flat Armpad ¹ – Left (14”L x 3.25”W)	\$90
ZFAPLR	Long Molded Flat Armpad ¹ – Right (14”L x 3.25”W)	\$90
ZTRML	Arm Trough ² - Medium Left (18.5” L x 5” W)	\$110
ZTRMR	Arm Trough ² - Medium Right (18.5” L x 5” W)	\$110
ZTRSL	Arm Trough ² - Large Left (21” L x 5.25” W) <i>NEW</i>	\$110
ZTRSR	Arm Trough ² - Large Right (21” L x 5.25” W) <i>NEW</i>	\$110
ZPAML	Pivoting Arm Mount ³ – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$165
ZPAMR	Pivoting Arm Mount ³ – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	\$165
ZAPSB_*	Anterior Posture Support Bar ⁴ * Enter frame width to complete code (22” / 24”)	\$245
ZAPSB_*	Anterior Posture Support Bar ⁴ * Enter frame width to complete code (26” / 28” / 30”)	\$275

¹May require ZPAM for Arm Support to flip up with Flat Arm Pad

² ZPAM is required for Arm Support to flip up with Arm Trough

³ Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered

⁴ Not compatible with ZFAP / ZTR

6. Leg Supports

ZFSHDL	Foot Support Left – 16”-20.5”	Standard
ZFSHDR	Foot Support Right – 16”-20.5”	Standard
ZFET2L	Footrest Extension Tube ¹ Left – 2” longer	\$25
ZFET2R	Footrest Extension Tube ¹ Right – 2” longer	\$25
ZADML	Adj Angle/Depth Molded Footplate - Left	\$60
ZADMR	Adj Angle/Depth Molded Footplate - Right	\$60
ZFET2AL	Footrest Extension Tube Left for ZADML – 2” longer	\$25
ZFET2AR	Footrest Extension Tube Right for ZADMR – 2” longer	\$25
ZFBT	Flip-Back Foot Support (1-piece) 16.5”-20” with calf pads	\$455
ZELRL	Elevating Leg Support – Left (includes ZADML)	\$495
ZELRR	Elevating Leg Support – Right (includes ZADMR)	\$495
ZRLSL	Universal Residual Limb Support – Left	\$195
ZRLSR	Universal Residual Limb Support – Right	\$195

¹ Not compatible with ZADM / ZVS / ZFBT

7. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$95
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (22" / 24")	\$105
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (26" / 28" / 30")	\$245
ZCSIC_*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (22" / 24" / 26" / 28" / 30")	\$125
ZHS2224	H-Strap – for 22" / 24"W frames – hook and loop attachment	\$110
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$60
ZNC	Neoprene Cover for Standard Footplate (pr)	\$90
ZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates	\$90
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$95
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$95
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$95

¹ Not compatible with ZFBT / ZVS

8. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

Head Supports		
ZHSP2	Head Support Pad – Large (12" x 6")	\$50
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$75
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate.		
ZLTS2LHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Left	\$225
ZLTS2RHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Right	\$225
Leg Adductors		
ZLGT2L	Leg Adductor w/ Large Pad – Left	\$225
ZLGT2R	Leg Adductor w/ Large Pad – Right	\$225
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors		
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
ZLAEPL	Lateral/Adductor Extension Plate - Left (ea)	\$28
ZLAEPR	Lateral/Adductor Extension Plate - Right (ea)	\$28
Lateral Thoracic Support	Left Right	
Leg Adductor	Left Right	
Hip Adductor	Left Right	
Total Quantity ___		
Medial Thigh Supports (Pommel) Large Pad Size – 6.75" x 4.3"		
ZMTST2	Medial Thigh Support w/ Large Pads	\$495

9. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBXLHDT	Pelvic Belt – XLarge (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$75	ZCBXL	Chest Belt – XLarge ² (back frame widths 24"-26", 1-piece w/ hook and loop closure)	\$55
ZPBXXLHDT	Pelvic Belt – XXL Large (seat frame widths 28"-30", 2-piece w/ side-release buckle)	\$110	ZCBXXL	Chest Belt – XXL Large ² (back frame widths 26"-30", 1-piece w/ hook and loop closure)	\$65
ZPBICXLHDT	Infection Control Pelvic Belt – XLarge ¹ (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$125	ZCBICXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 24"-26", 1-piece w/ side-release buckle)	\$95
ZPBICXXLHDT	Infection Control Pelvic Belt – XX-Large ¹ (seat frame widths 26"-30", 2-piece w/ side-release buckle)	\$135	ZCBICXXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 26"-30", 1-piece w/ side-release buckle)	\$105
ZPBBPXLHDT	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$235	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	\$195
Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPXLHDT	22"-30"W base frame		
	Chest belt	ZCBBPXL	22"-30"W back frame		
¹ Non-absorbent belts – recommended for infection-control environments. ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.					

10. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC23HD	Splash Curtain ¹ for 23"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$345
ZSC26HD	Splash Curtain ¹ for 26"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$375
ZSCCHD	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$365
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard urine drainage bag holder for 22" / 24" / 26" / 28" / 30"W frames	\$155
¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

11. Other

ZDCLT	Directional Caster Lock – fits on left rear caster (facilitates steering)	\$165	ZPSK1	Commode Pan Spacer Kit ² – 1” lower	\$200
ZATTRR	Anti-Tippers – Rear (pr)	\$345	ZPSK2	Commode Pan Spacer Kit ² – 2” lower	\$215
ZATTFR	Anti-Tippers – Front (pr)	\$345	ZPSK3	Commode Pan Spacer Kit ² – 3” lower	\$230
ZTLS	Tilt Limit Spacers – 1/2” spacers snap onto gas struts to limit tilt (1/2” reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C	ZPSK4	Commode Pan Spacer Kit ² – 4” lower	\$245
ZSPKTHD24	AT600 Self-Propel Kit – 24” wheels – includes rear anti-tippers	\$1,425	¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.		

12. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	-\$45	Z912	Omit Complete Head Support Assembly	-\$50

Some assembly of Raz products is required.

13. Notes:**Specifications****Raz-AT600 (Attendant Tilt)**

Weight Capacity	605 lbs
Seat Depth (with Standard Molded Seat)	18”-20”
Overall Width	Seat frame width plus 4.5” (with Standard Arm Support Pads)
Toilet Clearance Width	17”
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	42”
- at 25°	52”
Clearance Height (under seat)	with 6” casters
- with pan brackets	17.75”-18.75”
- without pan brackets	18.5”-19.5”
Arm Support Height - above molded seat	9.25”
Seat-to-Floor Height (with molded seat)	22.25”
Back Height (Molded Seat to top of upholstery)	20”
Product Weight (with 22”W seat frame)	84 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)