


TRANSFER BENCHES

PADDED KNOCK DOWN TRANSFER BENCH

Designed to provide assistance to those having difficulty stepping over bathtub walls.

Features

- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Side support arm aids in lateral transfer
- Seat height is adjustable in 1" increments
- Reversible to accommodate any bathroom
- Knock-down design makes transporting easy
- Quick and easy tool-free assembly
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED
- Three (3) month limited warranty



| | | overall Width | | overall Depth | | Seat Width | | Seat Depth | | minimum Seat Height | | maximum Seat Height | | Commode diameter | | Weight Capacity EVENLY DISTRIBUTED | | UOM |
|----------|-------------------------------|------------------|-------|------------------|-------|------------|-------|------------|-------|------------------------|-------|------------------------|-------|---------------------|-------|------------------------------------------|--------|------|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | Kg | |
| 7955KD-1 | padded transfer bench | 29 | 73.66 | 23 | 58.42 | 24 | 60.96 | 16 | 40.64 | 18 | 45.72 | 22 | 55.88 | - | - | 300 | 136.08 | 1 ea |
| 7955KD-2 | padded transfer bench | 29 | 73.66 | 23 | 58.42 | 24 | 60.96 | 16 | 40.64 | 18 | 45.72 | 22 | 55.88 | - | - | 300 | 136.08 | 2/cs |
| 7956KD-1 | padded commode transfer bench | 29 | 73.66 | 23 | 58.42 | 24 | 60.96 | 16 | 40.64 | 18 | 45.72 | 22 | 55.88 | 8 | 20.32 | 300 | 136.08 | 1 ea |

