


Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Technical Data Overview			
Product Width	14" to 20"	Base	Pre-contoured closed cell foam
Product Depths	14" to 20"	Insert	Multizone elastomer insert with 2 standard pressure adjustment pods
Product Height	2"	Cover	Double spacer mesh outer cover with durable side bands and non-slip base
Product Weight	2 lbs.		
User Weight Capacity	330 lbs.	JAY Your Way Available	Yes
HCPCS Code	E2622		

Measuring and Sizing Guide



JAY LitePro Cushion									
✓	OPTION	HCPCS	DESCRIPTION	RETAIL	✓	OPTION	HCPCS	DESCRIPTION	RETAIL
<input type="checkbox"/>	JLP1414	E2622	JAY LitePro 14" x 14"	\$525	<input type="checkbox"/>	JLP1719	E2622	JAY LitePro 17" x 19"	\$525
<input type="checkbox"/>	JLP1415	E2622	JAY LitePro 14" x 15"	\$525	<input type="checkbox"/>	JLP1816	E2622	JAY LitePro 18" x 16"	\$525
<input type="checkbox"/>	JLP1416	E2622	JAY LitePro 14" x 16"	\$525	<input type="checkbox"/>	JLP1817	E2622	JAY LitePro 18" x 17"	\$525
<input type="checkbox"/>	JLP1417	E2622	JAY LitePro 14" x 17"	\$525	<input type="checkbox"/>	JLP1818	E2622	JAY LitePro 18" x 18"	\$525
<input type="checkbox"/>	JLP1515	E2622	JAY LitePro 15" x 15"	\$525	<input type="checkbox"/>	JLP1819	E2622	JAY LitePro 18" x 19"	\$525
<input type="checkbox"/>	JLP1516	E2622	JAY LitePro 15" x 16"	\$525	<input type="checkbox"/>	JLP1820	E2622	JAY LitePro 18" x 20"	\$525
<input type="checkbox"/>	JLP1517	E2622	JAY LitePro 15" x 17"	\$525	<input type="checkbox"/>	JLP1916	E2622	JAY LitePro 19" x 16"	\$525
<input type="checkbox"/>	JLP1518	E2622	JAY LitePro 15" x 18"	\$525	<input type="checkbox"/>	JLP1917	E2622	JAY LitePro 19" x 17"	\$525
<input type="checkbox"/>	JLP1615	E2622	JAY LitePro 16" x 15"	\$525	<input type="checkbox"/>	JLP1918	E2622	JAY LitePro 19" x 18"	\$525
<input type="checkbox"/>	JLP1616	E2622	JAY LitePro 16" x 16"	\$525	<input type="checkbox"/>	JLP1919	E2622	JAY LitePro 19" x 19"	\$525
<input type="checkbox"/>	JLP1617	E2622	JAY LitePro 16" x 17"	\$525	<input type="checkbox"/>	JLP1920	E2622	JAY LitePro 19" x 20"	\$525
<input type="checkbox"/>	JLP1618	E2622	JAY LitePro 16" x 18"	\$525	<input type="checkbox"/>	JLP2017	E2622	JAY LitePro 20" x 17"	\$525
<input type="checkbox"/>	JLP1716	E2622	JAY LitePro 17" x 16"	\$525	<input type="checkbox"/>	JLP2018	E2622	JAY LitePro 20" x 18"	\$525
<input type="checkbox"/>	JLP1717	E2622	JAY LitePro 17" x 17"	\$525	<input type="checkbox"/>	JLP2019	E2622	JAY LitePro 20" x 19"	\$525
<input type="checkbox"/>	JLP1718	E2622	JAY LitePro 17" x 18"	\$525	<input type="checkbox"/>	JLP2020	E2622	JAY LitePro 20" x 20"	\$525

JAY Pressure Adjustment Pods

Standard pods (pair) included in base cushion selection above - only select JLPSTD if extra required. Firm pods (pair) not included in base cushion selection - select JLPFIRM if required.

✓	OPTION	HCPCS	DESCRIPTION	RETAIL	✓	OPTION	HCPCS	DESCRIPTION	RETAIL
<input type="checkbox"/>	JLPSTD		JAY LitePro Standard Pods (Pair)	\$50	<input type="checkbox"/>	JLPFIRM		JAY LitePro Firm Pods (Pair)	\$50

JAY Your Way Modifications

Add an "M" in front of the item number to order a modification through JAY Your Way.

OPTION	DESCRIPTION	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		RETAIL
IW	Increase Width	<input type="checkbox"/>	+1"	<input type="checkbox"/>	+2"	N/C
IL	Increase Length	<input type="checkbox"/>	+1"	<input type="checkbox"/>	+2"	N/C

<input checked="" type="checkbox"/>	OPTION	DESCRIPTION	RETAIL
<input type="checkbox"/>	CVR	Reverse Dartex Cover	N/C

JAY LitePro Product Images

