

Fusion**ONE**™

Power Hi-Lo, Manual Back Exam Chair with Foot Control

For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE Exam Chair offers the features healthcare providers need most — all in one chair!

This hi-lo exam chair enables simple ADA access and ease of examination for patients of all ages and abilities.

Model 3001



Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › Low-access ADA height 17"-19" (43.2cm-48.3cm)
- › Large, reversible side drawer and large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Durable, galvanized steel construction
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

Specifications

- › Height: min. 19" (±.5") 48.3 cm (±1.27 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front) 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side) 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International model: 3001-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76 cm x 163 cm x 79 cm)

Accessories

- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Model 3001
shown with
options

Colors

Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
True Graphite	Classic Black	Fine Wine	Blue Skies
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 06.20 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1370LT Rev A.