

COVID-19/FLU A&B ANTIGEN HOME TEST

The Consult COVID-19/Flu A&B Antigen Home Test is a lateral flow immunochromatographic assay intended for the qualitative detection and differentiation of influenza A and influenza B nucleoprotein antigens and SARS-CoV-2 nucleocapsid antigen directly in anterior nasal swab samples from individuals with signs and symptoms of respiratory tract infection. Symptoms of respiratory infections due to SARS-CoV-2 and influenza can be similar. This test is for non-prescription home use by individuals aged 14 years or older testing themselves, or adults testing individuals aged 2 years or older.

FEATURES AND BENEFITS

- Simultaneously detects 3 viruses (SARS-CoV-2, Flu A&B)
- Results in 15 minutes
- No repeat testing necessary
- Two result windows for easy readability
- Shallow nasal swab sample collection
- Suitable for 2+ years old*
- CLIA waived



*Note: Swabbing should be performed by an adult for children aged 2 to 13. User must be aged 14 + to perform self-test.

[Interpreting Test Results](#)

E1 #	MFR#	DESCRIPTION	CONTENTS	QTY
1267919	535-COVFLUOTC	COVID-19/Flu A&B Antigen Home Test	2 Sealed Test Cassettes, 2 Sterile Nasal Swabs, 2 Pre-filled Extraction Tubes, 2 Extraction Tube Tips, 1 Quick Reference Instructions (QRI)	2 test/bx 300 bx/cs

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.