

AeroWalk™ Ultra-lite Rollator

Elevate Every Step with the AeroWalk™ Ultra-lite Rollator



Refined Style, Designed to Move with You

Refined Style, Designed to Move with You
Elegant colors & woodgrain accents highlight the seamless, weld-free frame.



Ultra-Lightweight

Just 14.5 pounds with a compact side-to-side fold



Optimal Support

Extendable, adjustable-height backrest provides lumbar support at rest



Smooth Mobility

10" front wheels & 8" back wheels for indoor-outdoor use



Customized Comfort

Height-adjustable handles & ergonomic hand grips for all-day use



Convenient Storage

Stylish, removable bag keeps essentials close

Lightweight mobility meets modern design with the AeroWalk™ Ultra-lite Rollator. Combining a seamless, weld-free frame and innovative back support, it sets the new benchmark for portable performance at just 14.5 pounds!

Item #	Description	UPC	UOM
10277BL	AeroWalk Ultra-lite Rollator, Blue	822383027104	1/ea
10277GY	AeroWalk Ultra-lite Rollator, Gray	822383027098	1/ea
10277WT	AeroWalk Ultra-lite Rollator, White	822383027111	1/ea

Product Specs

Material	Aluminum	Product Weight	14.5 lbs
Wheel Size	10" front 8" back	Shipping Weight	19.5 lbs
Handle Height	31.5" - 37.25"	Weight Capacity	330 lbs
Unit Dimensions	23.25" (W) x 26.5" (D)	HCPC	E0143 + E0156
Seat Dimensions	17.75" (L) x 10" (W) x 23" (H)	Warranty	Limited Lifetime



Easy Side-to-Side Folding

Allows for effortless travel and storage