



N95 PARTICULATE FILTER RESPIRATOR AND SURGICAL MASK REGULAR

Item No. RP88020



SUBMICRON PARTICLE FILTRATION (PFE)	> 99% @ 0.1 micron
FLUID RESISTANCE	ASTM Level 3 (160 mm Hg)
NIOSH N95 APPROVAL #	TC-84A-5216
FDA 510(k) SURGICAL MASK APPROVAL #	K102092
MASK TYPE	N95 Respirator and Surgical Mask
MASK SIZE	Regular
MASK COLOR	White
MASK MATERIAL	Spunbond Polypropylene (SBPP)
HEADBAND MATERIAL	Nylon/Spandex Elastic
HYPOALLERGENIC STATUS	Non-irritant
LATEX CONTENT	Not made with natural rubber latex
ITEMS PER BOX	50 respirators
BOXES PER CASE	6 boxes
ITEMS PER CASE	300 respirators
CASE DIMENSIONS	19 ³ / ₈ " x 11 ¹ / ₄ " x 6 ⁵ / ₈ "
CASE WEIGHT	6 lbs
COUNTRY OF ORIGIN	USA