


CANDo®

Hi-Lo Treatment Table, 3 Section, Black Vinyl Top

- **Fully Adjustable for Maximum Comfort** - Head section tilts from -55° to +45° and leg section raises up to 75° to find your perfect position
- **Easy Electric Height Adjustment** - Convenient foot bar raises and lowers the table at the touch of your foot
- **Safe, Rounded Design** - Smooth, rounded frame legs eliminate sharp edges for added safety
- **Built-In Face Cradle** - Removable nose/face cushion and plug provide comfort when lying face-down
- **Clean, Modern Look** - Enclosed base hides all mechanical components for a sleek appearance
- **Stable on Any Surface** - Adjustable leveling feet keep the table steady, while quiet-glide 3" (7.6 cm) casters allow easy repositioning
- **Powerful, Quiet Motor** - Smoothly adjusts from lowest to highest position in just 30 seconds
- **Wide Height Range** - 20-inch (50.5 cm) adjustment range accommodates various treatment heights and client needs
- **Premium 2-Inch (5 cm) Cushioning** - High-density foam padding (79.4 lb/m³ / 36 kg/m³) provides lasting comfort and support
- **Durable, Easy-Care Upholstery** - Vinyl top resists stains, moisture, and everyday wear
- **Heavy-Duty Construction** - Safely supports a weight capacity of up to 550 lb (250 kg)
- **Quality Certified** - Meets FDA and CE safety standards for your peace of mind

Technical Data



Angle of head section	Angle of leg section	Wheel	Product weight	Packing size
-55°~+45°	0~+75°	φ3"	159lbs	74.8"x28"x19.7"
-55°~+45°	0~+75°	φ75mm	72kgs	190x71x50cm

15-4295


CANDo® Hi-Lo Treatment Table, 3 Section,
Black Vinyl Top, White Base, Nose Hole Cut Out



3" (7.6 cm)
silent casters,
adjustable feet



nose hole cut out



360° foot bar for
table height adjustments



FEI
FABRICATION
ENTERPRISES INC

Fabrication Enterprises Inc.