

Protekt® Aire 4000DX

Low Air Loss/Alternating Pressure Mattress System
with Digital Pump

HCPCS
E0277



FEATURES AND BENEFITS:

PROTEKT® AIRE 4000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates airflow may be compromised.

MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Static head section prevents unnecessary head/neck movement.
- Optional raised side rails.
- Weight Capacity: 400 lbs.
- Warranty: 2 Year Non-Prorated.

80040DX	Protekt® Aire 4000DX 36"x80"x8"
80040DXRR	Protekt® Aire 4000DX with Raised Rails

